



Interlaboratory Evaluation of an Off-Line SFE/IR Method for Determination of Petroleum Hydrocarbons in Solid Matrices

Project Report



INTERLABORATORY EVALUATION OF AN OFF-LINE SFE/IR METHOD FOR DETERMINATION OF PETROLEUM HYDROCARBONS IN SOLID MATRICES

By

Viorica Lopez-Avila, Richard Young, and Robert Kim
Midwest Research Institute
California Operations
Mountain View, California 94043

Contract Number 68-C1-0029

Project Officer: Werner F. Beckert
Quality Assurance and Methods Development Division
Environmental Monitoring Systems Laboratory
Las Vegas, Nevada 89119

ENVIRONMENTAL MONITORING SYSTEMS LABORATORY
OFFICE OF RESEARCH AND DEVELOPMENT
U.S. ENVIRONMENTAL PROTECTION AGENCY
LAS VEGAS, NEVADA 89119

NOTICE

The information in this document is the result of work funded wholly or in part by the United States Environmental Protection Agency under Contract #68-C1-0029 to Midwest Research Institute. It has been subjected to the Agency's peer and administrative review, and it has been approved for publication as an EPA document.

PREFACE

This is the final report for Work Assignment 2-7, EPA Contract 68-C1-0029, conducted at Midwest Research Institute, California Operations. The project was directed by Dr. Viorica Lopez-Avila.

This report was written by Dr. Viorica Lopez-Avila. Technical support for the project was provided by Richard Young and Robert Kim.

The following collaborators/laboratories (listed in alphabetical order) participated in this study:

Mark Bruce, Wadsworth/Alert Laboratories, North Canton, OH

Paul David, Amoco Corporation, Naperville, IL

Sally Eckert-Tilotta, Energy and Environmental Research Center, University of North Dakota, Grand Forks, ND

Anna P. Emergy, NIST, Gaithersburg, MD

Chuck Hecht, Chemical Waste Management, Riverdale, IL

Joe Hedrick, Hewlett-Packard Company, Avondale, PA

Joe Levy, Suprex Corporation, Pittsburgh, PA

Shirley Liebman, CCS Instruments, West Grove, PA

Mary Ellen McNally/Connie Deardorff, DuPont Agricultural Products, Wilmington, DE

Midwest Research Institute, California Operations, Mountain View, CA

Nathan Porter, Dionex-Lee Scientific, Salt Lake City, UT

Steve Pyle, EPA, EMSL-LV, Las Vegas, NV

John L. Snyder, Lancaster Laboratories, Lancaster, PA

Joe Tehrani, Isco, Inc., Lincoln, NE

We gratefully acknowledge their contributions.

ABSTRACT

A collaborative study was conducted, with Midwest Research Institute, California Operations (MRI-CO), as the lead laboratory, and with 13 additional volunteer laboratories participating, to determine the method accuracy and precision of EPA SW-846 draft Method 3560. This method describes the extraction of petroleum hydrocarbons from solid matrices with supercritical carbon dioxide at 340 atm and 80°C for 30 min (dynamic; flow rate 1 to 2 mL/min as liquid) and collection of the extracted materials in 3 mL tetrachloroethylene. The extracts generated by the participating laboratories were shipped to MRI-CO for analysis by infrared spectrometry (draft Method 8440). The collaborative study was based on the AOAC International blind replicate design with balanced replicates. Four soil samples were to be extracted in triplicate. Three of these samples were standard reference materials with TPH contents of 614, 2050, and 32,600 mg/kg, respectively. The fourth sample was a clay soil spiked with motor oil at 10,000 mg/kg. Each of the four samples was extracted by the laboratories in triplicate with supercritical carbon dioxide, and the extracts were analyzed by MRI-CO on two different infrared spectrometers. In addition, each of the participating laboratories extracted a sample of the unspiked clay soil, clay soil spiked with both corn oil and reference oil (a mixture of n-hexadecane, isooctane, and chlorobenzene) at 1,000 mg/kg each, and clay soil spiked with motor oil at 10,000 mg/kg and brought to a water content of 30 percent (the individual laboratories then had to add anhydrous sodium sulfate). These latter three samples were extracted only once. They were included in the study to address possible cross-contamination in the supercritical fluid extraction system (the laboratories were asked to extract the unspiked clay soil right after extraction of a clay soil sample spiked at 10,000 mg/kg with motor oil), the presence of interferences such as oily materials, and the effect of water content on extraction efficiency.

After outlier removal (using both the Cochran and the Grubbs tests), we calculated the mean concentration, repeatability standard deviation, reproducibility standard deviation, repeatability relative standard deviation, and reproducibility relative standard deviation for each of the four matrices extracted in triplicate. In addition, the relationship between the measured petroleum hydrocarbon contents of these samples and the true concentrations was established.

The results from the triplicate analyses of the four materials show that the method accuracy (percent recovery) for petroleum hydrocarbon concentrations ranging from 614 to 32,600 mg/kg was 83 percent; the mean recoveries of petroleum hydrocarbons from each of the four matrices ranged from 78 to 107 percent for the analyses performed with a Perkin-Elmer FTIR spectrometer, and from 76 to 101 percent for the analyses performed with a Buck-Scientific IR spectrometer. The differences of the results from the two instruments on a sample-by-sample basis were less than 17 percent for the total petroleum hydrocarbon determinations. The interlaboratory method precisions ranged from 17 to 45 percent for analyses on the Perkin-Elmer FTIR system and from 17 to 48 percent for analyses on the Buck-Scientific IR system; the intralaboratory method precisions ranged from 12 to 17 percent for analyses on the Perkin-Elmer FTIR system and from 11 to 18 percent for analyses on the Buck-Scientific IR system. Method accuracy and precision data are presented for the five laboratories that used the Isco supercritical fluid extraction systems and for the seven laboratories that used different supercritical fluid extraction systems, all having extraction vessels of 3.5-mL volume or less.

Analysis of variance (ANOVA) was used to separate (mathematically) the total variation of the experimental measurements into an intralaboratory portion and an interlaboratory portion, with a corresponding split of the total number of degrees of freedom. The ANOVA results indicate that, as expected, the variation from laboratory to laboratory was greater than that attributed to the analytical error displayed within laboratories. The matrix and operational parameters, such as flow rate, size of the extraction vessel, extraction vessel design and orientation, mode of collection of the extracted material, and temperature of the collection solvent/trap all seemed to be important.

The recoveries from the wet clay soil samples mixed with sodium sulfate were above 30 percent for only two laboratories, and eight laboratories had recoveries below 8 percent. Separate experiments conducted at MRI-CO showed that addition of anhydrous magnesium sulfate and Hydromatrix (diatomaceous earth) raised the recoveries to 72 percent or higher for samples containing not more than 20 percent water. When samples contained more than 20 percent water, they had to be mixed with anhydrous magnesium sulfate and allowed to equilibrate for several hours (preferably overnight in sealed containers and at 4°C to minimize losses of volatile petroleum hydrocarbons) before extraction with supercritical carbon dioxide.

CONTENTS

| | |
|---|-----|
| Notice | ii |
| Preface | iii |
| Abstract | iv |
| Tables | vii |
| 1. Introduction | 1 |
| 2. Conclusions | 2 |
| 3. Recommendations | 3 |
| 4. Experimental | 4 |
| Collaborative study | 4 |
| Treatment of data | 17 |
| Quality assurance/quality control | 18 |
| 5. Results and discussion | 20 |
| Interlaboratory method performance | 20 |
| Quality assurance/quality control | 55 |
| References | 72 |
| Appendices | |
| A. Proposed draft protocol for supercritical fluid extraction of petroleum hydrocarbons (Method 3560) | A-1 |
| B. Proposed draft protocol for determination of total recoverable petroleum hydrocarbons by infrared spectrophotometry (Method 8440) | B-1 |
| C. List of instructions for collaborators | C-1 |

TABLES

| Number | | Page |
|--------|--|------|
| 1 | Samples used in the interlaboratory study | 5 |
| 2 | Identification of the SFE systems used in the collaborative study | 6 |
| 3 | SFE operating conditions for the Isco SFE systems | 7 |
| 4 | SFE operating conditions for the Dionex-Lee Scientific SFE-703 systems | 8 |
| 5 | SFE operating conditions for the Suprex Prepmaster SFE systems | 9 |
| 6 | SFE operating conditions for the Suprex SFE-50 system | 10 |
| 7 | SFE operating conditions for the CCS Instruments SFE system | 11 |
| 8 | SFE operating conditions for the HP 7680A SFE systems | 12 |
| 9 | Flow rate of carbon dioxide used by the participating laboratories | 13 |
| 10 | Exact weights of the samples extracted in the study | 15 |
| 11 | Final volumes of the extracts before dilution | 16 |
| 12 | Concentrations of TPHs determined in TPH-1 soil extracts using the PE-FTIR spectrometer | 21 |
| 13 | Concentrations of TPHs determined in TPH-1 soil extracts using the BSci-IR spectrometer | 22 |
| 14 | Concentrations of TPHs determined in TPH-2 soil extracts using the PE-FTIR spectrometer | 23 |
| 15 | Concentrations of TPHs determined in TPH-2 soil extracts using the BSci-IR spectrometer | 24 |
| 16 | Concentrations of TPHs determined in SRS103-100 soil extracts using the PE-FTIR spectrometer | 25 |

TABLES (continued)

| Number | | Page |
|--------|---|------|
| 17 | Concentrations of TPHs determined in SRS103-100 soil extracts using the BSci-IR spectrometer | 26 |
| 18 | Concentrations of TPHs determined in spiked clay soil extracts using the PE-FTIR spectrometer | 27 |
| 19 | Concentrations of TPHs determined in spiked clay soil extracts using the BSci-IR spectrometer | 28 |
| 20 | Summary of results from all laboratories (before outlier removal) | 29 |
| 21 | Summary of results from all laboratories (after outlier removal) | 30 |
| 22 | Summary of results from laboratories using Isco systems (before outlier removal) | 31 |
| 23 | Summary of results from laboratories using Isco systems (after outlier removal) | 32 |
| 24 | Summary of results from laboratories using extraction vessels 3.5 mL or less in volume (before outlier removal) | 33 |
| 25 | Summary of results from laboratories using extraction vessels 3.5 mL or less in volume (after outlier removal) | 34 |
| 26 | Outlier testing—the Cochran statistic, the single Grubbs statistic, and the double Grubbs statistic | 35 |
| 27 | One-way ANOVA for the TPH-1 soil samples extracted by all laboratories | 38 |
| 28 | One-way ANOVA for the TPH-2 soil samples extracted by all laboratories | 39 |
| 29 | One-way ANOVA for the SRS103-100 soil samples extracted by all laboratories | 40 |
| 30 | One-way ANOVA for the spiked clay soil samples extracted by all laboratories | 41 |
| 31 | One-way ANOVA for the TPH-1 soil samples extracted by laboratories 05, 13, 14, 15, and 17 | 42 |
| 32 | One-way ANOVA for the TPH-2 soil samples extracted by laboratories 05, 13, 14, 15, and 17 | 43 |

TABLES (continued)

| Number | | Page |
|--------|---|------|
| 33 | One-way ANOVA for the SRS103-100 soil samples extracted by laboratories 05, 13, 14, 15, and 17 | 44 |
| 34 | One-way ANOVA for the spiked clay soil samples extracted by laboratories 05, 13, 14, 15, and 17 | 45 |
| 35 | One-way ANOVA for the TPH-1 soil samples extracted by laboratories 01, 04, 06, 08, 13, 14, and 17 | 46 |
| 36 | One-way ANOVA for the TPH-2 soil samples extracted by laboratories 01, 04, 06, 08, 13, 14, and 17 | 47 |
| 37 | One-way ANOVA for the SRS103-100 soil samples extracted by laboratories 01, 04, 06, 08, 13, and 14 | 48 |
| 38 | One-way ANOVA for the spiked clay soil samples extracted by laboratories 01, 04, 06, 08, 13, 14, and 17 | 49 |
| 39 | Recoveries of oil/grease and TPHs determined in extracts of clay soil spiked with corn oil and reference oil using the PE-FTIR spectrometer | 51 |
| 40 | Recoveries of oil/grease and TPHs determined in extracts of clay soil spiked with corn oil and reference oil using the BSci-IR spectrometer | 52 |
| 41 | Concentrations and percent recoveries of TPHs determined in extracts of wet clay soil spiked with motor oil using the PE-FTIR spectrometer | 53 |
| 42 | Concentrations and percent recoveries of TPHs determined in extracts of wet clay soil spiked with motor oil using the BSci-IR spectrometer | 54 |
| 43 | Percent recoveries of TPHs determined in extracts of wet clay soil spiked with motor oil using the PE-FTIR spectrometer | 56 |
| 44 | Percent recoveries of TPHs determined in extracts of wet clay soil spiked with motor oil or diesel oil using the PE-FTIR spectrometer | 56 |
| 45 | Correlation between the PE-FTIR data and the BSci-IR data | 57 |
| 46 | Calibration data for the PE-FTIR spectrometer using reference oil in PCE and the 10-mm path-length IR cell | 58 |

TABLES (concluded)

| Number | | Page |
|--------|--|------|
| 47 | Calibration data for the BSci-IR spectrometer using reference oil in PCE and the 10-mm path-length IR cell | 60 |
| 48 | Calibration data for the PE-FTIR spectrometer using motor oil in PCE and the 10-mm path-length IR cell | 62 |
| 49 | Calibration data for the BSci-IR spectrometer using motor oil in PCE and the 10-mm path-length IR cell | 63 |
| 50 | Calibration data for the PE-FTIR spectrometer using reference oil in PCE and the 10-mm path-length IR cell (no silica gel cleanup) | 64 |
| 51 | Calibration data for the BSci-IR spectrometer using reference oil in PCE and the 10-mm path-length IR cell (no silica gel cleanup) | 65 |
| 52 | Concentrations of daily standards analyzed to verify system reproducibility | 66 |
| 53 | Concentrations of TPHs in the extracts from the unspiked clay soil samples | 67 |
| 54 | Concentrations of TPHs in the PCE blanks analyzed during this study | 68 |
| 55 | Concentrations of TPHs in the extracts submitted as system blanks | 69 |
| 56 | Percent recoveries of TPHs from the spiked clay soil samples, stored at 4°C in the dark, as a function of time | 69 |

SECTION 1

INTRODUCTION

Widespread use of chlorofluorocarbons and other harmful organic solvents in environmental testing laboratories is of major concern to the U.S. Environmental Protection Agency (EPA). Therefore, EPA's Environmental Monitoring Systems Laboratory in Las Vegas, NV, started the development of supercritical fluid extraction (SFE) methods that replace those harmful organic solvents with supercritical carbon dioxide or carbon dioxide modified with small amounts of methanol or other innocuous organic solvents as extractants. The first SFE method, SW-846 Method 3560, was recently proposed (1, 2). The proposed method describes extraction of petroleum hydrocarbons with supercritical carbon dioxide at 340 atm and 80°C for 30 min (dynamic), with a carbon dioxide flow rate of 1 to 2 mL/min as liquid, or 500 to 1,000 mL/min as decompressed gas. Depending on the extraction system used, the stream of carbon dioxide containing the extracted materials passes either through a collection vial containing 3 mL tetrachloroethylene (also known as perchloroethylene, PCE), or through an adsorbent trap. The technique is simple, requiring approximately 30 min for the extraction of a 2- to 5-g solid sample, and 10 min for extract cleanup (to remove polar organic compounds) and analysis by infrared (IR) spectrometry. The results of a single-laboratory evaluation of this SFE method indicated that its performance is equivalent to the Soxhlet extraction method using Freon-113 as extractant, its accuracy is 80 percent or better, and its precision is ± 20 percent (1, 2).

A collaborative study of Method 3560 was conducted by 14 laboratories, with Midwest Research Institute, California Operations (MRI-CO), as the lead laboratory, and with 13 additional volunteer laboratories participating. The goal of this study was to determine the method accuracy and precision of draft Method 3560 (Supercritical Fluid Extraction of Petroleum Hydrocarbons) when used at different laboratories using instrumentation from different manufacturers.

The criteria for selecting the laboratories included availability to the laboratories of commercial SFE systems that could be clearly described by the participants, willingness to perform the extractions within a month after sample receipt, previous experience of the prospective collaborators with the SFE technique in general, and willingness to participate as a volunteer.

Subsequent sections of this report present the conclusions and recommendations from this study, the experimental details, and the results. The proposed draft protocols for SFE (Method 3560) and for the total recoverable hydrocarbon determination by IR (Method 8440) are included as Appendices A and B, respectively. The list of instructions sent to the collaborating laboratories is included as Appendix C.

SECTION 2

CONCLUSIONS

The results of this collaborative study, conducted with 14 laboratories, indicate that the proposed EPA SW-846 Method 3560 (when used with EPA Method 8440) has an accuracy of 82.9 percent (this is the overall method recovery of petroleum hydrocarbons from samples containing from 614 to 32,600 mg/kg TPHs). The interlaboratory method precisions ranged from 17.3 to 45.4 percent relative standard deviation for PE-FTIR analyses and from 16.7 to 47.9 percent for BSci-IR analyses; the intralaboratory method precisions ranged from 11.5 to 17 percent for PE-FTIR analyses and from 11.1 to 18.2 percent for BSci-IR analyses.

The matrices used in this study were three standard reference soils and clay soil spiked with TPHs; the TPH recoveries for these samples, which were homogeneous and dry (with the exception of the clay soil that had a water content of 10.6 percent), were greater than 75 percent. However, when the clay soil was mixed with additional water to bring its water content to 30 percent, followed by extraction by SFE, less than 38 percent recovery was achieved. Results from additional experiments performed with anhydrous magnesium sulfate and diatomaceous earth as drying agents indicated that samples containing 20 percent water or more need to be mixed with anhydrous magnesium sulfate and allowed to equilibrate for several hours (preferably overnight in sealed containers and at 4°C to minimize losses of volatile petroleum hydrocarbons) before extraction by SFE.

Contamination of the SFE system is likely, especially when high-contamination samples are extracted; the analyst must take precautions to minimize cross-contamination of extracts.

The participating laboratories generally stayed within the guidelines of the instructions provided by MRI-CO; however, variations in the carbon dioxide flow rate, the extraction vessel dimensions, design, and orientation, the mode of collection of the extracted material, and the temperature of the collection solvent/trap were unavoidable because the various SFE systems were not identical in design. Therefore, it is not surprising that the variation from laboratory to laboratory was greater than that attributed to the analytical error displayed within laboratories. Nonetheless, the participating laboratories did a good job; when the data were subjected to a software program from the Association of Official Analytical Chemists (AOAC), data from only one laboratory were rejected on two of the matrices.

SECTION 3

RECOMMENDATIONS

On the basis of this interlaboratory method validation study, the Office of Solid Waste may want to consider Method 3560 for incorporation in the SW-846 methods manual.

The specific recommendations from the participating laboratories are listed below:

- The temperature of the collection solvent should be kept at 0 to 5°C to minimize losses of volatile organics.
- Laboratory 02 reported that all glass wool in their laboratory was contaminated with petroleum hydrocarbons. Therefore, they recommended that glass wool be washed and dried in a muffle furnace before use.
- Laboratory 02 recommended that glass wool be replaced with a drying agent (magnesium sulfate, diatomaceous earth) to "protect" the frits and fill the void volume of the extraction vessel.
- Laboratory 13 recommended that precleaned sand be used to fill the void volume of the extraction vessel.
- Laboratories 05 and 13 recommended that the direction of carbon dioxide flow be recorded. Furthermore, they stated that the carbon dioxide flow from top to bottom of the extraction vessel is superior to flow from bottom to top when the extraction cell is not full.
- The forms used by the laboratories in reporting the operating conditions (Appendix C) were very helpful, and we recommend that these forms be used to record the SFE conditions.

The method presented here is not suitable for the extraction of hydrocarbons from gasoline-contaminated soil samples because of poor collection efficiencies of the volatile hydrocarbons present in gasoline. Improvements in the collection method, such as trapping onto an adsorbent trap held at -10°C, followed by rinsing of the trap with an organic solvent (PCE) and analysis by gas chromatography with infrared detection, is recommended for future studies.

SECTION 4

EXPERIMENTAL

COLLABORATIVE STUDY

Study Design

The study was based on the AOAC's blind replicate design with balanced replicates for the collaborative evaluation of precision and accuracy of an analytical method (3,4). In our study, four soil samples were to be extracted in triplicate. Samples 1, 2, and 3 were standard reference materials that had been certified for TPHs by the modified Method 418.1(5); their TPH contents were 614, 2050, and 32,600 mg/kg, respectively. The fourth sample was a clay soil that we spiked with motor oil at 10,000 mg/kg. The samples, including their code numbers and TPH levels, are described in Table 1.

Each laboratory received 10-g portions from matrices 1, 2, and 3, and three 3-g portions of clay soil spiked with motor oil at 10,000 mg/kg. The spiked clay soil samples were identified as samples 4, 6, and 7. In addition, they received the unspiked clay soil, identified as sample 5, and clay soil samples spiked with corn oil/reference oil and motor oil (identified as samples 8 and 9 respectively). The corn oil was used to simulate lipids and similar materials; the reference oil is a mixture of n-hexadecane, isooctane, and chlorobenzene. To sample 9 we added water (after spiking with the motor oil) to bring its water content to 30 percent. A total of 15 extractions were to be performed by each laboratory. The laboratories were instructed to extract three 3-g portions from each of matrices 1, 2, and 3, and all of the material they received for samples 4 through 9.

Sets of test samples were sent to 17 laboratories; 14 laboratories submitted extracts within the time frame specified by MRI-CO. The laboratories were to perform the extractions according to the instructions given in Appendix C and then ship the extracts to MRI-CO for IR analysis. To minimize errors due to reagent contamination, each laboratory was given tetrachloroethylene (PCE) for use as collection solvent, and anhydrous sodium sulfate that was to be added to the water-containing spiked clay sample. In addition, we made arrangements to provide each laboratory with the SFE-grade carbon dioxide from Scott Specialty Gases. However, due to delays in the arrival of the carbon dioxide, laboratories 05, 12, 15, and 17 used their own SFE-grade carbon dioxide (Air Products).

Apparatus

The SFE systems used by the laboratories in this study are identified in Table 2. The actual operating conditions, as reported by the laboratories, are given in Tables 3 through 8. The flow rates of the carbon dioxide reported by the participating laboratories are given in Table 9. The exact

TABLE 1. SAMPLES USED IN THE INTERLABORATORY STUDY^a

| Sample no. | Sample identification | TPH concentration (mg/kg) | Source |
|------------|--|---------------------------|---|
| 1A, B, C | TPH-1 soil (Lot 91017) without fatty acids | 614 ^b | Environmental Research Associates, Arvada, CO |
| 2A, B, C | TPH-2 soil (Lot 91017) with fatty acids | 2,050 ^b | Environmental Research Associates, Arvada, CO |
| 3A, B, C | SRS103-100 soil | 32,600 ^c | Fisher Scientific, Pittsburgh, PA |
| 4 | Clay soil spiked with motor oil | 10,000 ^d | MRI-CO ^e |
| 5 | Clay soil (unspiked) | -- | MRI-CO |
| 6 | Clay soil spiked with motor oil | 10,000 ^d | MRI-CO |
| 7 | Clay soil spiked with motor oil | 10,000 ^d | MRI-CO |
| 8 | Clay soil spiked with corn oil and reference oil | 1,000 ^d | MRI-CO |
| 9 | Clay soil spiked with motor oil and water (approximately 30 percent by weight) | 10,000 ^d | MRI-CO |

^a Samples 1, 2, and 3 were extracted in triplicate by each laboratory. The three replicates were identified as A, B, and C. Only one extraction was performed for samples 4 through 9. Approximately 3 g of sample was extracted in each case. The exact amounts extracted by each laboratory are given in Table 10. The final volumes of the extracts before dilution are given in Table 11.

^b Certified value.

^c Determined by Soxhlet extraction with Freon-113 and analysis by IR spectrometry. Average of duplicate determinations.

^d Spike value.

^e Midwest Research Institute, California Operations.

TABLE 2. IDENTIFICATION OF THE SFE SYSTEMS USED IN THE COLLABORATIVE STUDY

| Laboratory code | SFE system identification |
|-----------------|-------------------------------|
| 01 | Suprex SFE-50 |
| 02 | Suprex Prepmaster |
| 03 | CCS Instruments SFE |
| 04 | Dionex-Lee Scientific SFE 703 |
| 05 | Isco SFE System 1200 |
| 06 | Suprex Prepmaster |
| 08 | Dionex-Lee Scientific SFE 703 |
| 10 | HP 7680A |
| 11 | HP 7680A |
| 12 | HP 7680A |
| 13 | Isco SFE System 1200 |
| 14 | Isco SFE System 1200 |
| 15 | Isco SFE System 1200 |
| 17 | Isco SFE System 1200 |

TABLE 3. SFE OPERATING CONDITIONS FOR THE ISCO SFE SYSTEMS

| Parameter | 05 | 13 | 14 | 15 | 17 |
|------------------------------------|---------------------|---------------------|---------------------|------------------------------|--------------------------|
| CO ₂ Pressure (atm) | 340 | 340 | 340 | 340 | 340 |
| CO ₂ Density (g/mL) | 0.785 | 0.785 | 0.785 | 0.785 | 0.785 |
| CO ₂ Flow rate (mL/min) | Varied | Varied | Varied | Varied | Varied |
| (see Table 9) | (see Table 9) | (see Table 9) | (see Table 9) | (see Table 9) | (see Table 9) |
| Oven temperature (°C) | 80 | 80 | 80 | 80 | 80 |
| Extraction time (min) | 30 | 30 | 30 | 30 | 30 |
| Extraction vessel volume (mL) | 10 | 2.5 | 2.5 | 10 | 2.5 and 10 ^a |
| Extraction vessel dimensions | 15-mm ID x | 7.5-mm ID x | 7.5-mm ID x | 15-mm ID x | 7.5/15-mm ID x |
| 56-mm length | 56-mm length | 56-mm length | 56-mm length | 56-mm length | |
| Extraction vessel orientation | Vertical, down flow | Vertical, down flow | Vertical, down flow | Vertical, down flow | Vertical, down flow |
| Restrictor dimensions | 50-μm ID x | 32-μm ID x | 50-μm ID x | 50-μm ID x | 50-μm ID x |
| 35-cm length | 11-cm length | 50-cm length | 37.5-cm length | 60/70-cm length ^b | |
| Restrictor temperature (°C) | Not heated | 80 | Not heated | Not heated | Not heated |
| Collection solvent | PCE | PCE | PCE | PCE | PCE |
| Volume of solvent (mL) | 3 | 3 | 3 | 3 | 3 |
| Temperature of collection | Ambient | Ambient | Ambient | Ambient | Ambient |
| | temperature | temperature | temperature | temperature | temperature ^c |

^a The 10-mL vessel was used for samples 3A, 3B, 3C, and 9.

^b A 60-cm length restrictor was used with the 10-mL extraction vessel; a 70-cm length restrictor was used with the 2.5-mL extraction vessel.

^c Initial temperature was room temperature; however, no attempt was made to control the collection solvent temperature during SFE.

TABLE 4. SFE OPERATING CONDITIONS FOR THE DIONEX-LEE SCIENTIFIC SFE-703 SYSTEMS

| Parameter | 04 | 08 |
|-------------------------------------|--------------------------|--------------------------|
| CO ₂ Pressure (atm) | 340 | 340 |
| CO ₂ Density (g/mL) | 0.785 | 0.785 |
| CO ₂ Flow rate (mL/min) | Varied (see Table 9) | Varied (see Table 9) |
| Oven temperature (°C) | 80 | 80 |
| Extraction time (min) | 30 | 30 |
| Extraction vessel volume (mL) | 3.5 | 3.5 |
| Extraction vessel dimensions | 9.4-mm ID x 50-mm length | 9.4-mm ID x 50-mm length |
| Extraction vessel orientation | Horizontal | Horizontal |
| Restrictor dimensions | a | a |
| Restrictor temperature (°C) | 150 | 150 |
| Collection solvent | PCE | PCE |
| Volume of solvent (mL) | 5 | 5 |
| Temperature of collection vial (°C) | 2 | 2 |

^a Restrictor identified as "250" restrictor.

TABLE 5. SFE OPERATING CONDITIONS FOR THE SUPREX PREPMASTER SFE SYSTEMS

| Parameter | 02 | 06 |
|-------------------------------------|----------------------------------|-------------------------------|
| CO ₂ Pressure (atm) | 340 | 340 |
| CO ₂ Density (g/mL) | 0.785 | 0.785 |
| CO ₂ Flow rate (mL/min) | 1.2 | Varied (see Table 9) |
| Oven temperature (°C) | 80 | 80 |
| Extraction time (min) | 30 | 30 |
| Extraction vessel volume (mL) | 5 | 3 |
| Extraction vessel dimensions | 9.5-mm ID x 65-mm length | 10-mm ID x 38-mm length |
| Extraction vessel orientation | Vertical, up flow | Vertical, up flow |
| Restrictor dimensions | a | 50- μ m ID x 10-cm length |
| Restrictor temperature (°C) | 100 | Not heated |
| Collection solvent | PCE | PCE |
| Volume of solvent (mL) | 3 | 3 |
| Temperature of collection vial (°C) | Ambient temperature ^b | Ambient temperature |

^a A prototype VariFlow restrictor (variable) was used.

^b Initial temperature was ambient temperature; however, no attempt was made to control the collection solvent temperature during SFE.

TABLE 6. SFE OPERATING CONDITIONS FOR THE SUPREX SFE-50 SYSTEM

| Parameter | 01 |
|-------------------------------------|----------------------------------|
| CO ₂ Pressure (atm) | 340 |
| CO ₂ Density (g/mL) | 0.785 |
| CO ₂ Flow rate (mL/min) | Varied (see Table 9) |
| Oven temperature (°C) | 80 |
| Extraction time (min) | 30 |
| Extraction vessel volume (mL) | 3.5 |
| Extraction vessel dimensions | 9.4-mm ID x 50-mm length |
| Extraction vessel orientation | Horizontal |
| Restrictor dimensions | 50- μ m ID x 60-cm length |
| Restrictor temperature (°C) | Not heated |
| Collection solvent | PCE |
| Volume of solvent (mL) | 3 |
| Temperature of collection vial (°C) | Ambient temperature ^a |

^a Initial temperature was ambient temperature; however, no attempt was made to control the collection solvent temperature during SFE.

TABLE 7. SFE OPERATING CONDITIONS FOR THE CCS INSTRUMENTS SFE SYSTEM

| Parameter | 03 |
|-------------------------------------|-------------------------------------|
| CO ₂ Pressure (atm) | 340 |
| CO ₂ Density (g/mL) | 0.785 |
| CO ₂ Flow rate (mL/min) | 0.2 |
| Oven temperature (°C) | 80 |
| Extraction time (min) | 30 |
| Extraction vessel volume (mL) | 6.0 |
| Extraction vessel dimensions | 9.5-mm ID x 127-mm length |
| Extraction vessel orientation | Vertical, up flow |
| Restrictor dimensions | 20- μ m ID x 5-cm length |
| Restrictor temperature (°C) | 80 to 100 |
| Collection solvent | PCE |
| Volume of solvent (mL) | 3 |
| Temperature of collection vial (°C) | Ambient temperature ^a |

^a Initial temperature was ambient temperature; however, no attempt was made to control the collection solvent temperature during SFE.

TABLE 8. SFE OPERATING CONDITIONS FOR THE HP 7680A SFE SYSTEMS

| Parameter | 10 | 11 | 12 |
|---|----------------------------|----------------------------|----------------------------|
| CO ₂ Pressure (atm) | 340 | 370 | 370 |
| CO ₂ Density (g/mL) | 0.78 | 0.80 | 0.80 |
| CO ₂ Flow rate (mL/min) | 2.0 | 2.0 | 2.0 |
| Oven temperature (°C) | 80 | 80 | 80 |
| Extraction time (min) | 30 | 30 | 20 |
| Extraction vessel volume (mL) | 7.0 | 7.0 | 7.0 |
| Extraction vessel dimensions | 10-mm ID x 90-mm length | 10-mm ID x 90-mm length | 10-mm ID x 90-mm length |
| Volumes swept (mL) | 10.2 | 9.9 | 6.6 |
| Extraction vessel orientation | Vertical, up flow | Vertical, up flow | Vertical, up flow |
| Nozzle temperature (°C) | 50 | 45 | 45 |
| Trap temperature (°C) | -15 | 10 | 10 |
| Trap packing material | Stainless- steel beads | Stainless- steel beads | ODS ^a |
| Rinse solvent | PCE | PCE | PCE |
| Volume | 1.2 | 1.5 | 1.5 |
| Rate (mL/min) | 1.0 | 1.0 | 1.0 |
| Nozzle temperature during rinse (°C) | 40 | 45 | 45 |
| Trap temperature during rinse (°C) | 40 | 60 | 40 |
| Number of rinses | 2 | 2 | 2 |

^a ODS - octadecylsilyl-bonded silica.

TABLE 9. FLOW RATE (mL/min) OF CARBON DIOXIDE (AS LIQUID) USED BY THE PARTICIPATING LABORATORIES

| Sample ID | 01 | 02 | 03 | 04 | 05 | 06 | 08 | 10 | 11 | 12 | 13 | 14 ^a | 15 | 17 |
|-----------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----------------|-----|-----|
| 1A | 1.2 | 1.2 | 0.2 | 0.7 | 1.9 | 0.9 | 0.8 | 2.0 | 2.0 | 2.0 | 0.8 | | 2.4 | 1.6 |
| 1B | 1.2 | 1.2 | 0.2 | 0.7 | 1.7 | 0.7 | 0.7 | 2.0 | 2.0 | 2.0 | 0.8 | | 2.3 | 1.6 |
| 1C | 1.1 | 1.2 | 0.2 | 0.6 | 2.0 | 0.5 | 0.6 | 2.0 | 2.0 | 2.0 | 0.9 | | 2.3 | 1.9 |
| 2A | 1.0 | 1.2 | 0.2 | 0.7 | 1.9 | 1.2 | 0.7 | 2.0 | 2.0 | 2.0 | 0.6 | | 2.2 | 1.6 |
| 2B | 1.2 | 1.2 | 0.2 | 0.6 | 1.6 | 1.2 | 0.7 | 2.0 | 2.0 | 2.0 | 0.5 | | 2.2 | 1.6 |
| 2C | 1.0 | 1.2 | 0.2 | 0.6 | 1.9 | 0.6 | 0.6 | 2.0 | 2.0 | 2.0 | 0.6 | | 2.0 | 1.2 |
| 3A | 0.6 | 1.2 | 0.2 | 0.5 | 1.5 | 2.2 | 0.7 | 2.0 | 2.0 | 2.0 | 0.6 | | 2.0 | 0.8 |
| 3B | 0.5 | 1.2 | 0.2 | 0.6 | 1.5 | 1.2 | 0.6 | 2.0 | 2.0 | 2.0 | 0.6 | | 2.0 | 0.8 |
| 3C | 0.9 | 1.2 | 0.2 | 0.6 | 1.7 | 0.9 | 0.6 | 2.0 | 2.0 | 2.0 | 0.6 | | 1.5 | 0.8 |
| 4 | 0.9 | 1.2 | 0.2 | 0.6 | 1.7 | 1.5 | 0.4 | 2.0 | 2.0 | 2.0 | 0.7 | | 2.1 | 1.2 |
| 5 | 0.8 | 1.2 | 0.2 | 0.7 | 1.7 | 1.5 | 0.3 | 2.0 | 2.0 | 2.0 | 0.6 | | 2.0 | 1.8 |
| 6 | 1.1 | 1.2 | 0.2 | 0.6 | 1.2 | 1.0 | 0.6 | 2.0 | 2.0 | 2.0 | 0.7 | | 2.2 | 1.8 |
| 7 | 1.0 | 1.2 | 0.2 | 0.5 | 1.5 | 2.2 | 0.8 | 2.0 | 2.0 | 2.0 | 0.9 | | 2.4 | 1.0 |
| 8 | 0.6 | 1.2 | 0.2 | 0.7 | 1.3 | 2.4 | 0.6 | 2.0 | 2.0 | 2.0 | 0.8 | | 2.1 | 1.0 |
| 9 | 0.9 | 1.2 | 0.2 | 0.6 | 1.4 | 2.4 | 0.3 | 2.0 | 2.0 | 2.0 | 0.8 | | 2.2 | 0.8 |

^a Reported as 1 to 2 mL/min but not measured for each extraction.

weights of the samples extracted by the individual laboratories are listed in Table 10, and the volumes of the extracts delivered to MRI-CO are given in Table 11.

Two infrared spectrometers were used: a Perkin-Elmer Corporation (Norwalk, CT) Model 1605 FTIR, interfaced with a Digital Equipment Corporation (Cupertino, CA) DEC 386SX computer (scanning from 3200 to 2700 cm^{-1}), and a Buck-Scientific IR filter spectrometer Model 404, centered on 2924 cm^{-1} with a band of 15 cm^{-1} on each side (E. Norwalk, CT).

Materials

The reference oil standards used in the IR determination of TPHs were prepared as follows: Pipet 15 mL n-hexadecane, 15 mL isooctane, and 10 mL chlorobenzene into a 50-mL glass-stoppered flask. Mix by swirling the contents. Pipet 0.5 mL of this mixture into a tared 100-mL volumetric flask, weigh to the nearest milligram, and dilute to volume with PCE. Pipet appropriate volumes of this stock standard solution into 100-mL volumetric flasks and dilute to volume with PCE to make 10- $\mu\text{g/mL}$, 25- $\mu\text{g/mL}$, 50- $\mu\text{g/mL}$, 100- $\mu\text{g/mL}$, 250- $\mu\text{g/mL}$, and 500- $\mu\text{g/mL}$ calibration standards. The motor oil calibration standards at 10 $\mu\text{g/mL}$, 25 $\mu\text{g/mL}$, 50 $\mu\text{g/mL}$, 100 $\mu\text{g/mL}$, 250 $\mu\text{g/mL}$, and 500 $\mu\text{g/mL}$ were prepared by serial dilution of a stock solution made at 200 $\mu\text{g}/\mu\text{L}$ in PCE. Corn oil at 40 mg/mL in Freon-113 was used to spike the clay soil.

The following materials were used in the study:

n-Hexadecane, isooctane, chlorobenzene (Aldrich Chemical, Milwaukee, WI).

Corn oil (local supermarket).

Motor oil (local Shell gas station).

PCE, glass-distilled, HPLC grade, lot no. 06626TX (Aldrich Chemical).

Silica gel, 70 to 230 mesh, ASTM (Baxter Scientific Products, McGaw Park, IL).

Anhydrous sodium sulfate, analytical reagent (Mallinckrodt, St. Louis, MO).

Anhydrous magnesium sulfate, analytical reagent (Mallinckrodt).

Hydromatrix (diatomaceous earth) (Isco, Inc., Lincoln, NB).

Clay soil (Sandoz Experimental Station, San Jose, CA) (33.6% sand, 35.4% silt, 31% clay; organic carbon 1.8%, water 10.6%).

Carbon dioxide, SFC/SFE-grade (Air Products, Allentown, PA).

Carbon dioxide, SFE-grade (Scott Specialty Gases, Inc., Plumsteadville, PA).

The standard reference materials used are listed in Table 1.

TABLE 10. WEIGHTS (g) OF THE SAMPLES EXTRACTED IN THE STUDY

| Laboratory code | | | | | | | | | | | | | | | |
|--------------------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|
| | 1A | 1B | 1C | 2A | 2B | 2C | 3A | 3B | 3C | 4 | 5 | 6 | 7 | 8 | 9 |
| 01 | 3.03 | 3.00 | 3.03 | 3.00 | 3.03 | 3.04 | 2.98 | 3.07 | 3.05 | 3.02 | 3.00 | 3.00 | 3.02 | 3.01 | 1.62 |
| 02 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 2.40 |
| 03 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | a | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 |
| 04 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.02 | 2.98 | 3.00 | 3.01 | 2.99 | 4.29 |
| 05 | 3.03 | 3.00 | 3.01 | 3.00 | 3.01 | 3.01 | 3.06 | 3.01 | 3.05 | 3.01 | 2.99 | 3.01 | 3.01 | 3.00 | 3.01 |
| 06 | 3.00 | 3.00 | 3.00 | 3.01 | 3.00 | 3.01 | 3.00 | 3.00 | 3.00 | 3.00 | 2.99 | 3.00 | 3.00 | 3.07 | 2.91 |
| 08 | 3.09 | 3.02 | 3.03 | 3.02 | 3.01 | 3.07 | 3.11 | 2.98 | 3.02 | 3.01 | 3.00 | 2.99 | 3.00 | 2.99 | 2.78 |
| 10 | 3.12 | 3.07 | 2.99 | 3.01 | 3.00 | 3.00 | 3.02 | 3.00 | 3.03 | 2.99 | 2.99 | 2.98 | 3.00 | 2.99 | 2.52 |
| 11 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 2.99 | 3.00 | 3.00 | 3.01 | 3.03 | 2.34 |
| 12 | 3.09 | 3.07 | 3.10 | 3.03 | 3.10 | 3.06 | 3.03 | 3.10 | b | 3.01 | 2.98 | 3.00 | 3.00 | 2.99 | 2.48 |
| 13 | 3.01 | 2.98 | 3.00 | 3.01 | 3.00 | 3.00 | 2.99 | 3.03 | 3.07 | 3.01 | 2.99 | 3.01 | 3.02 | 3.00 | 2.86 |
| 14 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 1.24 | 3.00 | 3.00 | 2.70 | 3.00 | 2.50 | 3.00 |
| 15 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 3.00 | 2.55 |
| 17 | 3.00 | 3.00 | 3.00 | 3.01 | 3.00 | 3.00 | 3.02 | 3.03 | 3.01 | 3.00 | 2.99 | 3.00 | 3.00 | 3.04 | 2.95 |

^a The extraction was not performed.

^b The extract was not submitted because the SFE system developed a leak.

TABLE 11. FINAL VOLUMES OF THE EXTRACTS (mL) BEFORE DILUTION

| Laboratory code | 1A | 1B | 1C | 2A | 2B | 2C | 3A | 3B | 3C | 4 | 5 | 6 | 7 | 8 | 9 |
|--------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 01 | 3.8 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |
| 02 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |
| 03 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | a | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.1 |
| 04 | 4.5 | 4.5 | 4.5 | 4.5 | 4.5 | 4.5 | 4.5 | 4.6 | 4.6 | 4.5 | 4.4 | 4.5 | 4.6 | 4.5 | 4.4 |
| 05 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.1 | 3.0 | 3.0 | 3.0 | 3.4 | 3.0 |
| 06 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |
| 08 | 4.6 | 4.7 | 4.7 | 4.7 | 4.7 | 4.7 | 4.8 | 4.7 | 4.7 | 4.6 | a | 4.6 | 4.6 | 4.6 | 4.8 |
| 10 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.9 | 3.8 | 3.0 | 3.0 | 3.0 | 3.2 | 3.1 | 3.0 | 3.0 |
| 11 | 3.4 | 3.4 | 3.4 | 3.3 | 3.3 | 3.4 | 3.6 | 3.6 | 3.6 | 3.5 | 3.0 | 3.5 | 3.5 | 3.3 | 3.4 |
| 12 | 3.2 | 3.2 | 3.2 | 3.2 | 3.2 | 3.2 | 3.3 | 3.4 | a | 3.3 | 3.3 | 3.3 | 3.4 | 3.3 | 3.4 |
| 13 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |
| 14 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |
| 15 | 3.0 | 3.0 | 3.0 | 3.3 | 3.0 | 3.2 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |
| 17 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |

^a The extract was not submitted for analysis.

Sample Spiking Procedure

For spiking of the clay soil samples with motor oil, 3.0-g portions of the soil were weighed into 7-mL glass vials, and portions (150 μL) of concentrated stock solution containing motor oil in Freon-113 at 200 $\mu\text{g}/\mu\text{L}$ were added with a syringe while ensuring that the solution did not contact the walls of the vials. Mixing was performed with the tip of a disposable glass pipet. After the solvent had completely evaporated (approximately 15 min), the vials with the spiked samples were sealed with Teflon-lined caps and stored at 4°C for 5 days prior to shipping to the laboratories. Spiking with corn oil and reference oil was done as described above for the motor oil, except that the concentrations of the stock solutions were 20 $\mu\text{g}/\mu\text{L}$. The clay soil samples that required adjustment of their water content to 30 percent were also spiked individually; a 2.4-g portion of the clay soil sample was weighed into a glass vial, spiked, and then 0.6 mL water was added.

Sample Extraction

The methods used to extract the samples and to analyze the extracts are included in Appendices A and B to this report. Each laboratory used only one analyst to perform sample extractions and allowed one block of time for all sample extractions.

TREATMENT OF DATA

Outlier Testing

Outlier testing was done using both the Cochran test and the Grubbs test (3, 4). The former is used for the removal of results from laboratories that show significantly greater variability among replicate (intralaboratory) analyses than the other laboratories for a given material. To calculate the Cochran test statistic, the intralaboratory variance for each laboratory was computed, and the largest of these values was divided by the sum of all these variances. When a laboratory was rejected on the basis of these tests, its results were removed from the set of data for a particular matrix, and the test was repeated using the remaining data in the subset.

Grubbs tests were performed to remove results from laboratories with extreme averages. The single-value test (2-tail; $P = 0.01$) was run first; then, when no outlier was found, the pair value test (two values at the highest end, two values at the lowest end, and two values, one at each end, at an overall $P = 0.01$) was applied.

To calculate the single Grubbs test statistic, the average for each of the L laboratories was computed, and the standard deviation of these L averages (designated as the original s) was calculated. The standard deviation of the set of averages with the highest average removed (s_H) was calculated, then the standard deviation of the set of averages with the lowest average removed (s_L) was calculated. Finally, the percentage decrease in standard deviation was calculated as follows:

$$100 \times [1 - (s_L/s)] \text{ and } 100 \times [1 - (s_H/s)]$$

The higher of these two percentage decreases was the single Grubbs statistic, which signals the presence of an outlier to be omitted if it *exceeds* the critical value listed in the single Grubbs tables, at the $P = 0.01$ level, 2-tail, for L laboratories.

To calculate the Grubbs pair statistic, we proceeded in an analogous fashion, except that we calculated the standard deviations, from the original set of averages. The smallest of these three standard deviation values was taken and the corresponding percentage in standard deviation from the original s was calculated. A Grubbs outlier pair was present if the selected value for the percentage decrease from the original s *exceeded* the critical value listed in the Grubbs pair value table, at the $P = 0.01$ level, for L laboratories.

Statistical Analysis

Statistical analysis was performed using the AOAC Lotus spreadsheet program developed for the analysis of data from interlaboratory studies (4).

The summary statistics that are reported for each matrix include the mean concentration of TPHs and the method precision (percent relative standard deviation). The repeatability relative standard deviation (RSD_r), which was determined from the repeatability standard deviation (s_r) and the mean concentration of a particular matrix, is an indication of the intralaboratory precision. The reproducibility relative standard deviation (RSD_R), which was determined from the reproducibility standard deviation (s_R) and the mean concentration of a particular matrix, is an indication of the interlaboratory method precision.

QUALITY ASSURANCE/QUALITY CONTROL

The extractions of the study samples were performed according to the instructions provided by MRI-CO to all collaborators (Appendix C). To minimize errors due to reagent contamination, we provided each laboratory with the collection solvent, the extractant (carbon dioxide), and the drying agent (anhydrous sodium sulfate). We made arrangements with Scott Specialty Gases to supply carbon dioxide from the same lot to all participating laboratories. All but four laboratories used SFE-grade carbon dioxide from Scott Specialty Gases to extract the test samples. Laboratories 05, 12, 15, and 17 experienced delay in delivery of the carbon dioxide and therefore used SFC/SFE-grade carbon dioxide from Air Products. Each laboratory was instructed to report the SFE operating conditions on special forms (provided by MRI-CO), and to record the exact mass of the study sample extracted by SFE.

The IR analyses were performed on two different IR systems according to the proposed Method 8440 included in Appendix B. The following quality control procedures were implemented:

- A six-level calibration (at 10, 25, 50, 100, 250, and 500 $\mu\text{g/mL}$) was performed every day, on each instrument, during the period the IR analyses were performed. The multilevel calibration was then verified after every 10 analyses by analyzing a reference oil standard or a motor oil standard at 100 $\mu\text{g/mL}$. The reference oil standard was used to quantify TPHs in samples 1, 2, 3, 5, and 8. The motor oil standard was used to quantify TPHs in samples 4, 6, 7, and 9. Corn oil was quantified against the reference oil standard.
- PCE blanks were analyzed daily on each instrument during the period the IR analyses were performed.
- All system blanks received from the participating laboratories were analyzed for TPHs.
- Sample 5 was a blind QC sample (unspiked clay soil).

At the MRI-CO, six clay soil samples were spiked with motor oil at 10,000 mg/kg and stored at 4°C in the dark; two were extracted with supercritical carbon dioxide after 22 days, two after 34 days, and two after 40 days of storage. These six samples were part of the batch of samples prepared for the collaborative study and were expressly set aside for extraction at a later time but still within the time frame in which the collaborating laboratories would carry out their extractions (six laboratories submitted extracts within 30 days of sample spiking, five within 40 days, and three within 60 days).

SECTION 5

RESULTS AND DISCUSSION

INTERLABORATORY METHOD PERFORMANCE

The results of the interlaboratory study have been summarized by matrix and by the IR spectrometer used in the analysis. As mentioned in Section 1, all IR determinations were performed at MRI-CO on two different IR spectrometers. Tables 12 through 19 present the concentration data for each of the four matrices that were extracted in triplicate (TPH-1 soil, TPH-2 soil, SRS103-100 soil, and the clay soil spiked with motor oil). The mean concentrations, the repeatability standard deviations (s_r), the reproducibility standard deviations (s_R), the repeatability relative standard deviations (RSD_r), the reproducibility relative standard deviations (RSD_R) and the mean recoveries have been summarized in Tables 20 through 25.

Rejection of Outliers

For the entire study, the AOAC software program rejected only one laboratory (laboratory 03) on two matrices. This laboratory achieved very low recoveries on all test samples because the carbon dioxide flow rate was too low (0.2 mL/min); however, the outlier test only rejected values on the TPH-2 soil extract and on the extract from the clay soil spiked with motor oil. Table 26 summarizes the outlier testing results. When data from all laboratories were pooled, and the outlier testing was performed by matrix, then laboratory 03 data for TPH-2 and spiked clay soil samples were rejected on the basis of the single Grubbs test (lowest average). When data from five laboratories using Isco systems were pooled, and outlier testing was performed by matrix, then laboratory 14 data were rejected on the basis of the Cochran test (maximum intralaboratory variance). When data from laboratories using extraction vessels of 3.5 mL or less in volume were pooled, and outlier testing was performed by matrix, then data from laboratory 17 were rejected on the basis of the Cochran test (maximum intralaboratory variance).

Method Recovery

The summary statistics for the recoveries from each of the four matrices that were extracted by all laboratories, calculated after outlier removal, are presented in Table 21. The mean recoveries of TPHs from these four matrices ranged from 77.9 to 107 percent for analyses performed with the PE-FTIR spectrometer, and from 75.9 to 101 percent for analyses performed with the BSci-IR spectrometer (a discussion of the correlation between the data generated with the two IR spectrometers is given later in this section). These recoveries are in agreement with data that we reported previously for the single-laboratory evaluation of this method (1, 2); the percent differences between the mean recoveries for the interlaboratory study and the single-laboratory study were ranging from 1.8 to 28.7 percent.

TABLE 12. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN TPH-1 SOIL EXTRACTS USING THE PE-FTIR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|-----------------|----------------|-------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 618 | 808 | 1,000 | 809 | 23.6 |
| 02 | 449 | 453 | 487 | 463 | 4.5 |
| 03 | 67 | ND ^b | 77 | 51.3 | 70.5 |
| 04 | 352 | 391 | 420 | 388 | 8.8 |
| 05 | 943 | 859 | 854 | 885 | 5.7 |
| 06 | 250 | 478 | 44 | 257 | 84.5 |
| 08 | 737 | 799 | 933 | 823 | 12.2 |
| 10 | 693 | 1,070 | 789 | 851 | 23.0 |
| 11 | 972 | 1,210 | 1,100 | 1,090 | 10.9 |
| 12 | 691 | 701 | 656 | 683 | 23.6 |
| 13 | 673 | 604 | 620 | 632 | 5.7 |
| 14 | 639 | 605 | 580 | 608 | 4.9 |
| 15 | 883 | 963 | 885 | 910 | 5.0 |
| 17 | 575 | 712 | 832 | 706 | 18.2 |

^a The certified concentration of TPHs in this sample is 614 mg/kg.

^b ND - not detected; the estimated detection limit is 10 mg/kg.

TABLE 13. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN TPH-1 SOIL EXTRACTS USING THE BSci-IR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|-----------------|----------------|-------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 634 | 793 | 987 | 805 | 22.0 |
| 02 | 402 | 393 | 433 | 409 | 5.1 |
| 03 | 62.0 | ND ^b | 80.5 | 50.8 | 72.0 |
| 04 | 339 | 386 | 404 | 376 | 8.9 |
| 05 | 922 | 831 | 825 | 859 | 6.3 |
| 06 | 230 | 435 | 52.5 | 239 | 80.0 |
| 08 | 751 | 821 | 940 | 837 | 11.4 |
| 10 | 590 | 1,060 | 690 | 780 | 31.7 |
| 11 | 959 | 1,190 | 1,080 | 1,080 | 10.7 |
| 12 | 621 | 625 | 593 | 613 | 2.8 |
| 13 | 594 | 535 | 544 | 558 | 5.7 |
| 14 | 563 | 532 | 510 | 535 | 5.0 |
| 15 | 859 | 944 | 856 | 886 | 5.6 |
| 17 | 497 | 639 | 747 | 628 | 20.0 |

^a The certified concentration of TPHs in this sample is 614 mg/kg.

^b ND - not detected; the estimated detection limit is 10 mg/kg.

TABLE 14. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN TPH-2 SOIL EXTRACTS USING THE PE-FTIR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|----------------|----------------|-------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 1,840 | 2,160 | 1,880 | 1,960 | 8.9 |
| 02 | 1,940 | 1,990 | 1,730 | 1,890 | 7.3 |
| 03 | 257 | 410 | 356 | 341 | 22.8 |
| 04 | 1,350 | 1,160 | 1,520 | 1,340 | 13.4 |
| 05 | 1,490 | 2,000 | 2,230 | 1,910 | 19.8 |
| 06 | 1,850 | 1,730 | 1,380 | 1,650 | 14.8 |
| 08 | 1,480 | 1,100 | 1,500 | 1,360 | 16.7 |
| 10 | 1,880 | 2,240 | 1,980 | 2,030 | 9.2 |
| 11 | 2,140 | 2,200 | 2,050 | 2,130 | 3.5 |
| 12 | 2,320 | 2,440 | 2,110 | 2,290 | 7.3 |
| 13 | 1,900 | 2,010 | 1,940 | 1,950 | 2.9 |
| 14 | 2,130 | 1,900 | 1,340 | 1,790 | 22.7 |
| 15 | 1,850 | 1,900 | 1,990 | 1,910 | 3.7 |
| 17 | 1,780 | 1,920 | 1,980 | 1,890 | 5.4 |

^a The certified concentration of TPHs in this sample is 2,050 mg/kg.

TABLE 15. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN TPH-2 SOIL EXTRACTS USING THE BSci-IR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|----------------|----------------|-------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 1,620 | 1,950 | 1,660 | 1,740 | 10.4 |
| 02 | 1,760 | 1,790 | 1,540 | 1,700 | 8.0 |
| 03 | 250 | 374 | 380 | 335 | 21.9 |
| 04 | 1,290 | 1,120 | 1,460 | 1,290 | 13.2 |
| 05 | 1,370 | 1,740 | 1,980 | 1,700 | 18.1 |
| 06 | 1,670 | 1,550 | 1,240 | 1,490 | 14.9 |
| 08 | 1,350 | 1,050 | 1,330 | 1,240 | 13.5 |
| 10 | 1,660 | 2,010 | 1,770 | 1,810 | 9.9 |
| 11 | 1,960 | 2,040 | 1,910 | 1,970 | 3.3 |
| 12 | 2,110 | 2,210 | 1,780 | 2,030 | 11.1 |
| 13 | 1,670 | 1,760 | 1,710 | 1,710 | 2.6 |
| 14 | 1,900 | 1,680 | 1,180 | 1,590 | 23.2 |
| 15 | 1,700 | 1,730 | 1,820 | 1,750 | 3.6 |
| 17 | 1,530 | 1,680 | 1,730 | 1,650 | 6.3 |

^a The certified concentration of TPHs in this sample is 2,050 mg/kg.

TABLE 16. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN SRS103-100 SOIL EXTRACTS USING THE PE-FTIR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|----------------|----------------|--------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 32,600 | 35,600 | 33,700 | 34,000 | 4.5 |
| 02 | 28,900 | 29,600 | 29,500 | 29,300 | 12.9 |
| 03 | 4,130 | 11,100 | b | 7,620 | 64.7 |
| 04 | 23,200 | 20,800 | 25,000 | 23,000 | 9.2 |
| 05 | 25,200 | 29,500 | 10,800 | 21,800 | 44.9 |
| 06 | 46,400 | 45,100 | 54,600 | 48,700 | 10.6 |
| 08 | 25,500 | 24,400 | 24,400 | 24,800 | 2.6 |
| 10 | 19,500 | 14,800 | 17,500 | 17,300 | 13.6 |
| 11 | 28,900 | 27,800 | 28,300 | 28,300 | 1.9 |
| 12 | 33,200 | 32,700 | b | 32,900 | 1.1 |
| 13 | 30,900 | 30,900 | 31,500 | 31,100 | 1.1 |
| 14 | 26,500 | 29,200 | 11,400 | 22,400 | 42.8 |
| 15 | 31,100 | 32,200 | 30,800 | 31,400 | 2.3 |
| 17 | 21,600 | 24,300 | 22,600 | 22,800 | 6.0 |

^a The concentration of TPHs measured in our laboratory for this sample is 32,600 mg/kg.

^b This laboratory did not submit an extract for this sample.

TABLE 17. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN SRS103-100 SOIL EXTRACTS USING THE BSci-IR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|----------------|----------------|--------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 29,600 | 32,200 | 30,500 | 30,800 | 4.3 |
| 02 | 27,000 | 27,800 | 27,900 | 27,600 | 1.8 |
| 03 | 3,890 | 10,200 | b | 7,050 | 63.3 |
| 04 | 22,700 | 20,600 | 24,400 | 22,600 | 8.4 |
| 05 | 22,400 | 25,800 | 9,790 | 19,300 | 43.7 |
| 06 | 42,400 | 41,100 | 49,100 | 44,200 | 9.7 |
| 08 | 23,600 | 23,200 | 23,800 | 23,500 | 1.3 |
| 10 | 18,300 | 14,000 | 16,200 | 16,200 | 13.3 |
| 11 | 26,600 | 25,100 | 25,700 | 25,800 | 2.9 |
| 12 | 30,500 | 30,300 | b | 30,400 | 0.5 |
| 13 | 27,900 | 27,600 | 28,600 | 28,000 | 1.8 |
| 14 | 24,500 | 26,800 | 11,300 | 20,900 | 40.0 |
| 15 | 28,300 | 29,400 | 28,200 | 28,600 | 2.3 |
| 17 | 20,600 | 22,800 | 21,300 | 21,600 | 5.2 |

^a The concentration of TPHs measured in our laboratory for this sample is 32,600 mg/kg.

^b This laboratory did not submit an extract for this sample.

TABLE 18. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN SPIKED CLAY SOIL EXTRACTS USING THE PE-FTIR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|----------------|----------------|-------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 8,760 | 10,400 | 8,680 | 9,280 | 10.5 |
| 02 | 7,850 | 8,740 | 9,740 | 8,780 | 10.8 |
| 03 | 580 | 368 | 963 | 637 | 47.3 |
| 04 | 7,000 | 7,590 | 7,270 | 7,290 | 4.1 |
| 05 | 8,160 | 8,790 | 9,290 | 8,750 | 6.5 |
| 06 | 6,240 | 6,810 | 4,100 | 5,720 | 25.0 |
| 08 | 7,110 | 6,130 | 5,780 | 6,340 | 10.9 |
| 10 | 6,470 | 4,310 | 3,200 | 4,660 | 35.7 |
| 11 | 9,690 | 8,380 | 9,400 | 9,160 | 7.5 |
| 12 | 8,660 | 8,810 | 9,270 | 8,910 | 3.6 |
| 13 | 8,420 | 8,210 | 8,830 | 8,490 | 3.7 |
| 14 | 5,000 | 7,910 | 8,970 | 7,290 | 28.2 |
| 15 | 8,560 | 7,660 | 9,020 | 8,410 | 8.2 |
| 17 | 7,990 | 7,940 | 8,590 | 8,170 | 4.4 |

^a The clay soil samples were spiked with motor oil at 10,000 mg/kg.

TABLE 19. CONCENTRATIONS (mg/kg) OF TPHs DETERMINED IN SPIKED CLAY SOIL EXTRACTS USING THE BSci-IR SPECTROMETER^a

| Laboratory code | Concentration (mg/kg) | | | | Percent RSD |
|--------------------|-----------------------|----------------|----------------|-------|----------------|
| | X ₁ | X ₂ | X ₃ | Mean | |
| 01 | 9,220 | 10,800 | 8,870 | 9,630 | 10.7 |
| 02 | 8,020 | 9,380 | 9,980 | 9,130 | 11.0 |
| 03 | 586 | 411 | 954 | 650 | 42.6 |
| 04 | 7,180 | 7,670 | 7,270 | 7,370 | 3.5 |
| 05 | 8,440 | 9,120 | 9,480 | 9,010 | 5.9 |
| 06 | 6,040 | 6,610 | 6,760 | 6,470 | 5.9 |
| 08 | 7,460 | 6,650 | 6,150 | 6,750 | 9.8 |
| 10 | 7,310 | 5,230 | 4,130 | 5,560 | 29.0 |
| 11 | 9,880 | 9,000 | 9,770 | 9,550 | 5.0 |
| 12 | 9,110 | 9,210 | 9,310 | 9,210 | 1.1 |
| 13 | 8,720 | 8,670 | 9,020 | 8,800 | 2.2 |
| 14 | 5,400 | 8,360 | 9,220 | 7,660 | 26.2 |
| 15 | 8,770 | 8,120 | 9,680 | 8,860 | 8.8 |
| 17 | 7,970 | 8,120 | 8,920 | 8,340 | 6.1 |

^a The clay soil samples were spiked with motor oil at 10,000 mg/kg.

TABLE 20. SUMMARY OF RESULTS FROM ALL LABORATORIES (BEFORE OUTLIER REMOVAL)

| Matrix | Certified or spike value (mg/kg) | Mean concentration ^a (mg/kg) | s _r (mg/kg) | s _R (mg/kg) | Percent RSD _r | Percent RSD _R | Percent mean recovery | Number of laboratories in the study |
|------------------------------------|--|---|---------------------------|---------------------------|-----------------------------|-----------------------------|-----------------------------|---|
| PE-FTIR | | | | | | | | |
| ERA TPH-1 soil | 614 | 654 | 111 | 297 | 17.0 | 45.4 | 107 | 14 |
| ERA TPH-2 soil | 2,050 | 1,750 | 206 | 509 | 11.8 | 29.2 | 85.4 | 14 |
| SRS103-100 soil | 32,600 | 26,820 | 4,320 | 9,720 | 16.1 | 36.2 | 82.3 | 14 |
| Clay soil spiked with motor oil | 10,000 | 7,280 | 968 | 2,490 | 13.3 | 34.2 | 72.8 | 14 |
| BSci-IR | | | | | | | | |
| ERA TPH-1 soil | 614 | 618 | 113 | 296 | 18.2 | 47.9 | 101 | 14 |
| ERA TPH-2 soil | 2,050 | 1,570 | 188 | 446 | 12.0 | 28.4 | 76.7 | 14 |
| SRS103-100 soil | 32,600 | 24,750 | 3,740 | 8,650 | 15.1 | 35.0 | 75.9 | 14 |
| Clay soil spiked with motor oil | 10,000 | 7,640 | 88 | 2,470 | 11.5 | 32.4 | 76.4 | 14 |

^a The number of replicates per laboratory was three.

TABLE 21. SUMMARY OF RESULTS FROM ALL LABORATORIES (AFTER OUTLIER REMOVAL)

| Matrix | Certified or spike value (mg/kg) | Mean concentration ^a (mg/kg) | s _r (mg/kg) | s _R (mg/kg) | Percent RSD _r | Percent RSD _R | Percent mean recovery | Number of laboratories retained |
|------------------------------------|--|---|---------------------------|---------------------------|-----------------------------|-----------------------------|-----------------------------|---------------------------------------|
| PE-FTIR | | | | | | | | |
| ERA TPH-1 soil | 614 | 654 | 111 | 297 | 17.0 | 45.4 | 107 | 14 |
| ERA TPH-2 soil | 2,050 | 1,850 | 213 | 321 | 11.5 | 17.3 | 90.2 | 13 |
| SRS103-100 soil | 32,600 | 26,820 | 4,320 | 9,720 | 16.1 | 36.2 | 82.3 | 14 |
| Clay soil spiked with motor oil | 10,000 | 7,790 | 1,000 | 1,660 | 12.9 | 21.3 | 77.9 | 13 |
| BSci-IR | | | | | | | | |
| ERA TPH-1 soil | 614 | 618 | 113 | 296 | 18.2 | 47.9 | 101 | 14 |
| ERA TPH-2 soil | 2,050 | 1,670 | 194 | 278 | 11.7 | 16.7 | 81.5 | 13 |
| SRS103-100 soil | 32,600 | 24,750 | 3,740 | 8,650 | 15.1 | 35.0 | 75.9 | 14 |
| Clay soil spiked with motor oil | 10,000 | 8,180 | 910 | 1,500 | 11.1 | 18.3 | 81.8 | 13 |

^a The number of replicates per laboratory was three.

TABLE 22. SUMMARY OF RESULTS FROM LABORATORIES USING ISCO SYSTEMS (BEFORE OUTLIER REMOVAL)

| Matrix | Certified or spike value (mg/kg) | Mean concentration ^a (mg/kg) | S _r (mg/kg) | S _R (mg/kg) | Percent RSD _r | Percent RSD _R | Percent mean recovery | Number of laboratories retained |
|---------------------------------|----------------------------------|---|------------------------|------------------------|--------------------------|--------------------------|-----------------------|---------------------------------|
| PE-FTIR | | | | | | | | |
| ERA TPH-1 soil | 614 | 748 | 68.3 | 152 | 9.1 | 20.3 | 122 | 5 |
| ERA TPH-2 soil | 2,050 | 1,890 | 256 | 256 | 13.5 | 13.5 | 92.2 | 5 |
| SRS103-100 soil | 32,600 | 25,900 | 6,170 | 7,020 | 23.8 | 27.1 | 79.4 | 5 |
| Clay soil spiked with motor oil | 10,000 | 8,220 | 1,030 | 1,030 | 12.5 | 12.5 | 82.2 | 5 |
| BSci-IR | | | | | | | | |
| ERA TPH-1 soil | 614 | 693 | 67.7 | 177 | 9.8 | 25.5 | 113 | 5 |
| ERA TPH-2 soil | 2,050 | 1,680 | 222 | 222 | 13.3 | 13.3 | 82.0 | 5 |
| SRS103-100 soil | 32,600 | 23,700 | 5,350 | 6,150 | 22.6 | 26.0 | 72.7 | 5 |
| Clay soil spiked with motor oil | 10,000 | 8,530 | 1,020 | 1,020 | 12.0 | 12.0 | 85.3 | 5 |

^a The number of replicates per laboratory was three.

TABLE 23. SUMMARY OF RESULTS FROM LABORATORIES USING ISCO SYSTEMS (AFTER OUTLIER REMOVAL)

| Matrix | Certified or spike value (mg/kg) | Mean concentration ^a (mg/kg) | S _r (mg/kg) | S _R (mg/kg) | Percent RSD _r | Percent RSD _R | Percent mean recovery | Number of laboratories retained |
|------------------------------------|--|---|---------------------------|---------------------------|-----------------------------|-----------------------------|-----------------------------|---------------------------------------|
| PE-FTIR | | | | | | | | |
| ERA TPH-1 soil | 614 | 748 | 68.3 | 152 | 9.1 | 20.3 | 122 | 5 |
| ERA TPH-2 soil | 2,050 | 1,890 | 256 | 256 | 13.5 | 13.5 | 92.2 | 5 |
| SRS103-100 soil | 32,600 | 25,900 | 6,170 | 7,020 | 23.8 | 27.1 | 79.4 | 5 |
| Clay soil spiked with motor oil | 10,000 | 8,460 | 507 | 507 | 6.0 | 6.0 | 84.6 | 4 |
| BSci-IR | | | | | | | | |
| ERA TPH-1 soil | 614 | 693 | 67.7 | 177 | 9.8 | 25.5 | 113 | 5 |
| ERA TPH-2 soil | 2,050 | 1,680 | 222 | 222 | 13.3 | 13.3 | 82.0 | 5 |
| SRS103-100 soil | 32,600 | 23,700 | 5,350 | 6,150 | 22.6 | 26.0 | 72.7 | 5 |
| Clay soil spiked with motor oil | 10,000 | 8,530 | 1,020 | 1,020 | 12.0 | 12.0 | 85.3 | 5 |

^a The number of replicates per laboratory was three.

TABLE 24. SUMMARY OF RESULTS FROM LABORATORIES USING EXTRACTION VESSELS 3.5 mL OR LESS IN VOLUME (BEFORE OUTLIER REMOVAL)

| Matrix | Certified or spike value (mg/kg) | Mean concentration ^a (mg/kg) | s _r (mg/kg) | s _R (mg/kg) | Percent RSD _r | Percent RSD _R | Percent mean recovery | Number of laboratories retained |
|------------------------------------|--|---|---------------------------|---------------------------|-----------------------------|-----------------------------|-----------------------------|---------------------------------------|
| PE-FTIR | | | | | | | | |
| ERA TPH-1 soil | 614 | 603 | 127 | 236 | 21.1 | 39.1 | 98.2 | 7 |
| ERA TPH-2 soil | 2,050 | 1,710 | 224 | 322 | 13.1 | 18.9 | 83.4 | 7 |
| SRS103-100 soil | 32,600 | 30,650 | 4,580 | 10,660 | 14.9 | 34.8 | 94.0 | 6 |
| Clay soil spiked with motor oil | 10,000 | 7,510 | 1,070 | 1,520 | 14.2 | 20.2 | 75.1 | 7 |
| BSci-IR | | | | | | | | |
| ERA TPH-1 soil | 614 | 568 | 117 | 236 | 20.6 | 41.4 | 92.5 | 7 |
| ERA TPH-2 soil | 2,050 | 1,530 | 203 | 259 | 13.3 | 16.9 | 74.6 | 7 |
| SRS103-100 soil | 32,600 | 28,300 | 3,960 | 9,190 | 14.0 | 32.4 | 86.8 | 6 |
| Clay soil spiked with motor oil | 10,000 | 7,860 | 927 | 1,360 | 11.8 | 17.3 | 78.6 | 7 |

^a The number of replicates per laboratory was three.

TABLE 25. SUMMARY OF RESULTS FROM LABORATORIES USING EXTRACTION VESSELS 3.5 mL OR LESS IN VOLUME (AFTER OUTLIER REMOVAL)

| Matrix | Certified or spike value (mg/kg) | Mean concentration ^a (mg/kg) | S _r (mg/kg) | S _R (mg/kg) | Percent RSD _r | Percent RSD _R | Percent mean recovery | Number of laboratories retained |
|------------------------------------|--|---|---------------------------|---------------------------|-----------------------------|-----------------------------|-----------------------------|---------------------------------------|
| PE-FTIR | | | | | | | | |
| ERA TPH-1 soil | 614 | 603 | 127 | 236 | 21.1 | 39.1 | 98.2 | 7 |
| ERA TPH-2 soil | 2,050 | 1,710 | 224 | 322 | 13.1 | 18.9 | 83.4 | 7 |
| SRS103-100 soil | 32,600 | 32,300 | 2,600 | 10,420 | 8.1 | 32.3 | 99.1 | 5 |
| Clay soil spiked with motor oil | 10,000 | 7,510 | 1,070 | 1,520 | 14.2 | 20.2 | 75.1 | 7 |
| BSci-IR | | | | | | | | |
| ERA TPH-1 soil | 614 | 568 | 117 | 236 | 20.6 | 41.4 | 92.5 | 7 |
| ERA TPH-2 soil | 2,050 | 1,530 | 203 | 259 | 13.3 | 16.9 | 74.6 | 7 |
| SRS103-100 soil | 32,600 | 29,820 | 2,200 | 8,890 | 7.4 | 29.8 | 91.5 | 5 |
| Clay soil spiked with motor oil | 10,000 | 7,860 | 927 | 1,360 | 11.8 | 17.3 | 78.6 | 7 |

^a The number of replicates per laboratory was three.

TABLE 26. OUTLIER TESTING—THE COCHRAN STATISTIC, THE SINGLE GRUBBS STATISTIC, AND THE DOUBLE GRUBBS STATISTIC^a

| Instrument ^b | Matrix | Number of laboratories | Cochran test | Single Grubbs test | Double Grubbs test |
|---|---------------------------------|------------------------|---------------------|--------------------|--------------------|
| All laboratories | | | | | |
| 1 | TPH-1 soil | 14 | 0.2725 | 17.86 | 29.75 |
| 2 | TPH-1 soil | 14 | 0.3448 | 15.27 | 25.12 |
| 1 | TPH-2 soil | 14 | 0.2781 | 43.77 ^c | 51.67 ^c |
| | | 13 | 0.2810 | 14.06 | 38.22 |
| 2 | TPH-2 soil | 14 | 0.2739 | 45.42 ^c | 52.68 ^c |
| | | 13 | 0.2769 | 13.30 | 31.38 |
| 1 | SRS103-100 soil | 14 | 0.3768 | 22.58 | 44.00 |
| 2 | SRS103-100 soil | 14 | 0.3699 | 22.01 | 45.07 |
| 1 | Clay soil spiked with motor oil | 14 | 0.3219 | 38.91 ^c | 51.60 ^c |
| | | 13 | 0.3242 | 20.78 | 36.35 |
| 2 | Clay soil spiked with motor oil | 14 | 0.3703 | 45.13 ^c | 54.43 ^c |
| | | 13 | 0.3729 | 16.96 | 28.01 |
| Laboratories using Isco SFE system | | | | | |
| 1 | TPH-1 soil | 5 | 0.7097 | 11.29 | 63.77 |
| 2 | TPH-1 soil | 5 | 0.6868 | 11.62 | 71.21 |
| 1 | TPH-2 soil | 5 | 0.5046 | 59.63 | 83.00 |
| 2 | TPH-2 soil | 5 | 0.5504 | 32.38 | 57.00 |
| 1 | STS103-100 soil | 5 | 0.5036 | 9.95 | 89.75 |
| 2 | SRS103-100 soil | 5 | 0.4973 | 11.24 | 73.54 |
| 1 | Clay soil spiked with motor oil | 5 | 0.8041 ^c | 57.72 | 70.64 |
| | | 4 | 0.4649 | 30.57 | 78.03 |
| 2 | Clay soil spiked with motor oil | 5 | 0.7715 | 47.04 | 80.15 |
| Laboratories using extraction vessels 3.5 mL or less in volume | | | | | |
| 1 | TPH-1 soil | 7 | 0.4151 | 24.18 | 53.41 |
| 2 | TPH-1 soil | 7 | 0.3833 | 19.07 | 34.84 |
| 1 | TPH-2 soil | 7 | 0.4690 | 12.88 | 51.37 |
| 2 | TPH-2 soil | 7 | 0.4732 | 15.40 | 48.30 |
| 1 | SRS103-100 soil | 6 | 0.7314 ^c | 48.02 | 60.03 |
| | | 5 | 0.7850 | 49.23 | 58.26 |
| 2 | SRS103-100 soil | 6 | 0.7434 ^c | 52.24 | 64.41 |
| | | 5 | 0.7630 | 55.74 | 66.48 |
| 1 | Clay soil spiked with motor oil | 7 | 0.5280 | 15.71 | 31.75 |
| 2 | Clay soil spiked with motor oil | 7 | 0.6672 | 20.64 | 34.53 |

^a The critical values of the Cochran statistic, the single Grubbs statistic, and the double Grubbs statistic are given in Reference 3.

^b Instrument 1 is the PE-FTIR spectrometer; instrument 2 is the BSci-IR spectrometer.

^c This value is an outlier.

We pooled all data generated with the five Isco SFE systems used in this study. Tables 22 and 23 present the summary statistics for these data before and after outlier removal. The method recoveries ranged from 79.4 to 122 percent for PE-FTIR analyses and from 72.7 to 113 percent for BSci-IR analyses. When the results from all laboratories using extraction vessels with a volume of 3.5 mL or less were pooled (seven laboratories), the mean recoveries ranged from 75.1 to 99.1 percent for PE-FTIR analyses and from 74.6 to 92.5 percent for BSci-IR analyses (Tables 24 and 25).

We performed a 2-sample t-test to determine whether the data in Tables 21, 23, and 25 (specifically, percent recoveries for each of the four matrices) are different from each other at the 5-percent significance level. When comparing the data in Tables 21 and 23, we find significant differences only for the clay soil matrix; when comparing the data in Tables 21 and 25, we find significant differences only for the SRS103-100 soil matrix; and finally, when comparing the data in Tables 23 and 25, we find significant differences for all matrices.

To determine the overall method recovery, we performed a linear regression of the data presented in Tables 12, 14, 16, and 18 (after removal of outliers). The measured concentrations for the PE-FTIR analyses were plotted on the y-axis, and the true concentrations of TPHs in the four matrices were plotted on the x-axis. The slope of the regression equation was 0.8291, and the intercept was -18.87. The correlation coefficient was 0.9121 (160 data points). These data indicate that the overall method recovery (for levels ranging from 614 mg/kg to 32,600 mg/kg) is 82.9 percent. Since the true concentrations for the three reference materials were obtained by extracting them with Freon-113 and analyzing the extract by IR, and the fourth true value was the spike concentration, we concluded that the performance of Method 3560/8440 (SFE/IR) was comparable to that of Method 9071A/8440 (Freon-113 extraction/IR analysis).

Method Precision

The interlaboratory standard deviation (s_R , reproducibility) is the precision associated with measurements generated by a group of laboratories; the single-analyst standard deviation (s_r , repeatability) is the precision associated with performance in an individual laboratory. The values for s_r and s_R are given in Tables 20 through 25; they were used to calculate the method precision, given as the repeatability relative standard deviation (RSD_r) and the reproducibility relative standard deviation (RSD_R).

When data from all laboratories were pooled and outliers removed, the RSD_r ranged from 11.5 to 17.0 percent for PE-FTIR analyses and from 11.1 to 18.2 percent for BSci-IR analyses. The RSD_R ranged from 17.3 to 45.4 percent for PE-FTIR analyses and from 16.7 to 47.9 percent for BSci-IR analyses.

When the data from the Isco SFE systems were pooled and outliers removed, the percent RSD_r ranged from 6.0 to 23.8 percent for the PE-FTIR analyses and from 9.8 to 22.6 percent for BSci-IR analyses. The RSD_R ranged from 6.0 to 27.1 percent for PE-FTIR analyses and from 12 to 26 percent for BSci-IR analyses. The interlaboratory method precision for the Isco SFE systems alone (RSD_R in Table 23) was better than that for all SFE systems (RSD_R in Table 25), because the Isco RSD_R values were lower for each of the four matrices.

The method precision for the seven laboratories using extraction vessels with a volume of 3.5 mL or less was almost similar to that achieved for all laboratories. For example, the RSD_r in Table 25 ranged from 8.1 to 21.1 percent for PE-FTIR analyses and from 7.4 to 20.6 percent for BSci-IR analyses, and the RSD_R ranged from 18.9 to 39.1 percent for PE-FTIR analyses and from 16.9 to 41.4 percent for BSci-IR analyses. This seems to indicate that the vessel size within the range used in this study was of little importance.

The precision estimates for the SFE/IR method are ± 20 percent for the intralaboratory performance and ± 45 percent for the interlaboratory performance.

Analysis of Variance

Analysis of variance (ANOVA) was used to separate (mathematically) the total variation of the experimental measurements into an intralaboratory portion and an interlaboratory portion, with a corresponding split of the total number of degrees of freedom. The mean squares for the among laboratories and within laboratories were calculated from the sum of squares and the number of degrees of freedom. The ratio of the two mean squares (the mean of squares among laboratories and the mean of squares within laboratories) is distributed as F. Thus, we calculated the F and compared it with the critical values of F for the upper 95 percent point of distribution at the corresponding degrees of freedom. For example, in Table 27, we are performing an ANOVA for data submitted by the 14 laboratories for the TPH-1 sample; in each case, we have three replicates per laboratory. The among-laboratories degrees of freedom is the number of laboratories (14 minus 1). The within-laboratories degrees of freedom is the number of laboratories multiplied by 2; the number 2 in this case is the number of replicates minus 1. The ratio of the mean squares is 19.4. Since the critical value of F is 2.09, we see that the variation from laboratory to laboratory was greater than that attributed to the analytical error displayed within laboratories. The matrix and operational parameters such as flow rate, extraction vessel design and orientation, mode of collection of the extracted material, and temperature of the collection solvent/trap seemed to be important.

The ANOVA results for the other three matrices are presented in Tables 28 through 30. In all cases, we find that the interlaboratory variance is greater than intralaboratory variance for the reasons stated above.

To further substantiate that the interlaboratory variance can be attributed to matrix, we pooled the data obtained with the Isco SFE systems (laboratories 05, 13, 14, 15, and 17). In this case, all laboratories used an extraction vessel in vertical position, and the collection vial was kept in a beaker with water at room temperature. The vessel dimensions, however, varied slightly, the restrictor dimensions varied slightly, and the restrictor temperature was different for one laboratory (Table 3). Nonetheless, the ANOVA results (Tables 31 through 34) indicate that the interlaboratory variance is significantly less than that attributed to analytical error displayed within laboratories for three of the matrices (TPH-2 soil, SRS103-100 soil, and spiked clay soil), but significantly greater than that attributed to analytical error within laboratories for the TPH-1 soil.

One-way ANOVA was also performed for the laboratories using vessels of 3.5 mL or less (Tables 35 through 38). In all four cases, we find that the interlaboratory variance is greater than that attributed to analytical error displayed within laboratories.

TABLE 27. ONE-WAY ANOVA FOR THE TPH-1 SOIL SAMPLES EXTRACTED BY ALL LABORATORIES

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|--------------|
| | X ₁ | X ₂ | X ₃ | |
| 01 | 618 | 808 | 1,000 | 2,430 |
| 02 | 449 | 453 | 487 | 1,390 |
| 03 | 67 | 10 | 77 | 154 |
| 04 | 352 | 391 | 420 | 1,160 |
| 05 | 943 | 859 | 854 | 2,660 |
| 06 | 250 | 478 | 44 | 772 |
| 08 | 737 | 799 | 933 | 2,470 |
| 10 | 693 | 1,070 | 789 | 2,550 |
| 11 | 972 | 1,210 | 1,100 | 3,280 |
| 12 | 691 | 701 | 656 | 2,050 |
| 13 | 673 | 604 | 620 | 1,900 |
| 14 | 639 | 605 | 580 | 1,820 |
| 15 | 883 | 963 | 885 | 2,730 |
| 17 | 575 | 712 | 832 | <u>2,120</u> |
| Grand sum | | | | 27,480 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F _{crit} ^a |
|------------------------|-----------------------|-------------------|----------------|------|--------------------------------|
| Among laboratories | 13 | 3,114,300 | 239,560 | 19.4 | 2.09 |
| Within laboratories | <u>28</u> | <u>345,970</u> | 12,360 | | |
| Total | 41 | 3,460,270 | | | |

^a F_{crit} is F_{0.05, 13, 28}

TABLE 28. ONE-WAY ANOVA FOR THE TPH-2 SOIL SAMPLES EXTRACTED BY ALL LABORATORIES

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|--------------|
| | X ₁ | X ₂ | X ₃ | |
| 01 | 1,840 | 2,160 | 1,880 | 5,880 |
| 02 | 1,940 | 1,990 | 1,730 | 5,660 |
| 03 | 257 | 410 | 356 | 1,020 |
| 04 | 1,350 | 1,160 | 1,520 | 4,030 |
| 05 | 1,490 | 2,000 | 2,230 | 5,720 |
| 06 | 1,850 | 1,730 | 1,380 | 4,960 |
| 08 | 1,480 | 1,100 | 1,500 | 4,080 |
| 10 | 1,880 | 2,240 | 1,980 | 6,100 |
| 11 | 2,140 | 2,200 | 2,050 | 6,390 |
| 12 | 2,320 | 2,440 | 2,110 | 6,870 |
| 13 | 1,900 | 2,010 | 1,940 | 5,850 |
| 14 | 2,130 | 1,900 | 1,340 | 5,370 |
| 15 | 1,850 | 1,900 | 1,990 | 5,740 |
| 17 | 1,780 | 1,920 | 1,980 | <u>5,680</u> |
| Grand sum | | | | 73,350 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square (variance) | F | F _{crit} |
|-------------------------------|-----------------------|-------------------|------------------------------|------|-------------------|
| Before outlier removal | | | | | |
| Among laboratories | 13 | 9,013,060 | 693,310 | 16.4 | 2.09 ^a |
| Within laboratories | <u>28</u> | <u>1,187,300</u> | <u>42,400</u> | | |
| Total | 41 | 10,200,360 | 248,790 | | |
| After outlier removal | | | | | |
| Among laboratories | 12 | 2,630,900 | 219,240 | 4.6 | 2.15 ^b |
| Within laboratories | <u>26</u> | <u>1,175,300</u> | 45,200 | | |
| Total | 38 | 3,806,200 | | | |

^a F_{crit} is F_{0.05, 13, 28}

^b F_{crit} is F_{0.05, 12, 26}

TABLE 29. ONE-WAY ANOVA FOR THE SRS103-100 SOIL SAMPLES EXTRACTED BY ALL LABORATORIES

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|---------------|
| | X ₁ | X ₂ | X ₃ | |
| 01 | 32,600 | 35,600 | 33,700 | 101,900 |
| 02 | 28,900 | 29,600 | 29,500 | 88,000 |
| 03 | 4,130 | 11,100 | — | 15,230 |
| 04 | 23,200 | 20,800 | 25,000 | 69,000 |
| 05 | 25,200 | 29,500 | 10,800 | 69,500 |
| 06 | 46,400 | 45,100 | 54,600 | 146,100 |
| 08 | 25,500 | 24,400 | 24,400 | 74,300 |
| 10 | 19,500 | 14,800 | 17,500 | 51,800 |
| 11 | 28,900 | 27,800 | 28,300 | 85,000 |
| 12 | 33,200 | 32,700 | — | 65,900 |
| 13 | 30,900 | 30,900 | 31,500 | 93,300 |
| 14 | 26,500 | 29,200 | 11,400 | 67,100 |
| 15 | 31,100 | 32,200 | 30,800 | 94,100 |
| 17 | 21,600 | 24,300 | 22,600 | <u>68,500</u> |
| | | | Grand sum | 1,085,730 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F _{crit} ^a |
|------------------------|-----------------------|--------------------|----------------|------|--------------------------------|
| Among laboratories | 13 | 3,052,860,000 | 234,835,400 | 13.6 | 2.09 |
| Within laboratories | <u>28</u> | <u>484,740,000</u> | 17,312,140 | | |
| Total | 41 | 3,537,600,000 | | | |

^a F_{crit} is F_{0.05, 13, 28}

TABLE 30. ONE-WAY ANOVA FOR THE SPIKED CLAY SOIL SAMPLES EXTRACTED BY ALL LABORATORIES

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|---------|
| | X ₁ | X ₂ | X ₃ | |
| 01 | 8,760 | 10,400 | 8,680 | 27,840 |
| 02 | 7,850 | 8,740 | 9,740 | 26,330 |
| 03 | 580 | 368 | 963 | 1,910 |
| 04 | 7,000 | 7,590 | 7,270 | 21,860 |
| 05 | 8,160 | 8,790 | 9,290 | 26,240 |
| 06 | 6,240 | 6,810 | 4,100 | 17,150 |
| 08 | 7,110 | 6,130 | 5,780 | 19,020 |
| 10 | 6,470 | 4,310 | 3,200 | 13,980 |
| 11 | 9,690 | 8,380 | 9,400 | 27,470 |
| 12 | 8,660 | 8,810 | 9,270 | 26,740 |
| 13 | 8,420 | 8,210 | 8,830 | 25,460 |
| 14 | 5,000 | 7,910 | 8,970 | 21,880 |
| 15 | 8,560 | 7,660 | 9,020 | 25,240 |
| 17 | 7,990 | 7,940 | 8,590 | 24,520 |
| Grand sum | | | | 305,640 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F _{crit} ^a |
|-------------------------------|-----------------------|-------------------|----------------|------|--------------------------------|
| Before outlier removal | | | | | |
| Among laboratories | 13 | 217,320,000 | 12,783,500 | 13.6 | 2.09 ^a |
| Within laboratories | <u>28</u> | <u>26,250,000</u> | 937,500 | | |
| Total | 41 | 243,570,000 | | | |
| After outlier removal | | | | | |
| Among laboratories | 12 | 74,870,000 | 6,239,200 | 6.2 | 2.15 ^a |
| Within laboratories | <u>26</u> | <u>26,070,000</u> | 1,002,700 | | |
| Total | 38 | 100,940,000 | | | |

^a F_{crit} is F_{0.05, 13, 28}

^b F_{crit} is F_{0.05, 12, 26}

TABLE 31. ONE-WAY ANOVA FOR THE TPH-1 SOIL SAMPLES EXTRACTED BY LABORATORIES 05, 13, 14, 15, AND 17

| Laboratory code | Concentration (mg/kg) | | | Sum | |
|--------------------|-----------------------|-------|-----------|--------|--|
| | X_1 | X_2 | X_3 | | |
| 05 | 943 | 859 | 854 | 2,660 | |
| 13 | 673 | 604 | 620 | 1,900 | |
| 14 | 639 | 605 | 580 | 1,820 | |
| 15 | 883 | 963 | 885 | 2,730 | |
| 17 | 575 | 712 | 832 | 2,120 | |
| | | | Grand sum | 11,230 | |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F_{crit}^a |
|------------------------|-----------------------|-------------------|----------------|------|--------------|
| Among laboratories | 4 | 239,780 | 59,950 | 12.9 | 3.48 |
| Within laboratories | <u>10</u> | <u>46,600</u> | 4,660 | | |
| Total | 14 | 286,380 | | | |

^a F_{crit} is $F_{0.05, 4, 10}$

TABLE 32. ONE-WAY ANOVA FOR THE TPH-2 SOIL SAMPLES EXTRACTED BY LABORATORIES 05, 13, 14, 15, AND 17

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|--------|
| | X ₁ | X ₂ | X ₃ | |
| 05 | 1,490 | 2,000 | 2,230 | 5,720 |
| 13 | 1,900 | 2,010 | 1,940 | 5,850 |
| 14 | 2,130 | 1,900 | 1,340 | 5,370 |
| 15 | 1,850 | 1,900 | 1,990 | 5,740 |
| 17 | 1,780 | 1,920 | 1,980 | 5,680 |
| | | | Grand sum | 28,360 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F _{crit} ^a |
|------------------------|-----------------------|-------------------|----------------|------|--------------------------------|
| Among laboratories | 4 | 43,290 | 10,820 | 0.16 | 3.48 |
| Within laboratories | <u>10</u> | <u>654,400</u> | 65,440 | | |
| Total | 14 | 697,690 | | | |

^a F_{crit} is F_{0.05, 4, 10}

TABLE 33. ONE-WAY ANOVA FOR THE SRS103-100 SOIL SAMPLES EXTRACTED BY LABORATORIES 05, 13, 14, 15, AND 17

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|---------|
| | X ₁ | X ₂ | X ₃ | |
| 05 | 25,200 | 29,500 | 10,800 | 65,500 |
| 13 | 30,900 | 30,900 | 31,500 | 93,300 |
| 14 | 26,500 | 29,200 | 11,400 | 67,100 |
| 15 | 31,100 | 32,200 | 30,800 | 94,100 |
| 17 | 21,600 | 24,300 | 22,600 | 68,500 |
| Grand sum | | | | 388,500 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F _{crit} ^a |
|------------------------|-----------------------|-------------------|----------------|-----|--------------------------------|
| Among laboratories | 4 | 286,050,000 | 71,512,500 | 1.9 | 3.48 |
| Within laboratories | 10 | 380,950,000 | 38,095,000 | | |
| Total | 14 | 667,000,000 | | | |

^a F_{crit} is F_{0.05, 4, 10}

TABLE 34. ONE-WAY ANOVA FOR THE SPIKED CLAY SOIL SAMPLES EXTRACTED BY LABORATORIES 05, 13, 14, 15, AND 17

| Laboratory code | Concentration (mg/kg) | | | Sum |
|-----------------|-----------------------|-------|-------|---------|
| | X_1 | X_2 | X_3 | |
| 05 | 8,160 | 8,790 | 9,290 | 26,240 |
| 13 | 8,420 | 8,210 | 8,830 | 25,460 |
| 14 | 5,000 | 7,910 | 8,970 | 21,880 |
| 15 | 8,560 | 7,660 | 9,020 | 25,240 |
| 17 | 7,990 | 7,940 | 8,590 | 24,520 |
| Grand sum | | | | 123,340 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F_{crit}^a |
|---------------------|--------------------|----------------|-------------|-----|--------------|
| Among laboratories | 4 | 3,740,160 | 940,720 | 0.9 | 3.48 |
| Within laboratories | 10 | 10,509,700 | 1,050,970 | | |
| Total | 14 | 14,249,860 | | | |

^a F_{crit} is $F_{0.05, 4, 10}$

TABLE 35. ONE-WAY ANOVA FOR THE TPH-1 SOIL SAMPLES EXTRACTED BY LABORATORIES 01, 04, 06, 08, 13, 14, AND 17

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|--------|
| | X ₁ | X ₂ | X ₃ | |
| 01 | 618 | 808 | 1,000 | 2,430 |
| 04 | 352 | 391 | 420 | 1,160 |
| 06 | 250 | 478 | 44 | 772 |
| 08 | 737 | 799 | 933 | 2,470 |
| 13 | 673 | 604 | 620 | 1,900 |
| 14 | 639 | 605 | 580 | 1,820 |
| 17 | 575 | 712 | 832 | 2,120 |
| | | | Grand sum | 12,670 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F _{crit} ^a |
|------------------------|-----------------------|-------------------|----------------|-----|--------------------------------|
| Among laboratories | 6 | 804,350 | 134,060 | 8.3 | 2.85 |
| Within laboratories | 14 | <u>227,060</u> | 16,220 | | |
| Total | 20 | 1,031,410 | | | |

^a F_{crit} is F_{0.05, 6, 14}

TABLE 36. ONE-WAY ANOVA FOR THE TPH-2 SOIL SAMPLES EXTRACTED BY LABORATORIES 01, 04, 06, 08, 13, 14, AND 17

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|-------|-----------|--------|
| | X_1 | X_2 | X_3 | |
| 01 | 1,840 | 2,160 | 1,880 | 5,880 |
| 04 | 1,350 | 1,160 | 1,520 | 4,030 |
| 06 | 1,850 | 1,730 | 1,380 | 4,960 |
| 08 | 1,480 | 1,100 | 1,500 | 4,080 |
| 13 | 1,900 | 2,010 | 1,940 | 5,850 |
| 14 | 2,130 | 1,900 | 1,340 | 5,370 |
| 17 | 1,780 | 1,920 | 1,980 | 5,680 |
| | | | Grand sum | 35,850 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F_{crit}^a |
|------------------------|-----------------------|-------------------|----------------|-----|--------------|
| Among laboratories | 6 | 1,260,630 | 210,100 | 3.6 | 2.85 |
| Within laboratories | 14 | 704,000 | 58,670 | | |
| Total | 20 | 1,964,630 | | | |

^a F_{crit} is $F_{0.05, 6, 14}$

TABLE 37. ONE-WAY ANOVA FOR THE SRS103-100 SOIL SAMPLES EXTRACTED BY LABORATORIES 01, 04, 06, 08, 13, AND 14^a

| Laboratory code | Concentration (mg/kg) | | | Sum |
|--------------------|-----------------------|----------------|----------------|---------|
| | X ₁ | X ₂ | X ₃ | |
| 01 | 32,600 | 35,600 | 33,700 | 101,900 |
| 04 | 23,200 | 20,800 | 25,000 | 69,000 |
| 06 | 46,400 | 45,100 | 54,600 | 146,100 |
| 08 | 25,500 | 24,400 | 24,400 | 74,300 |
| 13 | 30,900 | 30,900 | 31,500 | 93,300 |
| 14 | 26,500 | 29,200 | 11,400 | 67,100 |
| | | | Grand sum | 552,200 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F _{crit} ^b |
|------------------------|-----------------------|--------------------|----------------|------|--------------------------------|
| Among laboratories | 5 | 1,496,270,000 | 299,254,000 | 14.3 | 3.11 |
| Within laboratories | <u>12</u> | <u>251,640,000</u> | 20,970,000 | | |
| Total | 17 | 1,747,910,000 | | | |

^a Laboratory 17 was not included because it used a 10-mL vessel for the extractions.

^b F_{crit} is F_{0.05, 5, 12}

TABLE 38. ONE-WAY ANOVA FOR THE SPIKED CLAY SOIL SAMPLES EXTRACTED BY LABORATORIES 01, 04, 06, 08, 13, 14, AND 17

| Laboratory code | Concentration (mg/kg) | | | Sum |
|-----------------|-----------------------|--------|-------|---------|
| | X_1 | X_2 | X_3 | |
| 01 | 8,760 | 10,400 | 8,680 | 27,840 |
| 04 | 7,000 | 7,590 | 7,270 | 21,860 |
| 06 | 6,240 | 6,810 | 4,100 | 17,150 |
| 08 | 7,110 | 6,130 | 5,780 | 19,020 |
| 13 | 8,420 | 8,210 | 8,830 | 25,460 |
| 14 | 5,000 | 7,910 | 8,970 | 21,880 |
| 17 | 7,990 | 7,940 | 8,590 | 24,520 |
| Grand sum | | | | 157,730 |

| Source of variation | Degrees of freedom | Sum of squares | Mean square | F | F_{crit}^a |
|---------------------|--------------------|----------------|-------------|------|--------------|
| Among laboratories | 6 | 27,630,000 | 4,605,000 | 4.03 | 2.85 |
| Within laboratories | 14 | 16,000,000 | 1,142,900 | | |
| Total | 20 | 43,630,000 | | | |

^a F_{crit} is $F_{0.05, 6, 14}$

Method Performance for the Clay Soil Samples Spiked with Corn Oil and Reference Oil

Table 39 shows the concentrations of oil/grease and TPHs determined in the extracts derived from spiked clay soil samples by SFE with carbon dioxide. Oil/grease is defined as the material that has been extracted from the soil sample with supercritical carbon dioxide and collected in PCE, but has not been subjected to silica gel cleanup. The spiking level was 1,000 mg/kg for each the corn oil and the reference oil. The oil/grease data indicate that six of the 14 laboratories achieved recoveries ranging from 69 to 84 percent, five laboratories had recoveries ranging from 34 to 42 percent, two had recoveries of 16 and 18 percent, and one laboratory did not submit an extract. These recoveries may be biased low because the quantification of corn oil was done against the reference oil standard. The TPH recoveries in these samples, after silica gel cleanup of the extract, ranged from 65 to 86 percent for seven laboratories, from 28 to 54 percent for four laboratories, 10 percent for two laboratories, with one laboratory not submitting an extract.

These data are difficult to interpret. We would have expected (based on the results reported earlier) much better and more consistent recoveries. It is possible that low recoveries of the reference oil are partly due to the volatilization of isooctane or chlorobenzene during extraction. We were not able to correlate these results with the flow rates of carbon dioxide in Table 9 for sample 8. For example, the laboratories using HP systems (laboratories 10, 11, and 12), with a carbon dioxide flow rate of 2 mL/min, recovered 53.6, 83, and 9.5 percent, respectively.

The data in Table 40 were obtained on the same extracts, but the IR analyses were performed with a BSci-IR spectrometer. The recoveries were slightly lower than those given in Table 39, but they followed the same trend.

Method Performance for the Wet Clay Soil Samples

Tables 41 and 42 show the method performance data for the wet clay soil samples spiked with motor oil at 10,000 mg/kg. These samples were mixed with equal portions of anhydrous sodium sulfate immediately prior to extraction. The recoveries were above 30 percent for only two laboratories (30.7 and 37.7 percent), and eight laboratories had recoveries below 8 percent.

These results indicated that additional experimental work was needed to identify a better drying method. Experiments were, therefore, performed at MRI-CO using anhydrous magnesium sulfate and Hydromatrix (diatomaceous earth) as drying agents with clay soil samples spiked with motor oil at 10,000 mg/kg. The water content of the samples was varied from 10.6 percent to 40 percent (Table 43); the extractions were performed at 340 atm/80°C/30 min (dynamic). We also extracted spiked clay soil samples containing water at 10.6, 20, 30, and 40 percent to which no drying agent was added, but the extractions were performed at 450 atm and 150°C for 30 min (dynamic). The recovery data in Table 43 indicate that at 10.6 percent water, the recoveries obtained with the two drying agents (e.g., 96 percent for anhydrous magnesium sulfate and 99 percent for Hydromatrix) were identical with those achieved at higher pressure and temperature but with no drying agent (e.g., 97 percent). As the water content increased, the recoveries decreased, and they became much lower for the experiments performed without a drying agent. Since we noticed that the clay particles tended to clump when we mixed them with anhydrous magnesium sulfate, we also used anhydrous magnesium sulfate without mixing it with the clay soil sample, but adding it as a plug in the extraction vessel such that the carbon dioxide would flow through the sample first and then through the bed of anhydrous magnesium sulfate. The recoveries achieved in those experiments were 75 percent at 20 percent

TABLE 39. RECOVERIES (PERCENT) OF OIL/GREASE AND TPHs DETERMINED IN EXTRACTS OF CLAY SOIL SPIKED WITH CORN OIL AND REFERENCE OIL (PE-FTIR SPECTROMETER)

| Laboratory code | Percent recovery | |
|--------------------|------------------|------|
| | Oil/grease | TPHs |
| 01 | 71.5 | 71.8 |
| 02 | 36.6 | 28.5 |
| 03 | 42.4 | 76.5 |
| 04 | 34.1 | 31.0 |
| 05 | b | b |
| 06 | 35.3 | 32.2 |
| 08 | 16.0 | 10.1 |
| 10 | 39.5 | 53.6 |
| 11 | 78.0 | 83.0 |
| 12 | 18.1 | 9.5 |
| 13 | 73.5 | 78.3 |
| 14 | 84.0 | 85.8 |
| 15 | 70.5 | 65.2 |
| 17 | 69.0 | 68.4 |

^a The spiking level was 1,000 mg/kg for reference corn oil and 1,000 mg/kg for reference oil. Single determinations.

^b This laboratory did not submit an extract for this sample.

TABLE 40. RECOVERIES (PERCENT) OF OIL/GREASE AND TPHs DETERMINED IN EXTRACTS OF CLAY SOIL SPIKED WITH CORN OIL AND REFERENCE OIL (BSci-IR SPECTROMETER)^a

| Laboratory code | Percent recovery | |
|--------------------|------------------|------|
| | Oil/grease | TPHs |
| 01 | 55.0 | 56.4 |
| 02 | 30.6 | 25.7 |
| 03 | 38.4 | 71.3 |
| 04 | 25.8 | 23.0 |
| 05 | b | b |
| 06 | 29.0 | 26.8 |
| 08 | 12.5 | 9.1 |
| 10 | 30.3 | 52.6 |
| 11 | 62.0 | 69.0 |
| 12 | 14.3 | 10.9 |
| 13 | 53.5 | 65.8 |
| 14 | 68.0 | 70.7 |
| 15 | 53.5 | 54.3 |
| 17 | 55.5 | 56.4 |

^a The spiking level was 1,000 mg/kg for reference corn oil and 1,000 mg/kg for reference oil. Single determinations.

^b This laboratory did not submit an extract for this sample.

TABLE 41. CONCENTRATIONS (mg/kg) AND PERCENT RECOVERIES OF TPHs DETERMINED IN EXTRACTS OF WET CLAY SOIL SPIKED WITH MOTOR OIL (PE-FTIR SPECTROMETER)^a

| Laboratory code | Concentration (mg/kg) | Percent recovery |
|-----------------|-----------------------|------------------|
| 01 | 3,770 | 37.7 |
| 02 | 1,920 | 19.2 |
| 03 | 135 | 1.4 |
| 04 | 2,210 | 22.1 |
| 05 | 489 | 4.9 |
| 06 | 3,070 | 30.7 |
| 08 | 76.0 | 0.8 |
| 10 | 69.0 | 0.7 |
| 11 | 1,780 | 17.8 |
| 12 | 533 | 5.3 |
| 13 | 289 | 2.9 |
| 14 | 608 | 6.1 |
| 15 | 2,330 | 23.3 |
| 17 | 775 | 7.8 |

^a The clay soil was spiked with motor oil at 10,000 mg/kg. The water content was 30 percent. Single determinations.

TABLE 42. CONCENTRATIONS (mg/kg) AND PERCENT RECOVERIES OF TPHs DETERMINED IN EXTRACTS OF WET CLAY SOIL SPIKED WITH MOTOR OIL (BSci-IR SPECTROMETER)^a

| Laboratory code | Concentration (mg/kg) | Percent recovery |
|-----------------|-----------------------|------------------|
| 01 | 3,730 | 37.3 |
| 02 | 1,850 | 18.5 |
| 03 | 125 | 1.3 |
| 04 | 2,150 | 21.5 |
| 05 | 489 | 4.9 |
| 06 | 2,930 | 29.3 |
| 08 | 86.0 | 0.9 |
| 10 | 71.8 | 0.7 |
| 11 | 1,740 | 17.4 |
| 12 | 561 | 5.6 |
| 13 | 316 | 3.2 |
| 14 | 612 | 6.1 |
| 15 | 2,320 | 23.2 |
| 17 | 754 | 7.5 |

^a The clay soil was spiked with motor oil at 10,000 mg/kg. The water content was 30 percent. Single determinations.

water, but they dropped to approximately 25 to 27 percent for samples with 30 and 40 percent water (Table 43). It is possible that the low recoveries were due to restrictor plugging and subsequently to reduced flow rate of carbon dioxide.

When spiked clay soil containing 40 percent water was mixed with anhydrous magnesium sulfate (equal weights) and allowed to equilibrate overnight, or for 5 days, at room temperature, we obtained much higher recoveries of TPHs when extracting at 340 atm/80°C for 30 min (dynamic). The average recovery \pm one standard deviation of eight determinations (two sets of four samples extracted in parallel) was 84.7 ± 3.4 percent for overnight equilibration (14 hours) and 76.9 ± 6.1 percent for 5-day equilibration for clay soil spiked with motor oil (Table 44), and 74.3 percent for overnight equilibration and 77.1 ± 4.6 percent for 5-day equilibration for clay soil spiked with diesel oil (Table 44). We used in both experiments a plug of anhydrous magnesium sulfate (1.5 g) and crushed the clumps that were formed upon adding magnesium sulfate. A disposable glass pipet was used for this purpose, and the crushing was done directly in the extraction vessel (to minimize losses of the more volatile petroleum hydrocarbons).

Correlation between the PE-FTIR Data and the BSci-IR Data

The development work for this method has been performed on an FTIR instrument, while the method specified the use of a filter or fixed-wavelength instrument. Questions were raised regarding how to compensate for the multiplex advantage of the FTIR and about the capability of the FTIR for spectral subtraction, which cannot be done with a filter instrument. To address these concerns, we analyzed all extracts for the interlaboratory study on two IR systems. The measurement was done first on the PE-FTIR, then the IR cuvette was placed in the BSci-IR system and the measurement was taken. The BSci-IR measurements and the PE-FTIR measurements (as TPH concentrations in mg/kg for each matrix) were then plotted on the y-axis and the x-axis, respectively. Table 45 shows the slope, intercepts, and the correlation coefficients of the FTIR and the IR data by matrix. Excellent correlation was demonstrated, as shown by the values for the correlation coefficients. The slopes indicate that the measurements obtained with the BSci-IR were always lower than those obtained with the PE-FTIR system. We cannot explain the differences; however, because these differences were 17 percent or less (the slopes of the linear regression equations ranged from 0.8267 to 0.9900 for four matrixes), we concluded that the two sets of data were comparable.

QUALITY ASSURANCE/QUALITY CONTROL

The quality control data generated in this study are presented in Tables 46 through 56. They include all calibration data generated during the analysis of extracts (by instrument, by material used in calibration, and by date) (Tables 46 through 52), the results from the analyses of the unspiked clay soil samples (Table 53), the results of the analyses of the PCE blanks (Table 54), the results from the analysis of all system blanks submitted by the participating laboratories (Table 55), and the results of the sample storage study (Table 56).

TABLE 43. PERCENT RECOVERIES OF TPHs DETERMINED IN EXTRACTS OF WET CLAY SOIL SPIKED WITH MOTOR OIL (PE-FTIR SPECTROMETER)^a

| Percent water | Percent recovery | | | |
|---------------|--|-------------------|--------------------------------------|------------------------------|
| | Addition of anhydrous MgSO ₄ ^b | | Addition of Hydromatrix ^b | No drying agent ^c |
| | Plug only | Mixed with sample | | |
| 10.6 | -- | 95.9 ± 1.2 | 98.8 | 96.8 |
| 20 | 74.9 | 86.4 ± 9.2 | 72.2 | 74.8 ± 12.2 |
| 30 | 24.8 | 55.8 | 53.1 | 52.0 |
| 40 | 26.9 | 45.3 | 47.1 | 29.2 |

^a The value given is the average recovery ± one standard deviation (for triplicate determinations) or the average recovery of duplicate determinations.

^b The extraction was performed at 340 atm/80°C/30 min (dynamic).

^c The extraction was performed at 450 atm/150°C/30 min (dynamic).

TABLE 44. PERCENT RECOVERIES OF TPHs DETERMINED IN EXTRACTS OF WET CLAY SOIL (40 PERCENT WATER) SPIKED WITH MOTOR OIL OR DIESEL OIL (PE-FTIR SPECTROMETER)^a

| Storage Condition | Motor oil | Diesel oil |
|-------------------|-------------------------|-------------------------|
| 14 hours at 22°C | 84.7 ± 3.4 ^b | -- |
| 120 hours at 22°C | 76.9 ± 6.1 ^c | -- |
| 14 hours at 4°C | -- | 74.3 ^d |
| 120 hours at 4°C | -- | 77.1 ± 4.6 ^e |

^a The extractions were performed at 340 atm/80°C/30 min (dynamic). Each spiked sample (3 g) was mixed with 3 g of anhydrous magnesium sulfate and stored as indicated. For extraction, a plug of 1.5 g of anhydrous sodium sulfate was put into each extraction vessel at the outlet of the vessel.

^b The number of determinations was eight.

^c The number of determinations was seven.

^d Duplicate determination.

^e The number of determinations was three.

TABLE 45. CORRELATION BETWEEN THE PE-FTIR DATA AND THE BSci-IR DATA^a

| Matrix | Sample identification | Slope | Intercept | Correlation coefficient | Number of data points |
|---|-----------------------|--------|-----------|-------------------------|-----------------------|
| TPH-1 soil | 1 | 0.9900 | -29.5 | 0.9929 | 42 |
| TPH-2 soil | 2 | 0.8711 | 50.3 | 0.9953 | 42 |
| SRS103-100 | 3 | 0.8880 | 945 | 0.9975 | 40 |
| Spiked clay soil | 4,6,7 | 0.9767 | 534 | 0.9831 | 42 |
| Spiked clay soil (oil/grease analysis) | 8 | 0.7689 | 21.8 | 0.9927 | 13 |
| Spiked clay soil (TPH analysis) | 8 | 0.8271 | 13.9 | 0.9891 | 13 |
| Wet clay soil | 9 | 0.9727 | 12.4 | 0.9997 | 14 |

^a The BSci-IR data were plotted on the y-axis and the PE-FTIR data were plotted on the x-axis.

**TABLE 46. CALIBRATION DATA FOR THE PE-FTIR SPECTROMETER USING
REFERENCE OIL IN PCE AND THE 10-mm PATH-LENGTH IR CELL^a**

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|------------|--------|-----------|----------------------------|
| 04/03/92 | 10 | 0.0560 | 0.0018 | -0.0016 | 0.9996 |
| | 25 | 0.0406 | | | |
| | 50 | 0.0856 | | | |
| | 100 | 0.1797 | | | |
| | 250 | 0.4525 | | | |
| | 500 | 0.8672 | | | |
| 04/06/92 | 10 | 0.0225 | 0.0017 | 0.0314 | 0.9987 |
| | 25 | 0.0725 | | | |
| | 50 | 0.1220 | | | |
| | 100 | 0.2222 | | | |
| | 250 | 0.4736 | | | |
| | 500 | 0.8768 | | | |
| 04/07/92 | 10 | 0.0361 | 0.0018 | 0.0288 | 0.9992 |
| | 25 | 0.0748 | | | |
| | 50 | 0.1085 | | | |
| | 100 | 0.2109 | | | |
| | 250 | 0.4921 | | | |
| | 500 | 0.8992 | | | |
| 04/10/92 | 10 | -0.0058 | 0.0017 | 0.0121 | 0.9982 |
| | 25 | 0.0785 | | | |
| | 50 | 0.1056 | | | |
| | 100 | 0.1899 | | | |
| | 250 | 0.4491 | | | |
| | 500 | 0.8722 | | | |
| 04/16/92 | 10 | 0.0242 | 0.0017 | 0.0247 | 0.9992 |
| | 25 | 0.0718 | | | |
| | 50 | 0.1082 | | | |
| | 100 | 0.2102 | | | |
| | 250 | 0.4749 | | | |
| | 500 | 0.8842 | | | |
| 04/17/92 | 10 | 0.0239 | 0.0017 | 0.0242 | 0.9992 |
| | 25 | 0.0715 | | | |
| | 50 | 0.1070 | | | |
| | 100 | 0.2092 | | | |
| | 250 | 0.4746 | | | |
| | 500 | 0.8828 | | | |

(continued)

^a Silica gel was added to the calibration standards.

TABLE 46. (concluded)^a

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|------------|--------|-----------|----------------------------|
| 04/21/92 | 10 | 0.0240 | 0.0017 | 0.0247 | 0.9992 |
| | 25 | 0.0713 | | | |
| | 50 | 0.1089 | | | |
| | 100 | 0.2115 | | | |
| | 250 | 0.4759 | | | |
| | 500 | 0.8869 | | | |
| 04/23/92 | 10 | 0.0240 | 0.0017 | 0.0242 | 0.9993 |
| | 25 | 0.0706 | | | |
| | 50 | 0.1070 | | | |
| | 100 | 0.2087 | | | |
| | 250 | 0.4738 | | | |
| | 500 | 0.8813 | | | |
| 04/24/92 | 10 | 0.0236 | 0.0017 | 0.0245 | 0.9993 |
| | 25 | 0.0715 | | | |
| | 50 | 0.1101 | | | |
| | 100 | 0.2097 | | | |
| | 250 | 0.4725 | | | |
| | 500 | 0.8848 | | | |
| 04/28/92 | 10 | 0.0237 | 0.0017 | 0.0234 | 0.9993 |
| | 25 | 0.0704 | | | |
| | 50 | 0.1074 | | | |
| | 100 | 0.2089 | | | |
| | 250 | 0.4731 | | | |
| | 500 | 0.8850 | | | |
| 05/01/92 | 10 | 0.0230 | 0.0017 | 0.0226 | 0.9994 |
| | 25 | 0.0704 | | | |
| | 50 | 0.1077 | | | |
| | 100 | 0.2086 | | | |
| | 250 | 0.4743 | | | |
| | 500 | 0.8903 | | | |
| 05/05/92 | 10 | 0.0231 | 0.0017 | 0.0234 | 0.9993 |
| | 25 | 0.0701 | | | |
| | 50 | 0.1066 | | | |
| | 100 | 0.2089 | | | |
| | 250 | 0.4722 | | | |
| | 500 | 0.8814 | | | |

^a Silica gel was added to the calibration standards.

**TABLE 47. CALIBRATION DATA FOR THE BSci-IR SPECTROMETER USING
REFERENCE OIL IN PCE AND THE 10-mm PATH-LENGTH IR CELL^a**

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|------------|--------|-----------|----------------------------|
| 04/03/92 | 10 | 0.006 | 0.0016 | -0.0074 | 1.000 |
| | 25 | 0.033 | | | |
| | 50 | 0.073 | | | |
| | 100 | 0.156 | | | |
| | 250 | 0.400 | | | |
| | 500 | 0.799 | | | |
| 04/06/92 | 10 | 0.020 | 0.0016 | 0.0260 | 0.9992 |
| | 25 | 0.072 | | | |
| | 50 | 0.109 | | | |
| | 100 | 0.198 | | | |
| | 250 | 0.422 | | | |
| | 500 | 0.813 | | | |
| 04/07/92 | 10 | 0.031 | 0.0016 | 0.0197 | 0.9998 |
| | 25 | 0.067 | | | |
| | 50 | 0.093 | | | |
| | 100 | 0.183 | | | |
| | 250 | 0.433 | | | |
| | 500 | 0.826 | | | |
| 04/10/92 | 10 | -0.003 | 0.0016 | 0.0055 | 0.9985 |
| | 25 | 0.069 | | | |
| | 50 | 0.092 | | | |
| | 100 | 0.165 | | | |
| | 250 | 0.393 | | | |
| | 500 | 0.805 | | | |
| 04/16/92 | 10 | 0.015 | 0.0016 | 0.0112 | 0.9997 |
| | 25 | 0.058 | | | |
| | 50 | 0.089 | | | |
| | 100 | 0.178 | | | |
| | 250 | 0.412 | | | |
| | 500 | 0.806 | | | |
| 04/17/92 | 10 | 0.024 | 0.0016 | 0.0162 | 0.9998 |
| | 25 | 0.061 | | | |
| | 50 | 0.093 | | | |
| | 100 | 0.182 | | | |
| | 250 | 0.416 | | | |
| | 500 | 0.811 | | | |

(continued)

^a Silica gel was added to the calibration standards.

TABLE 47. (concluded)^a

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|------------|--------|-----------|----------------------------|
| 04/21/92 | 10 | 0.019 | 0.0016 | 0.0143 | 0.9998 |
| | 25 | 0.060 | | | |
| | 50 | 0.093 | | | |
| | 100 | 0.182 | | | |
| | 250 | 0.417 | | | |
| | 500 | 0.813 | | | |
| 04/23/92 | 10 | 0.019 | 0.0016 | 0.0147 | 0.9997 |
| | 25 | 0.061 | | | |
| | 50 | 0.093 | | | |
| | 100 | 0.181 | | | |
| | 250 | 0.415 | | | |
| | 500 | 0.809 | | | |
| 04/24/92 | 10 | 0.022 | 0.0016 | 0.0184 | 0.9997 |
| | 25 | 0.064 | | | |
| | 50 | 0.098 | | | |
| | 100 | 0.185 | | | |
| | 250 | 0.419 | | | |
| | 500 | 0.813 | | | |
| 04/28/92 | 10 | 0.012 | 0.0016 | 0.0075 | 0.9997 |
| | 25 | 0.054 | | | |
| | 50 | 0.085 | | | |
| | 100 | 0.172 | | | |
| | 250 | 0.404 | | | |
| | 500 | 0.796 | | | |
| 05/01/92 | 10 | 0.022 | 0.0016 | 0.0179 | 0.9997 |
| | 25 | 0.063 | | | |
| | 50 | 0.096 | | | |
| | 100 | 0.187 | | | |
| | 250 | 0.420 | | | |
| | 500 | 0.815 | | | |
| 05/05/92 | 10 | 0.022 | 0.0016 | 0.0190 | 0.9996 |
| | 25 | 0.068 | | | |
| | 50 | 0.095 | | | |
| | 100 | 0.185 | | | |
| | 250 | 0.417 | | | |
| | 500 | 0.810 | | | |

^a Silica gel was added to the calibration standards.

TABLE 48. CALIBRATION DATA FOR THE PE-FTIR SPECTROMETER USING MOTOR OIL IN PCE AND THE 10-mm PATH-LENGTH IR CELL^a

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|------------|--------|-----------|----------------------------|
| 04/14/92 | 10 | 0.0306 | 0.0022 | 0.0263 | 0.9984 |
| | 25 | 0.0701 | | | |
| | 50 | 0.1323 | | | |
| | 100 | 0.2576 | | | |
| | 250 | 0.6152 | | | |
| | 500 | 1.0992 | | | |
| 04/17/92 | 10 | 0.0237 | 0.0022 | 0.0175 | 0.9987 |
| | 25 | 0.0629 | | | |
| | 50 | 0.1236 | | | |
| | 100 | 0.2499 | | | |
| | 250 | 0.6050 | | | |
| | 500 | 1.0987 | | | |
| 04/22/92 | 10 | 0.0252 | 0.0021 | 0.0281 | 0.9963 |
| | 25 | 0.0640 | | | |
| | 50 | 0.1257 | | | |
| | 100 | 0.2515 | | | |
| | 250 | 0.6159 | | | |
| | 500 | 1.0310 | | | |
| 04/30/92 | 10 | 0.0251 | 0.0022 | 0.0189 | 0.9987 |
| | 25 | 0.0640 | | | |
| | 50 | 0.1258 | | | |
| | 100 | 0.2514 | | | |
| | 250 | 0.6072 | | | |
| | 500 | 1.1016 | | | |
| 05/05/92 | 10 | 0.0255 | 0.0022 | 0.0182 | 0.9989 |
| | 25 | 0.0643 | | | |
| | 50 | 0.1249 | | | |
| | 100 | 0.2508 | | | |
| | 250 | 0.6020 | | | |
| | 500 | 1.1023 | | | |

^a Silica gel was added to the calibration standards.

TABLE 49. CALIBRATION DATA FOR THE BSci-IR SPECTROMETER USING MOTOR OIL IN PCE AND THE 10-mm PATH-LENGTH IR CELL^a

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|------------|--------|-----------|----------------------------|
| 04/14/92 | 10 | 0.030 | 0.0020 | 0.0115 | 0.9999 |
| | 25 | 0.057 | | | |
| | 50 | 0.109 | | | |
| | 100 | 0.215 | | | |
| | 250 | 0.516 | | | |
| | 500 | 1.002 | | | |
| | | | | | |
| 04/17/92 | 10 | 0.020 | 0.0020 | 0.0031 | 1.000 |
| | 25 | 0.052 | | | |
| | 50 | 0.103 | | | |
| | 100 | 0.206 | | | |
| | 250 | 0.506 | | | |
| | 500 | 1.000 | | | |
| | | | | | |
| 04/22/92 | 10 | 0.020 | 0.0020 | 0.0014 | 1.000 |
| | 25 | 0.050 | | | |
| | 50 | 0.101 | | | |
| | 100 | 0.205 | | | |
| | 250 | 0.502 | | | |
| | 500 | 1.001 | | | |
| | | | | | |
| 04/30/92 | 10 | 0.021 | 0.0020 | 0.0024 | 1.000 |
| | 25 | 0.051 | | | |
| | 50 | 0.101 | | | |
| | 100 | 0.206 | | | |
| | 250 | 0.505 | | | |
| | 500 | 1.001 | | | |
| | | | | | |
| 05/05/92 | 10 | 0.021 | 0.0020 | 0.0037 | 1.000 |
| | 25 | 0.053 | | | |
| | 50 | 0.103 | | | |
| | 100 | 0.207 | | | |
| | 250 | 0.506 | | | |
| | 500 | 1.002 | | | |
| | | | | | |

^a Silica gel was added to the calibration standards.

TABLE 50. CALIBRATION DATA FOR THE PE-FTIR SPECTROMETER USING REFERENCE OIL IN PCE AND THE 10-mm PATH-LENGTH IR CELL (NO SILICA GEL CLEANUP)^a

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|--------------|--------|-----------|----------------------------|
| 04/16/92 | 10 | 0.0305 | 0.0017 | 0.0192 | 0.9989 |
| | 25 | 0.0559 | | | |
| | 50 | 0.0977 | | | |
| | 100 | 0.1959 | | | |
| | 250 | 0.4784 | | | |
| | 500 | 0.8652 | | | |
| 04/23/92 | 10 | 0.0225 | 0.0017 | 0.0146 | 0.9992 |
| | 25 | 0.0497 | | | |
| | 50 | 0.0960 | | | |
| | 100 | 0.1936 | | | |
| | 250 | 0.4597 | | | |
| | 500 | 0.8494 | | | |
| 04/28/92 | 10 | 0.0202 | 0.0017 | 0.0121 | 0.9992 |
| | 25 | 0.0479 | | | |
| | 50 | 0.0933 | | | |
| | 100 | 0.1897 | | | |
| | 250 | 0.4569 | | | |
| | 500 | 0.8461 | | | |
| 05/01/92 | 10 | 0.0214 | 0.0017 | 0.0132 | 0.9992 |
| | 25 | 0.0484 | | | |
| | 50 | 0.0941 | | | |
| | 100 | 0.1920 | | | |
| | 250 | 0.4576 | | | |
| | 500 | 0.8470 | | | |
| 05/05/92 | 10 | 0.0208 | 0.0017 | 0.0101 | 0.9993 |
| | 25 | 0.0476 | | | |
| | 50 | 0.0944 | | | |
| | 100 | ^b | | | |
| | 250 | 0.4558 | | | |
| | 500 | 0.8468 | | | |

^a No silica gel was added to the calibration standards.

^b Calibration was not performed with the 100- $\mu\text{g/mL}$ standard.

TABLE 51. CALIBRATION DATA FOR THE BSci-IR SPECTROMETER USING REFERENCE OIL IN PCE AND THE 10-mm PATH-LENGTH IR CELL (NO SILICA GEL CLEANUP)^a

| Date | Concentration ($\mu\text{g/mL}$) | Absorbance | Slope | Intercept | Correlation coefficient |
|----------|---------------------------------------|------------|--------|-----------|----------------------------|
| 04/16/92 | 10 | 0.027 | 0.0016 | 0.0081 | 0.9998 |
| | 25 | 0.046 | | | |
| | 50 | 0.084 | | | |
| | 100 | 0.171 | | | |
| | 250 | 0.429 | | | |
| | 500 | 0.820 | | | |
| 04/23/92 | 10 | 0.019 | 0.0016 | 0.0046 | 0.9999 |
| | 25 | 0.042 | | | |
| | 50 | 0.083 | | | |
| | 100 | 0.169 | | | |
| | 250 | 0.411 | | | |
| | 500 | 0.804 | | | |
| 04/28/92 | 10 | 0.018 | 0.0016 | 0.0049 | 0.9999 |
| | 25 | 0.043 | | | |
| | 50 | 0.083 | | | |
| | 100 | 0.169 | | | |
| | 250 | 0.410 | | | |
| | 500 | 0.801 | | | |
| 05/01/92 | 10 | 0.019 | 0.0016 | 0.0051 | 0.9999 |
| | 25 | 0.043 | | | |
| | 50 | 0.083 | | | |
| | 100 | 0.171 | | | |
| | 250 | 0.412 | | | |
| | 500 | 0.806 | | | |
| 05/05/92 | 10 | 0.019 | 0.0016 | 0.0032 | 0.9999 |
| | 25 | 0.041 | | | |
| | 50 | 0.083 | | | |
| | 100 | b | | | |
| | 250 | 0.411 | | | |
| | 500 | 0.804 | | | |

^a No silica gel was added to the calibration standards.

^b Calibration was not performed with the 100- $\mu\text{g/mL}$ standard.

TABLE 52. CONCENTRATIONS OF DAILY STANDARDS ANALYZED TO VERIFY SYSTEM REPRODUCIBILITY

| Date of analysis | Concentration ($\mu\text{g/mL}$) ^a | |
|------------------|---|---------|
| | PE-FTIR | BSci-IR |
| 04/06/92 | 110 | 103 |
| 04/06/92 | 104 | 108 |
| 04/06/92 | 104 | 108 |
| 04/09/92 | 100 | 103 |
| 04/09/92 | 100 | 100 |
| 04/10/92 | 105 | 106 |
| 04/10/92 | 107 | 108 |
| 04/10/92 | 112 | 115 |
| 04/10/92 | 107 | 111 |
| 04/14/92 | 112 | 108 |
| 04/14/92 | 112 | 112 |
| 04/14/92 | 106 | 106 |
| 04/14/92 | 101 | 100 |
| 04/16/92 | 107 | 105 |
| 04/16/92 | 100 | 102 |
| 04/16/92 | 107 | 105 |
| 04/16/92 | 105 | 103 |
| 04/17/92 | 105 | 103 |
| 04/21/92 | 106 | 107 |
| 04/21/92 | 104 | 107 |
| 04/23/92 | 105 | 103 |
| 04/23/92 | 104 | 103 |
| 04/23/92 | 106 | 106 |
| 04/24/92 | 107 | 103 |
| 04/28/92 | 106 | 105 |
| 04/30/92 | 105 | 104 |
| 04/30/92 | 105 | 103 |
| 05/01/92 | 106 | 102 |
| 05/01/92 | 105 | 108 |
| 05/01/92 | 104 | 109 |
| 05/05/92 | 107 | 104 |

^a The true concentration of the calibration standard was 100 $\mu\text{g/mL}$.

TABLE 53. CONCENTRATIONS OF TPHs IN THE EXTRACTS FROM THE UNSPIKED CLAY SOIL SAMPLES

| Laboratory | Concentration (mg/kg) | |
|------------|-----------------------|---------|
| | PE-FTIR | BSci-IR |
| 01 | 207 | 204 |
| 02 | < 10 | < 10 |
| 03 | 492 | 469 |
| 04 | 43 | 60 |
| 05 | 306 | 305 |
| 06 | 103 | 88 |
| 08 | a | a |
| 10 | 423 | 422 |
| 11 | 204 | 196 |
| 12 | 58 | 54 |
| 13 | 325 | 327 |
| 14 | 398 | 400 |
| 15 | 328 | 332 |
| 17 | 251 | 250 |

^a Extract not submitted for analysis.

TABLE 54. CONCENTRATIONS OF TPHs IN THE PCE BLANKS ANALYZED DURING THIS STUDY

| Date of analysis | Number of blanks performed during day | Concentration ($\mu\text{g/mL}$ PCE) | |
|------------------|---------------------------------------|---------------------------------------|---------|
| | | PE-FTIR | BSci-IR |
| 04/03/92 | 3 | < 10 | < 10 |
| 04/06/92 | 5 | < 10 | < 10 |
| 04/09/92 | 3 | < 10 | < 10 |
| 04/10/92 | 4 | < 10 | < 10 |
| 04/14/92 | 5 | < 10 | < 10 |
| 04/16/92 | 2 | < 10 | < 10 |
| 04/17/92 | 1 | < 10 | < 10 |
| 04/21/92 | 2 | < 10 | < 10 |
| 04/23/92 | 2 | < 10 | < 10 |
| 04/24/92 | 1 | < 10 | < 10 |
| 04/28/92 | 1 | < 10 | < 10 |
| 04/30/92 | 2 | < 10 | < 10 |
| 05/01/92 | 3 | < 10 | < 10 |
| 05/05/92 | 3 | < 10 | < 10 |

TABLE 55. CONCENTRATIONS OF TPHs IN THE EXTRACTS SUBMITTED AS SYSTEM BLANKS

| Laboratory | Number of blanks submitted | Concentration (mg/mL) | |
|------------|-------------------------------|---------------------------------|----------------------------------|
| | | PE-FTIR | BSci-IR |
| 01 | 0 | a | a |
| 02 | 1 | < 10 | < 10 |
| 03 | 2 | < 10; < 10 | < 10; < 10 |
| 04 | 1 | < 10 | < 10 |
| 05 | 2 | < 10; < 10 | < 10; < 10 |
| 06 | 0 | a | a |
| 08 | 1 | < 10 | < 10 |
| 10 | 0 | a | a |
| 11 | 5 | 229; 18; 96; < 10; 37 | 238; 22; 98 27; 41 |
| 12 | 0 | a | a |
| 13 | 1 | 81 | 62 |
| 14 | 4 | 22; < 10; < 10; < 10 | 18; < 10; < 10; 11 |
| 15 | 4 | 10; < 10; < 10; < 10 | 14; < 10; < 10; < 10 |
| 17 | 5 | < 10; < 10; < 10; < 10; < 10 | < 10; < 10; < 10; < 10; < 10; |

^a Extract not submitted for analysis.

TABLE 56. PERCENT RECOVERIES OF TPHs FROM THE SPIKED CLAY SOIL SAMPLES, STORED AT 4°C IN THE DARK, AS A FUNCTION OF TIME

| Days of storage | PE-FTIR | | BSci-IR | |
|--------------------|---------|------|---------|------|
| 22 | 83.1 | 79.1 | 87.2 | 82.2 |
| 34 | 85.6 | 76.6 | 87.7 | 81.2 |
| 40 | 87.4 | 77.9 | 90.7 | 81.7 |

^a The clay soil samples were spiked with motor oil at 10,000 mg/kg. The values given are for duplicate determinations at each time.

The following conclusions can be drawn from these data:

- The PE-FTIR and the BSci-IR calibration data generated during this study agreed to within 10 percent. The slopes for the reference oil calibrations were 0.0017 for the PE-FTIR system and 0.0016 for the BSci-IR system, and the correlation coefficients for standards ranging from 10 to 500 $\mu\text{g/mL}$ were 0.999 or greater (Tables 46 and 47). The slopes for the motor oil calibrations were 0.0022 for the PE-FTIR system (Table 48) and 0.0020 for the BSci-IR system (Table 49).
- Both IR spectrometer systems gave reproducible results over the period of 1 month during which IR analyses were performed. The multilevel calibration data (slopes, intercepts, and correlation coefficients) from each instrument indicate deviations less than 10 percent for the IR determinations performed over a period of 1 month. A 100- $\mu\text{g/mL}$ reference oil or motor oil standard was analyzed after every 10 analyses. The largest deviation (+15 percent) was found on one day with the BSci-IR system. The other deviations averaged approximately 5 percent over a period of 1 month (Table 52).
- The extracts from the unspiked clay soil samples (sample 5) did not contain TPHs above 10 mg/kg. We analyzed this material in our laboratory on several occasions and found TPHs at levels ranging from 5 to 9 mg/kg. The sample-5 extracts submitted by the participating laboratories, however, contained TPHs at concentrations as high as 492 mg/kg (Table 53). This implies that the levels reported for sample 5 are due to cross-contamination from the previous extraction, since we instructed the laboratories to analyze the clay soil blank immediately after the spiked clay soil (concentration 10,000 mg/kg). Thus, contamination of the SFE system, especially when dealing with high-concentration samples, is likely, and the analyst must take the necessary precautions to minimize it (e.g., clean extraction vessel and frits, replace restrictor, perform system blanks).
- The Aldrich PCE (spectrophotometric grade) was acceptable for this study. The solvent blanks analyzed during the study indicated TPH levels of less than 10 $\mu\text{g/mL}$ (Table 54).
- The system blanks generated with the Aldrich PCE and the SFE-grade carbon dioxide did not show contamination, with the exception of laboratories 11 and 13 (Table 55). Since both laboratories used the Scott SFE-grade carbon dioxide, and since the PCE was provided by MRI-CO, it appears that the source of contamination was in the SFE system.
- The percent recoveries from the spiked clay soil samples, that had been stored at 4°C in the dark and extracted after 22, 34, and 44 days, were independent of the storage time.

REFERENCES

1. Lopez-Avila, V., N. S. Dodhiwala, J. Benedicto, and R. Young, "SFE/IR Method for the Determination of Petroleum Hydrocarbons in Soils and Sediments," EPA 600/X-92/046, April 1992 (W.F. Beckert, Project Officer, EMSL-LV).
2. Lopez-Avila, V., J. Benedicto, N. S. Dodhiwala, and W. F. Beckert, "Development of an Off-Line SFE/IR Method for Petroleum Hydrocarbons in Soils," J. Chromatogr. Sci. 30, 335-343, 1992.
3. Guidelines for Collaborative Study Procedure to Validate Characteristics of a Method of Analysis, J. Assoc. Off. Anal. Chem. 72, 694-704, 1989.
4. "Lotus Spreadsheet Program for the Calculation of Performance Parameters from Collaborative Study Data Including Outlier Analyses," revision 3/5/91, received from John G. Phillips, Chairman, AOAC Statistics Committee.
5. EPA Method 418.1, "Total Recoverable Petroleum Hydrocarbons," in "Methods for Chemical Analyses of Water and Wastewater," EPA 600/14-79/020, revised 1983.

APPENDIX A

PROPOSED DRAFT PROTOCOL FOR SUPERCRITICAL FLUID EXTRACTION OF PETROLEUM HYDROCARBONS (METHOD 3560)

METHOD 3560

SUPERCRITICAL FLUID EXTRACTION OF TOTAL RECOVERABLE PETROLEUM HYDROCARBONS

1.0 SCOPE AND APPLICATION

1.1 Method 3560 describes the extraction with supercritical fluids of total recoverable petroleum hydrocarbons (TPHs) from soils, sediments, fly ash, and other solid materials that are amenable to extraction with conventional solvents. The method is suitable for use with any supercritical fluid extraction (SFE) system that allows extraction conditions (e.g., pressure, temperature, flow rate) to be adjusted to achieve separation of the TPHs from the matrices of concern.

1.2 Method 3560 is not suitable for the extraction of low-boiling TPHs such as gasoline.

2.0 SUMMARY OF METHOD

2.1 A known amount of sample is transferred to the extraction vessel. The sample is then extracted in the dynamic mode for up to 30 min with supercritical carbon dioxide at 340 atm and 80°C and a gas flow rate of 500 to 1,000 mL/min. After depressurization of the carbon dioxide, the extracted TPHs are collected in 3 mL of tetrachloroethylene or other appropriate solvent (see Section 5.3), or on a sorbent material, depending on the SFE system used. In the latter case, the analytes are collected by rinsing the sorbent material with tetrachloroethylene or other suitable solvent. After collection, the TPHs are analyzed by the appropriate determinative method.

3.0 INTERFERENCES

3.1 The analyst must demonstrate through the analysis of reagent blanks (collection solvent treated as per Section 7.4) that the supercritical fluid extraction system is free of interferents. To do this, perform a simulated extraction using an empty extraction vessel and a known amount of carbon dioxide under the same conditions as those used for sample extraction, and determine the background contamination by analyzing the extract by the appropriate determinative method (e.g., Method 8015 or 8440). If glass wool and a drying agent are used with the sample, they should be included when performing a reagent blank check.

3.2 The extraction vessel(s), the frits, the restrictor(s), and the multiport valve may retain solutes whenever high-concentration samples are extracted. It is, therefore, good practice to clean the extraction system after each extraction. Replacement of the restrictor may be necessary when reagent blanks indicate carryover. At least one reagent blank should be prepared and analyzed daily when the instrument is in use. Furthermore, reagent blanks should be prepared and analyzed after each extraction of a high-concentration sample (concentration in the high ppm range). If reagent blanks continue to indicate contamination even after replacement of the extraction vessel and the restrictor, the multi-port valve must be cleaned.

4.0 APPARATUS AND MATERIALS

4.1 Supercritical fluid extractor and associated hardware.

WARNING - A safety feature to prevent overpressurization is required on the extractor. This feature should be designed to protect the laboratory personnel and the instrument from possible injuries or damage resulting from equipment failure under high pressure.

4.1.2 Extraction vessel -- Stainless-steel vessel with end fittings and 0.5- or 2- μ m frits. Use the extraction vessel recommended by the manufacturer of the SFE system being used. The volume of the extraction vessel should fit the sample size. PEEK (polyether ether ketone) extraction vessels are acceptable only for use with specifically designed instruments.

4.1.2.1 Fittings used for the extraction vessel must be capable of withstanding the required extraction pressures. The maximum operating pressure for most extractors is 500 atm; however, extractors with higher pressure ratings are available. Check with the manufacturer of the particular extraction system on the maximum operating pressure and temperature for that system. Make sure that the extraction vessels are rated for such pressures and temperatures.

4.1.3 Restrictor -- 50- μ m ID x 150- or 375- μ m OD x 25- to 60-cm length piece of uncoated fused-silica tubing (J&W Scientific or equivalent). Other restrictors may be used including tapered restrictors, static pinhole restrictors, frit restrictors, and variable orifice restrictors (manual and computer-controlled) or crimped metal tubing. Check with the manufacturer of the SFE system on the advantages and disadvantages of the various restrictor designs.

4.1.4 Collection device -- The extracted TPHs can be collected either in vials containing solvent, or they can be trapped on a sorbent material (e.g., octadecyl-bonded silica, stainless-steel beads).

4.1.4.1 When the analytes are collected in a solvent, install the restrictor through a hole made through the cap and septum of the vial, and position the restrictor end about 0.5 inch from the bottom of the vial. A syringe needle should also be inserted through the septum of the vial (with the tip positioned just below the septum) to prevent buildup of pressure in the vial. Use the type of vials appropriate for the SFE system used.

4.1.4.2 When the analytes are trapped on a sorbent material, it is important to ensure that breakthrough of the analytes from the trap does not occur. Desorption from the trapping medium can be accomplished by increasing the temperature of the trap and using a solvent to remove the analytes. To recover the analytes, use the conditions suggested by the manufacturer of the particular system used.

4.2 Carbon dioxide cylinder balance (optional) -- Balances from White Associates, catalog no. 30, Scott Specialty Gases Model 5588D, or equivalent, can be used to monitor the fluid usage. Such a device is useful because carbon dioxide tanks used for SFE are not equipped with regulators, and it is difficult to determine when the tank needs to be replaced.

4.3 Tools required include: screwdriver (flat-blade), adjustable wrench, pliers, tubing cutter, and various small open-end wrenches for small fittings.

4.4 Other materials -- Magnesium sulfate monohydrate can be used as received. The following materials require high-temperature treatment (muffling at 400°C for 2 to 4 hours) prior to use since they may be contaminated with petroleum hydrocarbons. Alternatively, a blank can be performed to determine whether these materials are sufficiently clean.

4.4.1 Silanized glass wool.

4.4.2 Drying agents such as anhydrous magnesium sulfate or diatomaceous earth.

5.0 REAGENTS

5.1 Carbon dioxide, CO₂ -- Either supercritical fluid chromatography (SFC-grade) or SFE-grade CO₂ is acceptable for use in SFE. Aluminum cylinders are preferred over steel cylinders. The cylinders are fitted with eductor tubes, and their contents are under 1500 psi of helium head pressure.

5.2 Carbon dioxide (CO₂) for cryogenic cooling -- Certain parts of some models of extractors (i.e., the high-pressure pump head and the analyte trap) must be cooled during use. The carbon dioxide used for this purpose must be dry (< 50 ppm water content), and it must be supplied in tanks with a full-length eductor tube.

5.3 Tetrachloroethylene, C₂Cl₄ (spectrophotometric grade) -- Used for the collection of TPHs for determination by IR. Analyze a blank to ensure no interferences are present at the TPH wavelengths. Chlorofluorocarbons are not suitable for use with this method because of the risk to the ozone layer.

5.4 Other appropriate pesticide-quality solvents may be used for the collection of TPHs for determination by GC (i.e., methylene chloride). Chlorofluorocarbons are not suitable for use with this method because of the risk to the ozone layer.

5.5 Copper filings -- Copper filings added to remove elemental sulfur must have a shiny bright appearance to be effective. To remove oxides from copper surfaces, treat with dilute nitric acid, rinse with reagent water to remove all traces of acid, rinse with acetone (copper will darken if acid is still present), and dry under a stream of nitrogen.

6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

6.1 Solid samples should be collected and stored in the same manner as any other solid samples containing semivolatile organics. See Chapter Four for guidance relating to semivolatile organics (including holding times).

7.0 PROCEDURE

7.1 Determination of sample percent dry weight -- In certain cases, sample results are desired based on dry-weight basis. When such data are desired, a portion of sample for this determination should be weighed out at the same time as the portion used for analytical determination.

WARNING: The drying oven should be contained in a hood or vented.. Significant laboratory contamination may result from a heavily contaminated hazardous waste sample.

7.1.1 Immediately after weighing the sample portion for extraction, weigh 5.00 to 10.00 g of the remaining sample into a tared crucible. Determine the percent dry weight of the sample by drying it overnight at 105°C. Allow it to cool in a desiccator before weighing. Calculate the percent dry weight as follows:

$$\% \text{ dry weight} = \frac{\text{g of dry sample} \times 100}{\text{g of sample}}$$

7.2 Safety considerations -- Read section 11.0 "Safety" before attempting this procedure.

7.3 Sample handling

7.3.1 Decant and discard any aqueous layer that has accumulated on a sediment sample. Mix the sample thoroughly, especially composited samples. Discard any foreign objects such as sticks, leaves, and rocks.

7.3.2 Weigh 3 g of sample into a precleaned aluminum dish. A drying agent (e.g., anhydrous magnesium sulfate or diatomaceous earth) may be added to samples that contain water in excess of 20% to increase porosity or to bind water. Alternatively, magnesium sulfate monohydrate is an excellent drying agent, and the amount of heat released (compared with anhydrous magnesium sulfate) is small, thereby minimizing the loss of volatile petroleum hydrocarbons. The amount of the drying agent will depend on the water content of the sample. Typically, a ratio of 1:1 works well for wet soils and sediment materials. However, a certain amount of water (up to 10%) in the sample has been shown to improve recoveries from certain matrices; therefore, if the sample is dry, water may optionally be added to bring the water content to approximately 20%.

7.3.2.1 If drying agent has been added to the sample, store the mixture of sample and drying agent for several hours (preferably overnight) at 4°C, with a minimum of headspace. This additional storage time is necessary to achieve acceptable analyte recovery.

7.3.3 Transfer the weighed sample to a clean extraction vessel. The volume of the extraction vessel should match the sample volume. Use two plugs of silanized glass wool to hold the sample in place and fill the void volume (the use of drying agent or clean sand after the second glass wool plug to fill the void volume is also recommended). Attach the end fittings, and install the extraction vessel in the oven. Always use clean frits for each extraction vessel.

7.4 Sample extraction

7.4.1 Fill the collection vessel with 3 mL of tetrachloroethylene or other appropriate collection solvent. Chlorofluorocarbons are not suitable for use with this method because of the risk to the ozone layer.

7.4.2 Set the pressure at 340 atm and the temperature at 80°C. Follow manufacturer's instructions in setting up the instrument. Extract for 30 minutes in the dynamic mode. Note the safety precautions in Section 11 on venting the instrument into a chemical fume hood.

7.4.3 After the extraction time has elapsed, the system should automatically go to the equilibrate mode. At this point, remove the collection vessel(s) containing the extract(s). Since the depressurization of the carbon dioxide at the end of the restrictor outlet results in a gas flow rate of about 500 to 1000 mL/min, part of the collection solvent will evaporate during the extraction. However, cooling caused by the rapid expansion of the carbon dioxide limits the loss of solvent, so that approximately 2 mL remains (when tetrachloroethylene is used) after a 30-min extraction. To prevent the collection solvent from freezing during the extraction, place the collection vial in a beaker with warm water (approximately 25°C). The extract is then brought to the desired volume, or concentrated further. See Method 3510 for concentration techniques by micro Kuderna-Danish or nitrogen blowdown. Concentration must be performed in a chemical fume hood to prevent contamination of the laboratory environment.

7.4.4 Record the volume of liquid carbon dioxide used for extraction. Calculate the average flow rate by dividing the volume of the carbon dioxide by the extraction time.

7.4.5 The extract is ready for analysis by Method 8015 -- Nonhalogenated Volatile Organics by Gas Chromatography, or by Method 8440 -- Total Recoverable Petroleum Hydrocarbons by Infrared Spectrophotometry.

7.5 SFE System Maintenance

7.5.1 Depressurize the system following manufacturer's instructions.

7.5.2 After extraction of an especially tarry sample, the frits may require replacement to ensure adequate extraction fluid flow through the restrictor. In addition, very fine particles contained in samples can clog the frits, necessitating replacement.

7.5.3 Clean the extraction vessel after each sample. The cleaning procedure depends on the type of sample. After removing the bulk of the extracted sample from the extraction vessel, the vessel should be scrubbed with an ionic detergent, water, and a bottle brush. After extraction of tarry materials, use solvent rinses or an ultrasonic bath to clean the extraction vessel.

7.5.4 For samples known to contain elemental sulfur, use copper filings to remove the dissolved sulfur from the fluid. The copper filings (1 to 2 g per sample) can be packed in a separate extraction vessel connected to the outlet end of the sample extraction vessel, or they can be mixed with the sample, and a plug of copper filings can be loaded in the extraction

vessel with the sample such that any sulfur extracted by the carbon dioxide can be removed before the stream of carbon dioxide containing the analytes reaches the restrictor.

7.5.5 The procedure to be followed in emptying the syringe pump depends upon the type of fluid being used. In the case of carbon dioxide, which is a gas at ambient temperature and pressure, it is only necessary to vent the gas to a fume hood by allowing it to expand across the purge valve. Follow manufacturer's instructions in emptying the syringe pump.

7.5.6 To change fluid supply cylinders, follow manufacturer's instructions.

7.5.7 Restrictor removal and installation -- Follow manufacturer's instructions. When using fused-silica restrictors, it may be necessary to replace the restrictor after each sample, especially when extracting samples contaminated with heavy oils.

8.0 QUALITY CONTROL

8.1 Reagent blanks or matrix-spiked samples must be subjected to the same analytical procedures (Section 7.4) as those used on actual samples.

8.2 Refer to Chapter One for specific Quality Control procedures and to Method 3500 for sample preparation quality control procedures.

8.3 All instrument operating conditions must be recorded.

9.0 METHOD PERFORMANCE

9.1 Refer to Method 8440 and 8015 for performance data.

9.2 Use standard reference materials to establish the performance of the method with contaminated samples.

10.0 REFERENCES

1. Lopez-Avila, V., N. S. Dodhiwala, J. Benedicto, and R. Young, "SFE/IR Method for the Determination of Petroleum Hydrocarbons in Soils and Sediments," EPA 600/X-92/046 (W.F. Beckert, Project Officer), US EPA, Environmental Monitoring Systems Laboratory, Las Vegas, April 1992.

2. Pyle, S. M., and M. M. Setty, "Supercritical Fluid Extraction of High-Sulfur Soils with Use of a Copper Scavenger," *Talanta*, 1991, 38 (10), 1125-1128.

11.0 SAFETY

11.1 When liquid carbon dioxide comes in contact with skin, it can cause burns because of its low temperature. Burns are especially severe when the carbon dioxide is modified with organic liquids.

11.2 The extraction fluid, which may contain a modifier, usually exhausts through an exhaust gas port on the rear of the panel of the extractor. This port must be connected to a chemical fume hood to prevent contamination of the laboratory atmosphere.

11.3 When liquid carbon dioxide is used for cryogenic cooling, typical coolant consumption is 5 L/min (as decompressed gas), which results in a carbon dioxide level of 900 ppm for a room of 4.5 m x 3.0 m x 2.5 m, assuming 10 air exchanges per hour. The NIOSH time-weighted average concentration is 9,000 ppm (American Conference of Governmental Industrial Hygienists, 1991-1992).

APPENDIX B

PROPOSED DRAFT PROTOCOL FOR DETERMINATION OF TOTAL RECOVERABLE PETROLEUM HYDROCARBONS BY INFRARED SPECTROPHOTOMETRY (METHOD 8440)

METHOD 8440

TOTAL RECOVERABLE PETROLEUM HYDROCARBONS BY INFRARED SPECTROPHOTOMETRY

1.0 SCOPE AND APPLICATION

1.1 Method 8440 (formerly Method 9073) is used for the measurement of total petroleum hydrocarbons (TPHs) extracted with supercritical carbon dioxide from sediment, soil, and sludge samples using Method 3560.

1.2 Method 8440 is not applicable to the measurement of gasoline.

1.3 Method 8440 can detect TPHs at concentrations of 10 $\mu\text{g/mL}$ in extracts. This translates to 10 mg/Kg in soils when a 3-g sample is extracted by SFE (assuming 100 percent extraction efficiency) and the final extract volume is 3 mL.

1.4 All organic solvents used in this method should be recovered and recycled.

2.0 SUMMARY OF METHOD

2.1 Soil samples are extracted with supercritical carbon dioxide using Method 3560. Interferences are removed with silica gel, either by shaking the extract with loose silica gel, or by passing it through a silica gel solid-phase extraction cartridge. After infrared (IR) analysis of the extract, TPHs are quantified by direct comparison with standards.

3.0 INTERFERENCES

3.1 The parameter being measured (TPHs) is defined within the context of this method. The measurement may be subject to interferences, and the results should be interpreted accordingly.

3.2 Determination of TPHs is a gross measure of mineral oils only, and does not include the biodegradable animal greases and vegetable oils captured in oil-and-grease measurements. These non-mineral-oil contaminants may cause positive interferences with the IR analysis if they are not completely removed by the silica gel cleanup.

4.0 APPARATUS AND MATERIALS

4.1 Infrared spectrophotometer -- Scanning or fixed wavelength, for measurement around 2950 cm^{-1} .

4.2 IR cells -- 10-mm, 50-mm, and 100-mm path length, made of sodium chloride or IR-grade glass.

4.3 Optional -- A vacuum manifold consisting of glass vacuum basin, collection rack and funnel, collection vials, replaceable stainless steel delivery tips, built-in vacuum bleed valve and gauge

is recommended for use when silica gel cartridges are used. The system is connected to a vacuum pump or water aspirator through a vacuum trap made from a 500-mL sidearm flask fitted with a one-hole stopper and glass tubing.

5.0 REAGENTS

5.1 Reagent-grade chemicals shall be used in all tests. Unless otherwise indicated, all reagents shall conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society, where such specifications are available. Other grades may be used, provided it is first ascertained that the reagents are of sufficiently high purity to permit their use without adversely affecting the accuracy of the determinations. For cleanup, use silica gel cartridges or loose silica gel.

5.2 Solvents -- Spectrophotometric grade, or equivalent.

5.2.1 Tetrachloroethylene, C_2Cl_4

5.3 Materials for the preparation of the reference oil mixture -- Spectrophotometric grade, or equivalent.

5.3.1 n-Hexadecane, $CH_3(CH_2)_{14}CH_3$

5.3.2 Isooctane, $(CH_3)_3CCH_2CH(CH_3)_2$

5.3.3 Chlorobenzene, C_6H_5Cl

5.4 Silica gel

5.4.1 Silica gel solid-phase extraction cartridges (40- μ m particles, 60-Å pores), 0.5 g, Supelco, J.T. Baker, or equivalent.

5.4.2 Silica gel, 60 to 200 mesh, Davidson Grade 950, or equivalent (deactivated with 1 to 2 percent water).

5.5 Calibration mixtures

5.5.1 The material of interest, if available, or the same type of petroleum fraction, if it is known and original sample material is unavailable, shall be used for the preparation of calibration standards. Reference oil is to be used only for unknowns. Whenever possible, a GC fingerprint should be run on unknowns to determine the petroleum fraction type.

5.5.2 Reference oil -- Pipet 15.0 mL n-hexadecane, 15.0 mL isooctane, and 10.0 mL chlorobenzene into a 50-mL glass-stoppered bottle. Maintain the integrity of the mixture by keeping the bottle stoppered, except when withdrawing aliquots. Refrigerate at 4°C when not in use.

5.5.3 Stock standard -- Pipet 0.5 mL calibration standard (Section 5.8.1 or 5.8.2) into a tared 100-mL volumetric flask and stopper immediately. Weigh and dilute to volume with tetrachloroethylene.

5.5.4 Working standards -- Pipet appropriate volumes of stock standard (Section 5.5.3) into 100-mL volumetric flasks according to the cell size to be used. Dilute to volume with tetrachloroethylene. Calculate the concentrations of the standards from the stock standard concentrations.

5.6 Calibration mixture for silica gel cleanup -- Prepare a stock solution of corn oil by placing about 1 mL (0.5 to 1 g) of corn oil into a tared 100-mL volumetric flask. Stopper the flask and weigh to the nearest milligram. Dilute to the mark with tetrachloroethylene, and shake the contents to effect dissolution. Prepare additional dilutions to cover the range of interest.

6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

6.1 Solid samples should be collected and stored as any other solid samples containing semivolatile analytes.

6.2 Samples should be analyzed with minimum delay, upon receipt in the laboratory, and must be kept refrigerated prior to analysis.

7.0 PROCEDURE

7.1 Prepare solid and sludge samples according to Method 3560.

7.2 Add 0.3 g of the loose silica gel to the 3-mL extract and shake the mixture for 5 minutes, or pass the extract through a 0.5-g silica gel solid-phase extraction cartridge (conditioned with 5 mL tetrachloroethylene). When working with loose silica gel, filter the extract through a plug of precleaned silanized glass wool in a disposable glass pipet.

7.3 After the silica-gel cleanup, fill a clean IR cell with the solution and determine the absorbance of the extract. If the absorbance exceeds the linear range of the IR spectrophotometer, prepare an appropriate dilution and reanalyze. The possibility that the absorptive capacity of the silica gel has been exceeded can be tested at this point by repeating the cleanup and determinative steps.

7.4 Select appropriate working standard concentrations and cell path lengths according to the following ranges:

| Path length (mm) | Concentration range ($\mu\text{g/mL}$ of extract) | Volume (mL) |
|------------------|---|-------------|
| 10 | 5 to 500 | 3 |
| 50 | 1 to 100 | 15 |
| 100 | 0.5 to 50 | 30 |

Calibrate the instrument for the appropriate cells using a series of working standards. It is not necessary to add silica gel to the standards. Determine absorbance directly for each solution at the absorbance maximum at about 2950 cm⁻¹. Prepare a calibration plot of absorbance versus concentration of petroleum hydrocarbons in the working standards.

7.5 Determine the concentration of petroleum hydrocarbons in the extract by comparing the response against the calibration plot.

7.6 Calculate the concentration of TPHs in the sample using the formula:

$$\text{Concentration (mg/Kg)} = \frac{R \times D \times V}{W}$$

where:

R = mg/mL of TPHs as determined from the calibration plot
V = volume of extract in milliliters
D = extract dilution factor, if used
W = weight of solid sample in kilograms.

7.7 Recover the tetrachloroethylene used in this method by distillation or other appropriate technique.

8.0 QUALITY CONTROL

8.1 Reagent blanks or matrix-spiked samples must be subjected to the same analytical procedures as those used with actual samples.

8.2 Refer to Chapter One for specific Quality Control procedures and to Method 3500 for sample preparation procedures.

9.0 PRECISION AND ACCURACY

9.1 Table 1 presents a comparison of certified values and the values obtained using Methods 3560 and 8440. Data are shown for both Freon-113 and tetrachloroethylene, since both solvents were found to be acceptable as collection solvent. However, only tetrachloroethylene is recommended as collection solvent for TPHs in Method 3560.

9.2 Tables 2 and 3 present accuracy and precision data from the single-laboratory evaluation and the interlaboratory evaluation of Methods 3560 and 8440, respectively. These data were obtained by extracting samples at 340 atm/80°C/60 min dynamic (Table 2) or at 340 atm/80°C/30 min dynamic (Table 3).

10. REFERENCES

1. Rohrbough, W. G.; et al. Reagent Chemicals, American Chemical Society Specifications, 7th ed.; American Chemical Society, Washington, DC, 1986.
2. Methods for Chemical Analysis of Water and Wastes; U.S. Environmental Protection Agency. Office of Research and Development, Environmental Monitoring and Support Laboratory. ORD Publication Offices of Center for Environmental Research Information, Cincinnati, OH, 1983; EPA-600/4-79-020.
3. Lopez-Avila, V., N. S. Dodhiwala, J. Benedicto, and R. Young, "SFE/IR Method for the Determination of Petroleum Hydrocarbons in Soils and Sediments," EPA 600/X-92/046, U. S. EPA, Environmental Monitoring Systems Laboratory, Las Vegas, NV, April 1992.
4. Lopez-Avila, V., R. Young, and R. Kim, "Interlaboratory Evaluation of an Off-Line SFE/IR Method for the Determination of Petroleum Hydrocarbons in Solid Matrices," EPA 600/X-93/XXX, U.S. EPA, Environmental Monitoring Systems Laboratory, Las Vegas, NV, January 1993.

**TABLE 1. CERTIFIED AND SPIKE VALUES COMPARED WITH RESULTS
OBTAINED BY METHODS 3560/8440**

| Reference Material | Spike level or certified level (mg/Kg) | Methods 3560/8440 (mg/Kg) |
|---|--|---------------------------------|
| Environmental Resource Assoc. TPH-1 soil (Lot 91012) | 1,830 | 1,920 \pm 126 ^a |
| Environmental Resource Assoc. TPH-2 soil (Lot 91012) | 2,230 | 2,150 \pm 380 ^a |
| Clay soil spiked with kerosene | 100 | 86.0; 93.0 ^b |
| Clay soil spiked with light gas oil | 100 | 84.0; 98.0 ^b |
| Clay soil spiked with heavy gas oil | 100 | 103; 108 ^b |
| Clay soil spiked with medium neutral oil | 100 | 125 \pm 19.4 ^c |
| Clay soil spiked with heavy neutral oil | 100 | 126 \pm 15.8 ^d |
| Clay soil spiked with heavy lube oil | 100 | 118; 155 ^b |
| Environmental Resource Assoc. TPH-1 soil (Lot 91017) | 614 | 562; 447 ^b |
| Environmental Resource Assoc. TPH-2 soil (Lot 91017) | 2,050 | 1,780; 1,780 ^b |
| SRS103-100 soil | 32,600 | 29,100 \pm 1,930 ^e |

^a Three 60-min extractions. The extracted materials were collected in Freon-113; the concentrations were determined against the reference oil standard.

^b Duplicate 30-min extractions. The extracted materials were collected in tetrachloroethylene; the concentrations were determined against standards made from the spiking materials, except TPH-1 and TPH-2 soils where reference oil was used to determine concentrations.

^c Six 30-min extractions. The extracted materials were collected in tetrachloroethylene; the concentrations were determined against standards made from the spiking material.

^d Four 30-min extractions. The extracted materials were collected in tetrachloroethylene; the concentrations were determined against standards made from the spiking material.

^e Three 30-min extractions. The extracted materials were collected in tetrachloroethylene; the concentrations were determined against the reference oil standard.

TABLE 2. SINGLE-LABORATORY METHOD ACCURACY AND PRECISION FOR METHODS 3560/8440 FOR SELECTED MATRICES

| Matrix | Certified or spike value (mg/Kg) | Spike Material | Method accuracy (% recovery) | Method precision (% RSD) |
|---------------------------------------|----------------------------------|----------------|------------------------------|--------------------------|
| Clay soil ^a | 2,500 | Motor oil | 104 | 8.5 |
| ERA TPH-1 ^a (Lot 91016) | 2,350 | Vacuum oil | 80.3 | 19.7 |
| ERA TPH-2 ^a (Lot 91016) | 1,450 | Vacuum oil | 88.6 | 19.6 |
| SRS103-100 ^b | 32,600 | c | 94.2 | 4.0 |

^a Eight determinations were made using two different supercritical fluid extraction systems. The extracted materials were collected in Freon-113.

^b Ten determinations were made using three different supercritical fluid extraction systems. The extracted materials were collected in Freon-113.

^c This is a standard reference soil certified for polynuclear aromatic hydrocarbons. No spike was added.

TABLE 3. INTERLABORATORY METHOD ACCURACY AND PRECISION FOR METHOD 3560/8440 FOR SELECTED MATRICES

| Compound name | True concentration (mg/Kg) | Mean concentration ^a (mg/Kg) | s _r ^b (mg/Kg) | s _R ^c (mg/Kg) | Percent RSD _r ^d | Percent RSD _R ^e | Percent mean recovery | Number of laboratories |
|------------------------------------|----------------------------------|---|--|--|--|--|-----------------------------|---------------------------|
| Perkin-Elmer FTIR | | | | | | | | |
| ERA TPH-1 soil | 614 | 654 | 111 | 297 | 17.0 | 45.4 | 107 | 14 |
| ERA TPH-2 soil | 2,050 | 1,850 | 213 | 321 | 11.5 | 17.3 | 90.2 | 13 |
| SRS103-100 soil | 32,600 | 26,820 | 4,320 | 9,720 | 16.1 | 36.2 | 82.3 | 14 |
| Clay soil spiked with motor oil | 10,000 | 7,790 | 1,000 | 1,660 | 12.9 | 21.3 | 77.9 | 13 |
| Buck-Scientific IR | | | | | | | | |
| ERA TPH-1 soil | 614 | 618 | 113 | 296 | 18.2 | 47.9 | 101 | 14 |
| ERA TPH-2 soil | 2,050 | 1,670 | 194 | 278 | 11.7 | 16.7 | 81.5 | 13 |
| SRS103-100 | 32,600 | 24,750 | 3,740 | 8,650 | 15.1 | 35.0 | 75.9 | 14 |
| Clay soil spiked with motor oil | 10,000 | 8,180 | 910 | 1,500 | 11.1 | 18.3 | 81.8 | 13 |

^a The number of replicates per laboratory was three.

^b s_r - repeatability standard deviation.

^c s_R - reproducibility standard deviation.

^d RSD_r - repeatability relative standard deviation.

^e RSD_R - reproducibility relative standard deviation.

8440 - 8

Revision WG 0
December 1992

APPENDIX C

LIST OF INSTRUCTIONS FOR COLLABORATORS

1 ~ 2 ~ 3 ~
4 ~
5 ~
6 ~

Dear 1 ~ 3 ~:

Thank you for participating in the EPA's interlaboratory study on extracting total petroleum hydrocarbons (TPHs) from environmental samples using supercritical fluid extraction (SFE). Its purpose is twofold: first, to verify SFE as a viable alternative to the currently approved extraction techniques for TPHs, and second, to establish a feasible non-CFC solvent for collecting the material that is extracted by the supercritical carbon dioxide.

Your package should contain the following:

1. One set of instructions.
2. Forms to record SFE conditions and observations.
3. Nine vials containing soil samples labeled Vials No. 1 through 9. Vials No. 1 through 3 contain at least 10 g soil each. Vials No. 4 through 9 contain 3.00 g soil each. Refrigerate all samples upon receipt.
4. Twenty empty vials with Teflon-lined screw caps for sending the extracts to MRI-California Operations.
5. An additional vial is included (labeled MRI# 000074) that contains anhydrous sodium sulfate (Mallinckrodt AR); anhydrous sodium sulfate is recommended for use with the sample contained in Vial No. 9.

If anything is missing or damaged, please contact me immediately for replacement. We made arrangements with Aldrich to send you tetrachloroethylene (HPLC-grade, Lot No. 06626TX) immediately. If the lot number differs from this number, do not use, but notify me so that the correct tetrachloroethylene can be sent to you. SFE-grade carbon dioxide from Scott Specialty Gases will also be sent to you. Please do not use other grades of carbon dioxide for this study.

After you have collected the 15 extracts, please send them to MRI-California Operations by March 27, 1992 for analysis by IR.

Mr. Young
Page 2

March 5, 1992

For questions or comments, do not hesitate to write, call, or fax a message to:

Richard Young
MRI-California Operations
625-B Clyde Avenue
Mountain View, CA 94043
Phone: (415) 694-7700
Fax: (415) 694-7983

For your information, there are 15 laboratories participating in the study. Our plans are to complete all IR analyses and the statistical analyses by the end of April 1992, and to present the results of the study at the 8th Annual Waste Testing and Quality Assurance Symposium in Washington, DC, in July 1992. Also, when the study is completed, we will finalize the method protocols (Method 3560 for the SFE and Method 8440 for the IR determination of TPHs) and send you a copy as soon as they are approved by EPA.

Sincerely,

MIDWEST RESEARCH INSTITUTE —
CALIFORNIA OPERATIONS

Richard Young
Senior Chemist

Enclosures

INSTRUCTIONS FOR THE EPA INTERLABORATORY STUDY ON SFE/IR METHOD FOR TPHs

Your Laboratory Code is: 15

REAGENTS

1. Tetrachloroethylene Aldrich Lot No. 06626TX. Please do not use any other grade of tetrachloroethylene for this study.
2. Sodium sulfate (anhydrous) known to be free of interferences in the C-H stretch band range 3200 to 2700 cm^{-1}
3. SFE-grade carbon dioxide from Scott Specialty Gases. Please do not use any other grade of carbon dioxide for this study.

APPARATUS

1. Supercritical fluid extraction system. Follow manufacturer's instructions on how to operate the system. Make sure that the extraction vessels and fittings are rated to at least 340 atm. Also make sure that the temperature is maintained during the extraction at 80°C.
2. Since tetrachloroethylene freezes when carbon dioxide is expanded into it, the collection vessel must be immersed in water of room temperature (20 to 25°C).

SAFETY

1. A safety feature to prevent overpressurization is required on the extraction system. This feature should be designed to protect the laboratory personnel and the instrument from possible injuries or damage resulting from equipment failure under high pressure.
2. Liquid carbon dioxide can cause "burns" because of its low temperature (-70°C).
3. The material safety data sheets (MSDS) for tetrachloroethylene and sodium sulfate should be reviewed before their use.

INTERFERENCES

1. The analyst must demonstrate through the analysis of system blanks that the SFE system is free from interferences. To do this, perform a simulated extraction, using an empty extraction vessel and a known amount of carbon dioxide under the same conditions as those that normally will be used for sample extraction, and analyze the extract by IR. The blank should have < 10 ppm TPHs (based on using approximately 30 mL of carbon dioxide and collecting the extract in 3 mL tetrachloroethylene). If you cannot perform the IR analysis of the system blank, send the extract for the system blank, labeled with your laboratory code, the term "blank", and the form recording SFE conditions and observations, to MRI—California Operations for analysis.

2. The extraction vessel(s), frits, restrictor(s), and selector valve may retain solutes. Therefore, it is good practice to clean the extraction system after each extraction. Replacement of the restrictor may be necessary when system blanks indicate carryover. At least one system blank should be performed daily when the instrument is in use. Furthermore, system blanks should be performed after each extraction of a high-concentration sample. If system blanks continue to indicate contamination even after replacement of the extraction vessel and the restrictor, the multiport valve must be cleaned.
3. Be aware of any organic solvents used to clean any parts of the extraction system or glassware since these may introduce interferents.

SAMPLE PREPARATION

Remove the sample vial from the refrigerator and follow the specific directions for each vial.

1. Vials 1 and 2: perform triplicate extractions using 3-g portions of this material. Weigh 3.00 g sample each time.
2. Vial 3: perform triplicate extractions using 3-g portions of this material. Due to the concentration and nature of this sample (may plug restrictors), it is suggested that it be extracted last.
3. Vials 4 though 8: extract the entire contents of the vial weighing the vial before and after emptying to verify that the vial contained 3.00 g sample. Record any differences from this amount.
4. Vial 9: extract the entire contents, weighing the vial before and after emptying to verify that the vial contained 3.00 g sample. Record any differences from this amount. Add 3 g sodium sulfate as a drying agent.

For each extraction, transfer the sample to a clean extraction vessel and use two plugs of silanized glass wool to hold the sample in place and to fill the void volume. Attach the end fittings and install the extraction vessel in the oven. Always use clean frits for each extraction vessel.

EXTRACTION

1. Fill the collection vessel with 3 mL tetrachloroethylene Aldrich Lot No. 06626TX. Do not use any other grade of tetrachloroethylene. Since tetrachloroethylene freezes when carbon dioxide is expanded into it, the collection vessel must be immersed in a beaker with water of room temperature (20 to 25°C).
2. Set the SFE pressure to 340 atm and the temperature to 80°C. Extract for 30 min at a flowrate of 1 to 2 mL/min (as liquid carbon dioxide).
3. Record extraction conditions and observations on the forms provided. Identify each extraction by vial number. When triplicate extractions are performed on the same sample, add the suffix A,

March 5, 1992

B, or C. That is, the extracts should be labeled 1A, 1B, 1C, 2A, 2B, 2C, 3A, 3B, 3C for matrices in vials 1, 2, and 3, respectively.

4. Transfer the tetrachloroethylene extract to one of the vials provided. Do not readjust the solvent volume back to 3 mL. The vial should be sealed tightly with the Teflon-lined cap. Label the vial in a legible and permanent manner in the following way:

[15] - [extraction identification]

For example, if your lab code is 15 and you have just extracted the third sample from Vial 2, the label should read "15-2C."

5. Store the extract in a refrigerator at 4°C until it can be shipped to MRI-California Operations.

EXTRACT SHIPMENT

After all of the extractions have been performed, ship the extracts in the vials provided to:

Richard Young
MRI-California Operations
625-B Clyde Avenue
Mountain View, CA 94043

Ship the extracts packed in dry ice via overnight carrier. Include the forms with the recorded SFE conditions and observations. Please ship so that receipt occurs during weekdays (i.e., do not send out on Friday). Fax a message indicating when you are sending the extracts. Please ship extracts back by March 27, 1992.

Your participation is highly appreciated, and we look forward to receiving your extracts to make this a successful study. We also welcome your comments and suggestions concerning this study.

**FORM I: INSTRUMENT OPERATING CONDITIONS FOR SUPREX, ISCO,
DIONEX-LEE SCIENTIFIC, AND CCS EXTRACTION SYSTEMS^a**

Instrument: _____ Operator: _____ Date: _____
Supercritical fluid: _____ Supplier: _____ Grade: _____
Modifier: _____ Percent modifier: _____
Extraction vessel volume (mL): _____ Number: _____
Extraction vessel dimensions: _____ ID x _____ length (mm)
Extraction vessel position: Horizontal _____ Vertical _____
Restrictor: _____ ID x _____ OD x _____ length (mm)
Collection system: _____ Vial size: _____
Manufacturer: _____
Catalog number: _____
Solvent: _____ Volume: _____ Initial _____
Final _____

Operating conditions:

| | |
|--------------------------------|-------|
| Dynamic or static | _____ |
| Pressure (atm) | _____ |
| Temperature (°C) | _____ |
| Flowrate (mL/min) | _____ |
| Density (g/mL) | _____ |
| Time (min) | _____ |
| Restrictor temp. (°C) | _____ |
| Volume of CO ₂ (mL) | _____ |

^a If conditions vary for different samples, submit separate form.

**FORM II: INSTRUMENT OPERATING CONDITIONS FOR HP MODEL 7680A
EXTRACTOR**

SAMPLE HISTORY

Original entry:

Last update:

SAMPLE INFORMATION

Sample name:

Sample ID:

Operator:

Sample amount:

Comments:

EXTRACTION STEP

FLUID DELIVERY

Density:

g/mL

Pressure:

bars (convert bars to atm)

Flowrate:

mL/min

Extraction fluid:

CO₂

EXTRACTION CHAMBER

Chamber temperature:

°C

Equilibration time:

min

Extraction time:

min

Thimble size:

mL

Thimble volumes swept:

ANALYTE TRAP

Analyte:

Intermediate volatile

Trap material:

Nozzle temperature:

°C

Trap temperature:

°C

FORM II: INSTRUMENT OPERATING CONDITIONS FOR HP MODEL 7680A
EXTRACTOR (continued)

FRACTION OUTPUT

| Rinse substep | Solvent name | Volume (mL) | Rate (mL/min) | Nozzle temperature | Trap temperature | Vial number |
|------------------|-----------------|----------------|------------------|-----------------------|---------------------|----------------|
|------------------|-----------------|----------------|------------------|-----------------------|---------------------|----------------|