

June 11, 1997

ICR Laboratory Quality Control (QC) Database System Release 1.1

The attached diskettes contain the Release 1.1 update to the ICR Laboratory QC Database System. Release 1.1 is an incremental update to the application. It enhances Release 1.0 but does not replace it.

NOTES

The following describes changes made to the ICR Laboratory QC Database System subsequent to delivery of Release 1.0 of the application that are incorporated in Release 1.1. Unless otherwise specified, the changes made to the application should have no effect on data previously recorded in the application.

Sample Analytical Results

The following changes and corrections have been made to the Sample Analytical Results data entry portion of the ICR Laboratory QC Database System:

- The Internal Standard for samples analyzed under Method 551.1 for analyte groups Chloral Hydrate and Trihalomethane/Chloral Hydrate is no longer required in Release 1.1.
- When rejecting a Cryptosporidium results, Release 1.0 required the user to enter a numerical result. This was inconsistent with other analytes in the application. This has been corrected in Release 1.1.
- When comments were added for Giardia and Cryptosporidium in Release 1.0, the comments were not being properly being copied to the diskette sent to EPA. This has been corrected in Release 1.1.

INSTALLATION INSTRUCTIONS

Release 1.1 must be installed in the same directory as Release 1.0. Install Release 1.1 according to the following instructions:

- Start Windows.
- Insert the diskette labeled "ICR Laboratory Quality Control (QC) System Release 1.1, Disk 1 of 2" into a floppy disk drive on your PC.
- Choose Run from the File menu in Program Manager.

- Type a:setup in the Run dialog box. (Type b:setup if your diskette is in the B:\drive.)
- Press Enter or choose OK on the dialog box with your mouse.
- Follow the setup instructions that appear on the screen.

READ ME FILE

Additional text has been added to the "Read Me" file for the application. When Release 1.1 is opened for the first time, a window is displayed containing a button that allows the user to view the Read Me file.

INSTALLATION FILES

The following installation files are included on Disk 1 and Disk 2:

Disk 1

<u>File Name</u>	<u>Size</u>	<u>Date/Time</u>
_mssetup.exe	14KB	6/6/97 9:01 AM
acmsetup.ex_	119KB	6/6/97 9:01 AM
acmsetup.hl_	13KB	6/6/97 9:01 AM
admin.inf	1KB	6/6/97 9:03 AM
decomp.exe	39KB	6/6/97 9:01 AM
msacahbb.dl	73KB	6/6/97 9:01 AM
mscopydis.dl_	10KB	6/6/97 9:01 AM
mssetup.dll	187KB	6/6/97 9:01 AM
readme.tx_	3KB	6/6/97 9:01 AM
setup.exe	39KB	6/6/97 9:01 AM
setup.inf	2KB	6/6/97 9:01 AM
setup.ini	1KB	6/6/97 9:03 AM
setup.lst	1KB	6/6/97 9:01 AM
setup.sstf	3KB	6/6/97 9:03 AM

Disk 2

<u>File Name</u>	<u>Size</u>	<u>Date/Time</u>
lab_app.md_	965KB	6/6/97 9:03 AM