

United States **Environmental Protection** Agency

National Center for Environmental Assessment Research Triangle Park, NC 27711 EPA/600/R-05/102 October 2005

Research and Development

Guidance Manual

for the

All Ages Lead Model (AALM)

Draft Version 1.05

Prepared by

National Center for Environmental Assessment Office of Research and Development U. S. Environmental Protection Agency Research Triangle Park, N C 27711



NOTICE

The U.S. Environmental Protection Agency, through its National Center for Environmental Assessment at Research Triangle Park, produced this report. This document is a preliminary draft. It has not been formally released by the U.S. Environmental Protection Agency and should not be construed to represent Agency policy. It is circulated as a guide for reviewers of the model.

TABLE OF CONTENTS

AUTHORS AND CONTRIBUTORS	5
LIST OF FIGURES	6
SECTION 1: INTRODUCTION	8
SECTION 1: INTRODUCTION	8
1.1 INSTALLATION	9
1.2 BACKGROUND AND BASIS FOR PARAMETER VALUES	9
1.3 RUNNING THE MODEL	10
1.4 A QUICK REVIEW OF THE MODEL	12
2: BUILDING THE MODEL RUN	. 17
2.1. EDITING PARAMETERS 2.1.1. MEDIA EXPOSURES	19 19
2.2 EDITING EXPOSURE ACTIVITY PATTERNS	28
2.3 SUMMARY AND CONCLUSIONS	28
SECTION 3: SPECIAL EXPOSURE FEATURES	. 32
3.1 PICA EXPOSURE	32
3.2 DERMAL EXPOSURE	32
3.3 HISTORICAL LEAD EXPOSURES	33
3.4 MATERNAL EXPOSURE (Conceptual Design)	33
SECTION 4: ABSORPTION MODULE	. 34
4.1 MODELING LEAD ABSORPTION	34
4.2 ABSORPTION OF INHALED LEAD	36
4.3 ABSORPTION OF INGESTED LEAD	37
SECTION 5: BIOKINETICS MODULE	. 39
5.1 MODEL SETTINGS	. 39
5.1.1 PARAMETERS WITH EXPLICIT VALUES	. 39
5.2 BIOKINETIC PARAMETERS	42
5.2.1 Pb DECAY RATE.	43
5.2.2 CORTICAL BONE TURNOVER.	43
5.2.3. TRABECULAR BONE TURNOVER RATE.	44
5.2.3 TRANSFER FROM CORTICAL SURFACE TO BLOOD PLASMA.	46
5.2.4. TRANSFER FROM TRABECULAR SURFACE TO BLOOD PLASMA.	47
5.2.6 CORTICAL SURFACE TO VOLUME TRANSFER.	48
5.2.7 TRABECULAR SURFACE TO TRABECULAR VOLUME TRANSFER.	49
5.2.8 TOTAL TRANSFER FROM EXCHANGE BONE VOLUME.	. 50
5.2.9. IKANSFEK FKUM THE EACHANGE TO THE NON-EACHANGE VOLUME. 5.2.10. TRANSFER FROM LIVER1	. 50
5.2.10. IKANSFEK FKUM LIVEKI	
5,2,11 IKANSFEK FKUM KIDNE II.	

5.2.12 TRANSFER FROM BLADDER TO URINE	
5.2.13 TRANSFER FROM LIVER2	
5.2.14 TRANSFER FROM KIDNEY2	
5.2.15 TRANSFER FROM FAST SOFT TISSUE.	54
5.2.16 TRANSFER FROM INTERMEDIATE SOFT TISSUE	55
5.2.17 TRANSFER FROM SLOW SOFT TISSUE	
5.2.18 TRANSFER RATES FROM BRAIN	57
5.2.19 DEPOSITION FRACTION IN URINE.	58
5.2.20 DEPOSITION FRACTION IN BONE.	59
5.2.21 TRABECULAR BONE DEPOSITION FRACTION	60
5.2.22 DEPOSITION FRACTION IN COMPARTMENTS WITHOUT AGE VARIABI	IL 61
5.2.23 DEPOSITION FRACTION FROM DIFFUSIBLE PLASMA TO FAST SOFT	
5.2.24 DEPOSITION FRACTION FROM DIFFUSIBLE PLASMA TO INTERMEDIA	TE. 63
5.2.25 DEPOSITION FRACTION FROM DIFFUSABLE PLASMA TO SLOW SOFT.	64
5.2.26 DEPOSITION FRACTION FROM DIFFUSIBLE PLASMA TO BRAIN	65
5.2.27 COMPONENTS OF THE CIRCULATORY SYSTEM	
5.2.28 TRANSFER RATE FROM RBC TO PLASMA	
5.4.29 AMOUNT OF BLOOD	67
5.2.30 CHELATION FACTORS 1,2,and 3	
SECTION 6. CONTROLLING THE MODEL OUTPUT	69
6.1. A SIMPLE CONFIDENCE TEST	69
6.2. BIOKINETIC OUTPUTS	
6.2.1 TABLES	71
6.1.2 PLOTTING	
6.1.3 EDITING THE CHART	72
SECTION 7.	

AUTHORS AND CONTRIBUTORS

Robert Elias, U.S. Environmental Protection Agency Brian Gulson, Macquire University, Sidney Australia Gary Diamond, Syracuse Research Corporation Ramon Olivera, Lockheed-Martin Information Technology Goeffrey Nonato, Lockheed Martin Information Technology

LIST OF FIGURES

Figure 1.	Opening screen with NCEA Logo and disclaimer	10
Figure 2:	Window with prompt to create a new study or open an existing study.	11
Figure 3.	Drop Down window that allows the Modeler to customize three features of the program setup	11
Figure 4:	Primary Window - Modeler defines the individual by sex, life span, and time frame	12
Figure 5:	Model component selection window. Click on exposure to run this portion of the model	13
Figure 6:	Model status window. Select options to use historical data	13
Figure 7:	Retrieve model output file.	13
Figure 8:	Output data file in spreadsheet format.	14
Figure 9:	Primary Window with Male gender, End Date computation, and one day time step selected	15
Figure 10.	Drop-down window for resetting the model date.	16
Figure 11.	Bar along top of screen with Tools option.	18
Figure 12.	Tool Selection Menu.	18
Figure 13.	Drop down window to select the output parameters to be displayed.	18
Figure 14:	Primary window with full specifications for model run, ready for entering exposure parameters	19
Figure 15.	Save Current Study	19
Figure 16.	Select Exposure Parameters to review and modify the exposure portion of the model	20
Figure 17:	Exposure options group with Air component ready for review and modification	21
Figure 18:	Options for Dietary exposure.	22
Figure 19.	Dietary exposure with garden fruits and vegetables.	23
Figure 20:	Exposure options for dust exposure	24
Figure 21:	Drop Down window for the Dust component of Media Exposure.	25
Figure 22:	Window for the Dust component of Media Exposure	
Figure 23:	Window for the dust component of Media Exposure	26
Figure 24:	Window for the Water component of Media Exposure with the default option of Fully Flushed Water.	27
Figure 25:	Window for the Water component of Media Exposure with the options for Other Water Sources.	27
Figure 26:	Dron down window for the Activity Patterns component of Media Exposure	28
Figure 27:	Following set of the Exposure Scenario, the model run is initiated by checking the Exposure	29
Figure 28.	Two files are available for review. * mod and * mpa.	30
Figure 29.	Line by line output for the exposure module. * mod	30
Figure 30.	Total exposure by age for each exposure category (* mpa)	30
Figure 31.	Exposure Model options with options to change internal exposure and growth parameters	31
Figure 32:	Window for the Pica component of Media Exposure	32
Figure 33:	Window for the Dermal component of Media Exposure	33
Figure 34	Schematic diagram of the Leggett Model	35
Figure 35	Dron down window for reviewing and revising	36
Figure 36	Window for editing Absorption Parameters	36
Figure 37:	Drop down window for gastrointestinal absorption showing default values	37
Figure 38:	Plot of the default values for the gastrointestinal absorption fraction	37
Figure 39:	Default values for the gastrointestinal rate of movement function	38
Figure 40:	Plot of the default values for the rate of astrointestinal movement function	38
Figure 41.	Schematic Diagram of O'Flaherty Model	
Figure 42:	Window for editing the model settings	41
Figure 43	Dron down window for editing the default Delta Sten lengths	42
Figure 44	Window for editing hiskinetic parameters	42
Figure 45	Schematic Diagram of Leggett Model	43
Figure 46	Values for cortical hone turnover rate which decreases with age	15 44
Figure 47.	Plot of default values for cortical hone turnover rate by app	
Figure 48	Age-specific values for Surface to Volume in Trabecular Rone	
Figure 40.	Plot of Trabecular Surface to Volume	
Figure 50.	Values for trabecular hone turnover rate by age	+5 ⊿5
Figure 51.	Plot of default values for trabecular hone turnover rate by age	<i>1</i> 5 46
1 15m C 21.	I tot of activity rained for inaccount cone innover raile by agenininininininininininininininininini	0

Figure 52:	Dropdown window for editing the parameter for transfer from cortical surface to blood.	46
Figure 53:	Plot of the default values for the parameter for transfer from cortical surface to blood	47
Figure 54:	Drop down window for editing the parameter for the transfer rate for trabecular surface to blood	47
Figure 55:	Plot of the default values for the transfer rate for trabecular surface to blood	48
Figure 56.	Window for editing the value for transfer from cortical surface to volume	48
Figure 57:	Plot of the default values for transfer from cortical surface to volume	49
Figure 58:	Window for editing the values for transfer from trabecular surface to volume	49
Figure 59:	Plot of the default values for transfer from trabecular surface to volume	50
Figure 60:	Window for editing values for total transfer from the exchange bone volume.	50
Figure 61:	Window for editing values for total transfer from the non-exchange bone volume	51
Figure 62:	Window for editing the values for the transfer of lead from the bladder to urine.	52
Figure 63:	Plot of the default values for the transfer of lead from the bladder to urine	52
Figure 64.	Age-related coefficients for transfer of lead from Liver2	53
Figure 65.	Plot of the transfer of lead from the liver 2 compartment	53
Figure 66:	Drop down window for editing the default values for the transfer of lead from the kidney	54
Figure 67:	Plot of the default values for the transfer of lead from Kidney 2 (From Kidney to Bladder)	54
Figure 68:	Drop down window for editing the default values for the deposition fraction of lead in fast soft tissue.	. 55
Figure 69:	Plot of the default values for the deposition fraction of lead in fast soft tissue	55
Figure 70:	Drop down window for editing the default values for the deposition in intermediate soft tissue	56
Figure 71:	Plot of the default values for the deposition fraction of lead in fast soft tissue	56
Figure 72:	Drop down window for editing the default values for the deposition fraction in slow soft tissue	57
Figure 73:	Plot of the default values for the deposition fraction of lead in slow soft tissue	57
Figure 74:	Drop down window for editing the default values for the deposition fraction of lead in the brain	58
Figure 75:	Plot of the default values for the deposition fraction of lead in the brain	58
Figure 76.	Deposition fraction in Bone.	59
Figure 77.	Plot of Deposition Fraction in Bone.	59
Figure 78.	Bottom of biokinetic Parameters Window	60
Figure 79.	Age-dependent Variables for Trabecular Bone Deposition Fraction.	60
Figure 80.	Trabecular bone deposition fraction	61
Figure 81:	Plot of the default values for the transfer of lead from the liver	61
Figure 82.	Age-dependent Variables for Diffusible Plasma to Fast Soft Tissue	62
Figure 83.	Plot of the default values for the Diffusible Plasma to Fast Soft Tissue parameter	63
Figure 84.	Age-dependent Variables for Diffusible Plasma to Intermediate Soft Tissue	63
Figure 85.	Plot of the default values for the diffusible Plasma to Intermediate Soft Tissue parameter	64
Figure 86.	Age-dependent Variables for Diffusible Plasma to Slow Soft Tissue	64
Figure 87.	Plot of the default values for the Slow Soft Tissue Parameter.	65
Figure 88.	Age-dependent Variables for the Deposition Fraction from Diffusible Plasma to Brain	65
Figure 89.	Plot of the default values for the Deposition Fraction from Diffusible Plasma to Brain parameter	66
Figure 90.	Age-dependent variable for the Transfer Rate from RBC.	67
Figure 91.	Plot of the default values for the Transfer Rate for RBC parameter	67
Figure 92.	Age-dependent variable for total volume of blood (deciliters	68
Figure 93.	Plot of the default values for the Amount of Blood parameter.	68
Figure 94.	Output with all default settings	69
Figure 95.	Total Blood lead with elevated dust lead exposure.	.70
Figure 96.	Using the Built in Graphic Feature.	72

SECTION 1: INTRODUCTION

The purpose of the All Ages Lead Model (AALM) is to make available to the risk assessor a sufficiently complex, multi-compartment biokinetic model that predicts with reasonable accuracy the tissue concentrations of lead in humans and informs the modeler of the resulting changes in tissue concentrations that might occur over the lifetime of the modeled individual. The original model developed for this purpose, the Integrated Exposure Uptake Biokinetic Model for Lead in Children (IEUBK) met some of these objectives in that it is a multi-compartmental model that predicts blood concentrations of lead in humans, but the IEUIBK model makes this prediction only through age six and provides no information on the tissue concentrations (White et al, 1998). The AALM accomplishes this objective by providing an Exposure and Absorption interface to a biokinetic model.

This interface guides the modeler toward a thorough, accurate description of the total daily exposure to lead, calculating the amounts of lead absorbed. The model then offers a choice of two biokinetic models to determine how the lead is distributed to the key body tissues¹. These two models are referred to as the Leggett Model (Leggett, 1993) and the O'Flaherty Model O'Flaherty, 1998). One advantage of the Leggett Model is that it characterizes the distribution of lead based on the function of the human circulatory system, including the extracellular fluids. Likewise, the Leggett model also recognizes the chemical relationship between Pb and Ca. especially in describing the complexity of the bone compartments. There are some deficiencies in the Leggett model that can perhaps be overcome by the use of the O'Flaherty model, and we expect to make this feature available in the near future.

This document takes you step-by-step through the process of getting the AALM up and running on your computer. As with all software packages, expertise is achieved by trial and error, and by aggressively investigating all model features. You will find that by starting from the simplest case and moving systematically to the more complex applications, you will soon develop expertise that matches your interests.

This Guide takes into account comments and suggestions from many competent reviewers, evaluators and experienced modelers. It is not intended as a primary source of information regarding the source code for the AALM model, which is documented elsewhere. That document identifies the inherent structure of the three components of the AALM: lead exposure (intake) model, lead absorption (uptake) model, and biokinetic distribution model. It documents how the software system was designed to implement these models and describes each of the data elements used by the system, the functional flow of data through the system and the interfaces between the components.

Model development is not an easy effort. In addition to the science behind the process, there are also concerns about the development process itself. During the AALM development, which occurred over a period of approximately eight years, several symposia and workshops were held, including one on Model Validation (Elias et al, 1998). At that workshop, two papers have

¹ The O'Flaherty model is not fully implemented at this time. We expect to complete this task in the next few weeks. Meanwhile, if the O'Flaherty option is selected, the model will switch to the Leggett model and notify the modeler thereof.

provided significant guidance to this project: Oreskes: Evaluation (Not Validation) of Quantitative Models; and Mickle: Structure, Use, and Validation of the IEUBK Model.

The development of a large modeling system such as this is an ongoing process. The current version will, hopefully, meet the needs of most modelers. But in the process of developing this model, the planners recognized several features that could be incorporated with minimal effort, given the successful completion of the basic model. Construction of these additional features has been initiated and is discussed in Chapter 4 of this document. Comments are most welcome.

1.1 INSTALLATION

AALM requires Windows 95 or later (Windows 98, or Windows NT, Windows 2000, Windows XP). If you already have a previous review version of the AALM installed, the new version includes a user option to automatically uninstall the previous version. AALM is distributed as a bundled package that, when activated, will place all required files in a single directory with a default designation of C:\Program Files\AALM. The AALM program is installed with InstallShield using the Add/Remove feature of the Windows Control Panel. Thus, a new subdirectory called will appear in this directory after installation, and the All Ages Lead Model icon should appear on your opening window. Unless you designate otherwise, your project files will be stored in a subdirectory of this directory. One additional file, DEMO1.xal has also been installed to assist you in becoming familiar with the AALM. You will have an opportunity to access this file as you work your way through this manual.

1.2 BACKGROUND AND BASIS FOR PARAMETER VALUES

The purpose of the All Ages Lead Model is to mathematically provide an exposure, absorption, and biokinetic infrastructure that allocates, by simulation, the simultaneous distribution of absorbed lead in several major body components and thereby predict at any point in time the concentration of lead in these components.

The All Ages Lead Model has about 190 parameters with assigned values that can be revised by the modeler to meet the description of a modeled individual designed by the modeler. All values for these parameters have default settings that permit the model to run to a conclusion, reaching an endpoint whereby the predicted outcome can be determined that is consistent with the modeled input. This means that the model will run in the default mode, and will produce reasonable results consistent with the default inputs. The modeler, therefore, has the option to change one or more of the input parameters and rerun the model to compare the new results with the default output, or to create and save another "default" model parameters consistent with the needs and requirements of his/her investigation. A companion document, the AALM Parameters and Equation Dictionary, describes these model components in greater detail.

The AALM default values for the biokinetic parameters are based on published experimental values or ranges of values used for the Leggett and O'Flaherty models. Many of these default values have no meaning or significance other than as a place holder for a value that may be assigned by the modeler. With few exceptions, there are no documented values that can be validated by reliable literature reports. Consequently, the modeler must always be aware of the need to review the results rigorously, and establish internal constraints that adequately bracket the desired objective.

1.3 RUNNING THE MODEL

After installation, click on the All Ages Lead Model icon to open the program. The NCEA logo (Figure 1) should appear with text disclaiming responsibility prior to official release. Click OK. If you highlight box to "close the screen automatically", the next time you run the model, this Logo screen will be displayed for about five seconds, then move on to the window in Figure 2.



Figure 1. Opening screen with NCEA Logo and disclaimer

This window gives you the option to create a new study or open an existing study. The default is to create a new study, but you may want to start your orientation with a demonstration study that will to provide you with a quick review of the model.

Select Action - AALM
What would you like to do?
Action
C Create a new study
 Open an existing study
Don't ask me this again You can create or open a study at any time from the File menu or from the Tool Bar. OK Cancel

Figure 2: Window with prompt to create a new study or open an existing study. Your new file will be saved in a subdirectory of the AALM Program directory. You can change this default directory by clicking on Tools at the top of the screen, then Options. Type the new directory and any subdirectories in the window labeled Data.

Options 🛛 🔀
Startup Preferences Close the opening screen automatically Display the Select Action dialog box after the opening screen
Default Directories Data: C:\PROGRAM FILES\AALM
Definition Files: C:\PROGRAM FILES\AALM\DefFiles
OK Cancel

Figure 3. Drop Down window that allows the Modeler to customize three features of the program setup.

The features of the Startup Preferences window are: The automatic closing of the opening screen after a 5 second delay, the choice to display or not to display the Select Action dialog box after the opening screen, and the opportunity to change the directory where the user-created output data files are stored. This last choice, which is, by default, the AALM subdirectory of the C:\Program directory, is especially important if the Modeler is accustomed to backing up only the files located in the C:\Documents and Settings directory, in which case, the AALM output files would not be saved. This Window may be accessed from the upper bar below the Screen Title by selecting the "Tools" option.

1.4 A QUICK REVIEW OF THE MODEL

If this is your first time using the model, perhaps it would be a good idea to begin with a demonstration file that will provide the required settings you will need to get started. Select 'Open an existing study' from the Select Action window, then choose **Demo1.xal** to get started.

Follow these steps to run through the model. (You can return to this place in the Guidance Manual to continue your review of the model features):

If you installed the program according to the instructions in Section 1.1, you should be able to select **Demo1.xal** from the File Open window. Click on the Open button to get the next screen, which is the hub of the model and is referred to in this manual as the **Primary window**. From here you can specify the input parameters for the model using the options on the screen. First note that this is a typical Windows® screen with drop down menus for many Windows® functions at the top. The model will use the default screen color format you have selected for your computer.

Cender Gender • Male C Female	Comp End [utation at Date	7/29/2005 💌	Simulation time step	Parameters Media Exposure	_ [] [
Age Range Category	-Age Rar	nge (days)	Age Range (yrs)	Date Ranges	Scenarios	Bup
Frenatal	-155	5 to 0	< 0		Maternal	
🔽 Postnatal	1	to 183	0.0 to 0.5	08/22/1915 - 02/20/1916	Residential	Insert Event
Toddler	184	to 1460	0.5 to 4.0	02/21/1916 - 08/20/1919	Residential	Defaults
Preschool	1461	to 2190	4.0 to 6.0	08/21/1919 - 08/19/1921	Residential, Daycare	<u>C</u> lear All
Vouth	2191	to 4745	6.0 to 13.0	08/20/1921 - 08/17/1928	Residential, Recreational, School	
Adolescent	4746	to 6935	13.0 to 19.0	08/18/1928 - 08/16/1934	Residential, Recreational, School	
Voung Adult	6936	to 14600	19.0 to 40.0	08/17/1934 - 08/11/1955	Residential, Recreational, Occupational	
Middle Age	14601	to 23725	40.0 to 65.0	08/12/1955 - 08/04/1980	Residential, Recreational, Occupational	
🔽 Early Retirement	23726	to 29930	65.0 to 82.0	08/05/1980 - 07/31/1997	Residential, Recreational	Close
🔽 Late Retirement	29931	to 32850	82.0 to 90.0	08/01/1997 - 07/29/2005	Residential	Help

Figure 4: Primary Window - Modeler defines the individual by sex, life span, and time frame, and accesses the parameter definition windows.

These settings will run the model using the default settings for the full lifetime of the hypothetical person that you define by entering or changing the appropriate values in the exposure section of the model. This file can be saved and used in your own research and studies, and you can create as many other study files as you need. In the **Demo1** file, the individual is tracked through all age ranges. Press Run and you will begin this simulation leading to the next screen:

12. Model Run	$\overline{\mathbf{X}}$
🔽 Exposure	
🔽 Absorption	
Biokinetics	Edit Model Settings
Co	ontinue Cancel

Figure 5: Model component selection window. Click on exposure to run this portion of the model. Select all three to run the full model.

To run the full model, click the small boxes next to Exposure and Absorption (the Biokinetics box will automatically insert a $\sqrt{}$). Press **Continue** to get the next screen.

12a. Exposure Model Processing	X
Processing Study File	Start
	Cancel
1	0%
Show Exposure Results	🔲 Use Historical Diet Pb Concentration
	🔲 Use Historical Air Pb Concentration
Press [Start] button to run process, [Cancel] to stop.	🔲 Use Historical Dust Pb Concentrarion

Figure 6: Model status window. Select options to use historical data.

Click **Start**, and the model will run for a few seconds. Then you will get two output files which can be viewed (Figure 7); these will be saved as a part of this study.

12b. Select	×
Select the exposure model output file to view.	
DEM01.mod	
DEMOT.mpa	
) (iour Continue	

Figure 7: Retrieve model output file. The *.mod file has day by day record of values for the list of selected output options. This list may take a few seconds to pop up on the screen, as it may contain over 32000 records, depending on your choice of age ranges. Press **Continue** to save these files using the same filename and to pass onto the next screen. In a few seconds, you will get a screen that looks like this:

	Days	Blood	Cortical	Trabecular	Liver	
10	00.00	59.57320730	620.7635935	155.1908984	84.05649215	
20	00.00	74.99385985	1014.870227	253.7175567	134.6291828	
30	00.00	87.52161683	1068.310926	267.0777316	159.2715834	
40	00.00	92.66780516	1102.618125	275.7683800	185.7980257	
50	00.00	93.95207832	1137.706059	285.4432296	192.2514812	
60	00.00	96.95696595	1178.948697	297.0188847	196.2516904	
70	00.00	100.5488671	1218.005640	308.1954584	199.0910815	
80	00.00	104.3684032	1253.168382	318.4629236	201.2287977	
90	00.00	108.3233116	1284.685231	327.8326323	202.8755518	
10	00.00	112.4094359	1313.288044	336.4502432	204.1488684	
11	00.00	116.6512371	1339.745199	344.4785205	205.1256620	
12	200.00	121.0825655	1364.743795	352.0649610	205.8612626	
13	300.00	125.7409657	1388.873538	359.3356916	206.3972082	
14	400.00	130.6665373	1412.639076	366.3969987	206.7650919	
15	500.00	120.7866609	1415.891046	367.6363811	196.6072323	
16	600.00	117.6172172	1380.196595	357.8301706	182.1601065	
17	700.00	119.0561251	1360.012449	352.0897615	174.0447650	
18	300.00	122.0980072	1349.402495	348.8150850	168.7021013	
19	900.00	123.7529231	1348.450382	348.1877407	165.2048646	
20	00.00	125.1638764	1357.137548	350.5088979	163.6887077	
21	00.00	127.2397138	1372.120817	354.8038135	163.7381440	
22	200.00	126.0974126	1390.575201	360.2104014	163.8271536	
23	300.00	115.9016073	1375.607222	355.7690429	153.6127219	
104	00.00	445 4500005	4070 474500	054 0000004	4.40.0500404	

Figure 8: Output data file in spreadsheet format. All output parameters may be selected, and the file may be exported as a Lotus 123 file or Microsoft Excel file.

This table shows the output for the model run, giving the lead concentrations $(\mu g/g)$ in each of six compartments every one hundred days for the lifetime of the modeled individual. You can plot these now if you wish.

At this point you should understand that the main purpose of the model is to estimate the concentration of lead in selected body tissues based on the exposure of the modeled individual over the desired age range (the model runs from birth to the highest selected age range). You will see later how you can use the plotting capabilities of the model to view the data in graphic form, or export the data to a spreadsheet and develop your own custom made plots and perform further statistical analyses to meet your specific data evaluation requirements.

Now let's get back to the screen-by-screen description of the model. Returning to the Primary Window, observe the following options:

🗒 1. Untitled1					_ 🗆 🔀
Gender Male C Female	Computation at End Date	7/28/2005 💌	Simulation time step	Parameters Media Exposure	Edit
Age Range Category	Age Range (days)	Age Range (yrs)	Date Ranges	Scenarios	Bun
Prenatal	-155 to 0	< 0		Maternal	
Postnatal	0 to 0	0.0 to 0.5		Residential	Insert Event
Toddler	0 to 0	0.5 to 4.0		Residential	Defaults
Preschool	0 to 0	4.0 to 6.0		Residential, Daycare	<u>C</u> lear All
C Youth	0 to 0	6.0 to 13.0		Residential, Recreational, School	
Adolescent	0 to 0	13.0 to 19.0		Residential, Recreational, School	
🔲 Young Adult	0 to 0	19.0 to 40.0		Residential, Recreational, Occupational	
🔲 Middle Age	0 to 0	40.0 to 65.0		Residential, Recreational, Occupational	
Early Retirement	0 to 0	65.0 to 82.0		Residential, Recreational	Close
Late Retirement	0 to 0	82.0 to 90.0		Residential	<u>H</u> elp

Figure 9: Primary Window with *Male* gender, *End Date* computation, and *one day* time step selected.

You can:

- 1) Select Male or Female
- 2) Compute with an end date of today
- 3) Change the date from the default date, which is your computer's date
- 3) Select a simulation timestep of one day or one hour²
- 4) Run the model through one to nine age ranges³.

The Calendar date <u>cannot</u> be changed by editing the date in the date window, because the program does not read this date back in to the study file. The date <u>can</u> be changed by clicking on this window, thereby calling up a calendar which is capable of recording your selected date in the study file (Figure 10). Click on the up/down arrow at the right of the year to change the year. Click on the month right or left arrow to change the month, and then click on the day to set the day of the month.

Default options are preselected for all choices on this screen except the Age Range, and you need select only the highest Age Range that you want for your simulation. The rest will be populated for you. Take time to note that when you select an age range, the model displays the ages and the dates automatically. Thus, if you are interested in an exposure during a specific time period in this person's life, you can make an upper age selection that will include this time period in the person's life. You can also manually adjust the number of days in each age range.

² Our experience to date has not revealed any situations where the assessment would be enhanced by a timestep more frequent than one day.

³ The Prenatal option is not fully implemented but is included for demonstration purposes. (See Section 4)

Associated with each age range are one or more exposure scenarios that can be manipulated to define the daily exposure for one to four physical locations (e.g. Residential, School, Occupational, Recreational) during each of the life stages⁴.

When you are through with this exploration phase, exit the simulation and return to the main window. If you wish to reset the parameters to their default value you may do so on any of the Media Exposure windows.

📅 1. Untitled1		
Gender C Male 💿 Female	Computation at Simulation time step Parameters End Date 1/20/2006 Addition time step Media Exposure	<u>E</u> dit
Age Range Category Prenatal Postnatal Toddler Preschool Youth Adolescent	Image: A state of the stat	Ru <u>n</u> Insert Event Defaults <u>C</u> lear All
 ✓ Young Adult ✓ Middle Age ✓ Early Retirement ✓ Late Retirement 	6936 to 14600 19.0 to 40.0 02/08/1935 - 02/02/1956 Residential, Recreational, Occupational 14601 to 23725 40.0 to 65.0 02/03/1956 - 01/26/1981 Residential, Recreational, Occupational 23726 to 29930 65.0 to 82.0 01/27/1981 - 01/22/1998 Residential, Recreational 29931 to 32850 82.0 to 90.0 01/23/1998 - 01/20/2006 Residential	Close Help

Figure 10. Drop-down window for resetting the model date. Clicking on the date window will give the calendar window as shown. The years is selected by toggling up or down on the small bars to the right of the year. The month is selected by clicking on the left or right arrow on the left or right of the top of the screen. The day is selected by clicking on the desired date.

⁴ The two youngest age groups do not have the Occupational option, and therefore have only three location options.

2: BUILDING THE MODEL RUN

Building the model run consists of five steps:

- 1) define the individual
- 2) construct a full exposure scenario
- 3) modify the uptake and biokinetic parameters, if necessary
- 4) run the model to produce and save the output file
- 5) patch the output file into an Excel or Lotus spreadsheet from which you can perform further graphic and statistical analyses.

In every case throughout the model, default input values have been inserted when such a value is required for continuation. Accepting these values during the initial exploration of the model may enhance the learning curve and give you a feel for the range of values that will meet the needs of your study. With the default values that are provided, the model should run as installed with minimum modeler input. They can easily be changed and reset.

- a. Age Range: Select the maximum age category you would like for this model run. This will automatically populate all younger age categories. Note that these age categories are for the purpose of creating the exposure scenario based on reasonable activity patterns. Later, you will see age-related distributions for absorption and biokinetic parameters that are independent of these exposure and activity pattern age groups.
- b. Select the output you desire, with respect to the tissue concentrations. This is done by clicking on the Tools option on the bar at the top of the screen (Figure 11), then selecting "Set Output." See Figure 12. The Window is shown in Figure 13. Left click to place an 'X' in any box for which you want output, and remove the "X" from any box for which you don't want output by left clicking on the "X." Left click on OK when you are done selecting the output options.

For the purposes of graphic presentation, it is best to group your output selections so that they are approximately on the same scale and will display over about the same range. This can best be accomplished by trial and error.

c. Although there are several other options that could be brought into play, it might be best to run the model in default mode, just to get a feel for the type of output you can expect. Click the RUN button at the right of the box. You must save this study to a file. Click OK when asked "Save current study changes?" Clicking Cancel will stop the RUN mode at this point and wait for your further input. You can enter a new file name (the .xal extension will be added automatically) or click on a previous file from the list shown (Figure 14). Click SAVE.

😹 All Ages Lead Model - DEMO1.xal	_ 0 🛛
File Edit View Parameters Run Tools Window Help	

Figure 11. Bar along top of screen with Tools option. Clicking on Tools, then Set Output, will give the drop down window in Figure 13.

🗟 🛛 All Ages Lead Model - C:\Pi	rogram Files\AALM\Revie	w 22.xat	_ 7 🗙
File Edit View Parameters Run	Tools Window Help		
	View Results Set Output Plot Multiple MOB Options Calculator (Days/Year)		

Figure 12. Tool Selection Menu. Click on 'Set Output'.

	Output	Selected?	^
1	Sigma		
2	Plasma		
3	RBC		
4	Blood	\boxtimes	
5	Skeleton		
6	Cortical	$\overline{\boxtimes}$	
7	Trabecular	$\overline{\boxtimes}$	
8	Liver	$\overline{\boxtimes}$	=
9	Kidney		-
10	Soft		
11	Brain		
12	Lung		
13	Urine		
14	Feces		
15	Total Body2		
16	Total Excretion		
17	Bone Fraction		
18	Brain Fraction		
19	Liver Fraction		
20	Blood Fraction		
21	Renal Fraction		
22	Other Fraction		
23	Blood Concentratio	n 🗖	_
24	RBC Concentration		4

Figure 13. Drop down window to select the output parameters to be displayed and for which data are needed for export.

NOTE: The Insert Event option on the Opening window has been temporarily disabled.

💬 Demonstration.xa	d					
Gender © Male © Female	Comp	utation at Date 💌	7/29/2005 💌	Simulation time step	▼ Parameters	Edit
Age Range Category	-Age Rar	nge (days)	Age Range (yrs)	Date Ranges	Scenarios	Bun
Prenatal	-158	5 to 0	< 0		Maternal	
Postnatal	1	to 183	0.0 to 0.5	08/22/1915 - 02/20/1916	Residential	Insert Event
Toddler	184	to 1460	0.5 to 4.0	02/21/1916 - 08/20/1919	Residential	Defaults
Preschool	1461	to 2190	4.0 to 6.0	08/21/1919 - 08/19/1921	Residential, Daycare	<u>C</u> lear All
Vouth	2191	to 4745	6.0 to 13.0	08/20/1921 - 08/17/1928	Residential, Recreational, School	
Adolescent	4746	to 6935	13.0 to 19.0	08/18/1928 - 08/16/1934	Residential, Recreational, School	
Voung Adult	6936	to 14600	19.0 to 40.0	08/17/1934 - 08/11/1955	Residential, Recreational, Occupational	
Middle Age	14601	to 23725	40.0 to 65.0	08/12/1955 - 08/04/1980	Residential, Recreational, Occupational	
Early Retirement	23726	to 29930	65.0 to 82.0	08/05/1980 - 07/31/1997	Residential, Recreational	Close
✓ Late Retirement	29931	to 32850	82.0 to 90.0	08/01/1997 - 07/29/2005	Residential	Help
				1	J	

Figure 14: Primary window with full specifications for model run, ready for entering exposure parameters.

AALM		
2	Save curi	rent study changes?
	ОК	Cancel

Figure 15. Save Current Study

2.1. EDITING PARAMETERS

2.1.1. MEDIA EXPOSURES

There is considerable flexibility in the AALM to change the assigned values of the media exposure parameters. Each exposure category (AIR, DIET, DUST, WATER, PICA, DERMAL) has a separate window, and you can reset the values for these six windows. Air, Diet, Dust, and Water ingestion values are set to typical U.S. values and are discussed in this section. Pica and Dermal exposure values are, by default, set to zero exposure because they apply to special situations not fully documented in the literature. Pica and Dermal exposure are available for research purposes and are discussed further in Section 4.

Set the parameter window to "Media Exposure" and click the EDIT button. There are two options in Window 1a: "Exposure Parameters" and "Historical Exposure Settings." At this time, select "Exposure Parameters" to gain access to Window 1a1: Media Exposure (Figure 17). [Historical Exposure Settings are under development and may not execute correctly for this evaluation]. The Exposure Parameter window has six tabs, one each for Air, Diet, Dust, Water, Pica, and Dermal Exposure. In each of these windows, you can set the values for the key parameters for each of the age range options that you have selected, shown in Figure 17 as Postnatal to Late Retirement.

2.1.1.1 Air

The table for air exposure, shown on this figure, has eight rows of exposure data populated with two exposure parameters (amount of lead in air ($\mu g Pb/m^3$) and amount of air inhaled (m^3/day) for each exposure component: outdoor, residential, school, and occupational. In this way, air exposure is calculated as the product of the two parameters $\mu g/m^3$ and m^3/day , giving $\mu g/day$. In the case of air exposure in the three indoor environments, the amount of lead in air is, for convenience, expressed as a percentage of the outdoor lead concentration. The default values are set to 30% but may exceed 100% if the modeler chooses to do so.

The exposure module uses a smoothing feature that assigns the value you select to the middle of the age range, then interpolates to the adjacent points of the younger and older age ranges when assigning values, day by day, to the intervening time steps.

1a. Select Exposure Pa	arameter	To	X
Exposure Parameters	-	Edit	
	Г	_	
	L	Cance	

Figure 16. Select Exposure Parameters to review and modify the exposure portion of the model.

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retirem
Start Date	08/21/1915	02/20/1916	08/20/1919	08/19/1921	08/17/1928	08/16/1934	08/11/1955	08/04/1980	07/31/1997
Duration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
Outdoor Air Lead Conc. (ug/m3)	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Residential Air Pb Conc. (% of outdoor)	30	99	65	30	30	30	30	30	30
School Air Pb Conc. (% of outdoor)	30	88	50	30	30	30	30	30	30
Occupational Air Pb Conc. (% of outdoor)	30	30	30	30	30	30	30	30	30
Residential Ventilation Rates (m3/h)	0.37	0.37	0.47	0.5	0.53	0.55	0.55	0.55	0.48
School Ventilation Rates (m3/h)	0.37	0.37	0.47	0.5	0.53	0.55	0.55	0.55	0.48
Occupational Ventilation Rates (m3/h)	0	0	0	0	0	3.4	3.4	2.5	1.33
Recreational Ventilation Rates (m3/h)	0.65	0.65	1.74	1.9	2.4	2.6	2.6	2.5	1.33

Figure 17: Exposure options group with Air component ready for review and modification.

The default value of outdoor air is $0.1 \ \mu g \ Pb/m^3$ (Figure 17). Other options depend on the age range selected. The Clear button will clear all values in the "Air" window. The Clear All button will clear *all six exposure categories* and should not be used unless you plan to enter new values for all six windows. Likewise, the Default button will repopulate this window with the default values as shown, and the All to Default will do so for all six windows. The Historical Lead Concentration check window will be discussed in Section 4.

2.1.1.2 Diet

The data for lead concentrations in food are taken from the Market Basket Food Survey conducted by the Food and Drug Administration. In this survey, specific food products in over two hundred categories are purchased at supermarkets periodically throughout the United States and sent to a central laboratory where they are prepared in a kitchen environment according to the manufacturer's instructions. These are analyzed for many organic and inorganic components according to a strict analytical protocol. These data are reported quarterly and are available at the Food and Drug Administration website:

Market Basket Survey: <u>http://vm.cfsan.fda.gov/~lrd/pestadd.html</u> Total Diet Study: <u>http://www.cfsan.fda.gov/~comm/tds-hist.html#fca</u>

For this report, the data have been compiled, with some adjustments, and have been incorporated into the AALM for the periods 1982 to the present. Over the years, the categories of food have changed according to national food preferences, and the lead concentrations in these foods have declined dramatically, for several reasons⁵. From the standpoint of lifetime exposures to lead,

⁵ Food lead concentrations have decline considerably since the period 1982-85 with the removal of cans with

this decline has been relatively recent (i.e. since 1985), and the historical food lead exposure of an individual born in the 1960s or earlier can be a significant component of the total body burden of lead.

Select the "Diet" tab to move to the Dietary exposure mode as shown on Figure 19. The common features (age categories, dates, option buttons) are the same, and the rows of data are similar (amounts ingested, lead concentration) as for the air window. The default setting is for a Market Basket Survey Diet subdivided into the main items of Fruit, Vegetables, Meat and Fish. To access these alternate sources of food exposure, click on the window tab at the lower left corner shown in Figure 19 with the "Market Diet" label. Drop this window down to reveal options for Home Garden, Recreational Fishing/Hunting and Subsistence Fishing/Hunting. Figure 20 shows the selection of the Home Garden option that includes garden fruit and garden vegetables, each with a line for intake percent of diet and lead concentration. The intake percent has a default of zero, so that no exposure will be added to the total diet until a number is added to this row. You can do this one column at a time by entering a number in each cell or, for all age groups by entering a value in the window on the lower left. Note that the default lead concentrations are the same as for the market fruit above, so that no change will be seen in the diet until new lead concentrations for garden fruits are added, based on situational information.

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retirement
tart Date	08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
uration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
larket Fruit Intake [g/day]	13	55	60	70	75	75	85	85	70
arket Fruit Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arket ∀egetable Intake [g/day]	70	70	92	113	150	160	170	170	160
arket Vegetable Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arket Meat Intake [g/day]	20.2	75.4	85.6	115	145	160	160	150	130
arket Meat Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arket Fish Intake [g/day]	6.1	8	12.2	18.1	20	22	22	20	20
larket Fish Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01

Figure 18: Options for Dietary exposure.

soldered side seams. Air lead concentrations have also decreased, reducing the amount of lead in food crops.

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retirement
tart Date	10/01/1915	04/01/1916	09/30/1919	09/29/1921	09/27/1928	09/26/1934	09/21/1955	09/14/1980	09/10/1997
uration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
arket Fruit Intake [g/day]	13	55	60	70	75	75	85	85	70
arket Fruit Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arket Vegetable Intake [g/day]	70	70	92	113	150	160	170	170	160
arket Vegetable Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arket Meat Intake [g/day]	20.2	75.4	85.6	115	145	160	160	150	130
arket Meat Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arket Fish Intake [g/day]	6.1	8	12.2	18.1	20	22	22	20	20
arket Fish Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arden Fruit Intake [% of diet]	0	0	0	0	0	0	0	0	0
arden Fruit Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
arden Vegetable Intake [% of diet]	0	0	0	0	0	0	0	0	0
arden Vegetable Pb Conc. [ug/g]	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01	0.01
ome Garden 💌]+			. [] U:	se Default Hist	orical Pb Conce	ntrations Ga	ardn Fruit Intake[% di	et 💌 🛛 Set D

Figure 19. Dietary exposure with garden fruits and vegetables.

The exposure information for garden fruits and vegetables may be entered in a similar manner, again with no change in the exposure unless values on the table are different from their default values.

Other options include modifications if: 1) produce is home grown, or 2) recreational fishing/hunting, or 3) subsistence fishing/hunting are significant contributors to the diet. The Recreational Fishing and Hunting option provides the opportunity to account for lead exposure from the occasional ingestion of foods in the fish and meat category. Enter the data for recreational hunting/fishing category through the same process as for home gardening, and select garden fruit and vegetable options as well, if they apply.

The third option is Subsistence Hunting and Fishing, and in this case the intent is to replace all market basket meats and fish with meat and fish from non-market sources. This option may also be used in conjunction with the home gardening option, in which case the assumption is that all food comes from non-market sources.

Finally, it is important to note that alternate food sources may not necessarily be lower in lead content than market food. Historically, the main source of lead in market food was lead from cans with a lead-based seal on the side seam. These were removed from the commercial food production process during the period 1982-1985, resulting in a dramatic decline in the lead concentration of canned food. Likewise, the decline in air lead concentrations during the 1980s, with the banning of tetraethyl lead from gasoline, also resulted in a substantial reduction of lead in crops.

2.1.1.3 Dust

Household dust (Figure 20) is assumed to be a mixture of deposited particles of atmospheric origin, tracked-in soil, deteriorating paint, dust from local commercial industries, and on occasion, certain dust generating activities related to the home environment. There are several options available with this exposure group. The first decision is to determine the manner in which the dust sample was taken: a wet wipe sample gives μg Pb per area sampled (μg Pb/cm²); a dry vacuum sample can give μg Pb per area sampled or per gram of dust collected. If the total mass of dust sampled is available, then this can be related to an ingestion rate of dust mass/day, usually about 100 mg dust/day. If the mass of dust is not known and only the area sampled is known, then the modeler must base the exposure estimate on the hand-to-surface area contact and the hand-to-mouth activity of the individual. Furthermore, the efficiency of transfer at the time of surface contact, which depends largely on the dampness of the hand, should also be estimated.

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retire
tart Date	08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
ouration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
otal Dust Ingested [mg/day]	40	85	135	85	50	40	40	40	30
esidential Dust Pb Conc [ug/g]	200	200	200	200	200	200	200	200	200

Figure 20: Exposure options for dust exposure.

The choice between the concentration method (Conc) (Figure 21) or the area method (Area) (Figure 22) can be made by dropping down the window in the lower part of the dust screen. The default option for Residential Dust using concentration measurements is shown below.

The second dust decision is to determine the important sources of dust other than the home that might be associated with school, occupational or recreational activities. If this information is known or can be determined, these options can be implemented by dropping the window in the lower left part of the dust screen, the same window as the choice of [Conc] vs. [Area]. The alternate locations are the same options as for air, and are zeroed out until specific information is entered by the modeler.

In determining how to distribute the dust exposure over several alternate sources, the modeler must make an estimate of the percent of total dust ingestion that occurs at each location, and the

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retire
Start Date	08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
)uration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
fotal Dust Ingested [mg/day]	40	85	135	85	50	40	40	40	30
Residential Dust Ingested [% Total]	100	100	100	50	50	50	50	50	50
Residential Dust Pb Conc [ug/g]	200	200	200	200	200	200	200	200	200
School Ingested Dust [%]	0	0	0	20	20	20	20	20	20
chool Ingested Dust Pb Conc [ug/g]	0	0	0	0	0	0	0	0	0
occupational Ingested Dust [%]	0	0	0	30	30	30	30	30	30
Occupational Ingested Dust_Pb Conc [ug/g]	0	0	0	0	0	0	0	0	0
Recreational Ingested Dust [%]	0	0	0	0	0	0	0	0	0
Recreational Ingested Dust Pb Conc [ug/g]	200	200	200	200	200	200	200	200	200
	- Historical	Dust Ph Con	centration						

Figure 21: Drop Down window for the Dust component of Media Exposure with the options for Other Dust Sources (as Concentration) shown in the lower left drop-down window.

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retire
Start Date	08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
Duration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
Fotal Contact Surface Area [m2/day]	0.06	0.13	0.2	0.13	0.07	0.06	0.06	0.06	0.04
Residential Dust Pb Loading [ug/m2]	135	135	135	135	135	135	135	135	135

Figure 22: Window for the Dust component of Media Exposure with the options for Residential Dust sources (Area) shown in the sub-window at the lower left.

concentration of lead in that dust. The contribution of lead from dust sources other than the home environment will remain at zero until the modeler makes the appropriate alternate dust source entries. Another option using the lead loading approach allows for exposure from school, occupational and recreational activities (Figure 22 and Figure 23).

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retire
Start Date	08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
Duration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
otal Contact Surface Area [m2/day]	0.06	0.13	0.2	0.13	0.07	0.06	0.06	0.06	0.04
Residential Contact Surface Area [% Total]	100	100	100	100	100	100	100	100	100
Residential Dust Pb Loading [ug/m2]	135	135	135	135	135	135	135	135	135
School Contact Surface Area [% Total]	0	0	0	0	0	0	0	0	0
School Dust Pb Loading [ug/m2]	0	0	0	0	0	0	0	0	0
Occupational Contact Surface Area [% Total]	0	0	0	0	0	0	0	0	0
Occupational Dust Pb Loading [ug/m2]	0	0	0	0	0	0	0	0	0
Recreational Contact Surface Area [% Total]	0	0	0	0	0	0	0	0	0
Recreational Dust Pb Loading [ug/m2]	0	0	0	0	0	0	0	0	0

Figure 23: Window for the dust component of Media Exposure with the options for several Dust Sources (as Area and Loading) shown in the sub-window at lower left.

2.1.1.4 Drinking Water

As with the previous routes of exposure, conservative estimates of drinking water consumption are set at default values of one liter of fully flushed drinking water per day (Figure 24). This would include tap water used to prepare coffee or tea, and to reconstitute beverages such as frozen juices. It would not include canned or bottled soft drinks, which are treated as food items. On this basis, it is likely that the mean consumption of drinking water is higher than one liter, possibly 5-6 liters, but this varies by age, climate, season and several other factors, so that a conservative estimate seems appropriate as a baseline starting point. Lead exposure from drinking water is calculated as $\mu g Pb/L x L/day = \mu g Pb/day$.

Other options allow for inputs from first draw water, fountain water and bottled water. These options may be accessed through the drop-down window at the lower left, by clicking the option "Other Water Sources" (Figure 25). This will add options for First Draw Tap Water, Drinking Fountain Tap Water, and Bottled Water. Unlike previous exposure sources, the Drinking Water options are fully implemented when the "Other Water Sources" button is selected. To change the default values, note that each source has a percent-of-total-volume option and this percent may be changed for each source and age group. Remember that bottled water and perhaps drinking fountains were not readily available in the early 1900s; likewise, there is little reliable site specific information on historical drinking water lead concentrations.

a1. Media Exposure									
Air Diet Dust Water Pica	Dermal								
Parameter	Postnatal T	Toddler Pi	reschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retirement
Start Date	08/22/1915 02	2/21/1916 08	8/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
Duration (days)	182 12	276 72	29	2554	2189	7664	9124	6204	2919
Total Vol. of Water Consumed [L/day]	1 1	1		1	1	1	1	1	1
Flushed Pb Conc [ug/L]	5 5	5		5	5	5	5	5	5
Flushed Water									
Clear All Default	All to De	efault						OK Can	cel Help

Figure 24: Window for the Water component of Media Exposure with the default option of fully Flushed Water in the sub-window at the lower left.

	Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retiremen
S	tart Date	08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
D	uration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
T	otal Vol. of Water Consumed [L/day]	1	1	1	1	1	1	1	1	1
-1	ushed Amt [% of Total Volume]	84	84	84	84	84	84	84	84	84
-1	ushed Pb Conc [ug/L]	5	5	5	5	5	5	5	5	5
i	irst Draw Tap Amt [% of Total Volume]	5	5	5	5	5	5	5	5	5
-i	irst Draw Tap Pb Conc. [ug/L]	4	4	4	4	4	4	4	4	4
F	ountain Amt [% of Total Volume]	10	10	10	10	10	10	10	10	10
F	ountain Pb Conc [ug/L]	7	7	7	7	7	7	7	7	7
B	ottled Amt [% of Total Volume]	1	1	1	1	1	1	1	1	1
Ð	ottled Pb Conc [ug/L]	4	4	4	4	4	4	4	4	4
	Ther Water Sources									

Figure 25: Window for the Water component of Media Exposure with the options for Other Water Sources shown in the sub-window at lower left.

[Pica and Dermal Exposure Factors Are Discussed in Chapter 3]

2.2 EDITING EXPOSURE ACTIVITY PATTERNS

Activity patterns determine the amount of time (hours) the modeled individual spends at each activity. There are 168 hours in a week. For each age group, these hours are distributed, as applicable, among four exposure locations (residential, school, occupational, and recreational) using the 168 hours, minus the designated sleep hours per day. Within each location, the modeler may designate a percentage of the time spent indoors, allowing for separate estimates of indoor/outdoor air and dust lead exposures.

Three factors determine exposure: 1) the amount of time spent at a particular activity; 2) the amount of "media" ingested during that time; and 3) the concentration of lead in that media. For food and drinking water lead concentrations, the location may not be a factor if supermarket food and municipal water are consumed, as these are usually monitored closely and tend to be the same throughout the United States. Air and dust lead concentrations are likely to be different by location and are treated as such in this model.

During the postnatal period, which lasts six months, the child is presumed to be in a highly protective environment with limited access to dust collecting surfaces and under continuous observation by a caregiver. Consequently, only one exposure environment, "Residential" is provided in the model. In the event there needs to be a second environment, this can be accomplished by reducing the time frame for Postnatal status and expanding the Toddler time frame.

F	Postnatal Todd	ler Preschool	Youth
	Residential	168.0	90.0
	School		0.0
	Occupational		0.0
	Recreational	0.0	0.0
	Weekly Total	168.0	
	Daily sleep hou	rs: 12.0	

Figure 26: Drop down window for the Activity Patterns component of Media Exposure showing default values for activity patterns for Postnatal.

2.3 Summary and Conclusions

This concludes the introduction to the General Guidance for Media Exposure. Most risk assessment problems for lead exposure can be handled with the development of an exposure scenario involving the model features discussed in this section followed by a statistical

extrapolation of the predicted blood lead concentration for an individual to a mean blood lead concentration for a population of similarly exposed individuals.

At this point, you may find it helpful to enter your own exposure settings and save the file as a useful starting point for further research or investigation. Section 4 provides guidance on special features and work-arounds for special problems and situations. If these features do not apply to your needs, you may proceed directly to Section 5 and begin adjusting the model for lead absorption.

When you have completed constructing the exposure scenario for your subject, it is wise to give it a brief reality check. To determine whether you have correctly modeled the exposure of the individual, imagine yourself in that person's place for one week and determine whether or not you would have captured all or most of your exposure for that time. Remember that, in some cases, you may account for exposure at two locations as they were one. That is, your exposure at the office may be, for all practical purposes, the same as at the library or at a concert hall, and could be treated as one location.

12. Model Run	
Exposure Absorption	
Biokinetics	Edit Model Settings
C	ontinue Cancel

Figure 27: Following set of the Exposure Scenario, the model run is initiated by checking the Exposure and Absorption boxes. The Biokinetics box will be automatically checked.

At this point you may choose to run just the exposure portion of the model by selecting that box in the Model Run Window (Figure 27). In this case, you will get two files, *.mod and *.mpa (Figure 28). The first (*.mod) is a complete list, day by day, of exposure and may contain up to 32841 lines of data for the six exposure categories (air, diet, water, dust, dermal and pica).

The second file (*.mpa) contains a summary of these six exposure categories for each age category selected. You will also have an opportunity to view these files after the absorption/biokinetics portion of the model has been executed.

Select the expos	ure model output	t file to view.
Review1.mod		
Review1.mpa		
	View	Continue

Figure 28. Two files are available for review, *.mod and *.mpa.

	Interval Start Date	Interval End Date	Time Step (Day)	Air(ug Pb/TS)	Diet(ug Pb/TS)	Water(ug Pb/TS)	Dust(ug Pb/TS)	Dermal[ug/T 📥
1	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
2	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
3	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
4	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
5	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
6	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
7	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
8	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
9	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
10	09121997	09092005	1	0.3859	3.8000	420.0000	6.0000	0.00(
č	100101007	0000005			0.0000	100,0000		N

Figure 29. Line by line output for the exposure module. *.mod

	outp					
Age Category	Air(ug Pb)	Diet(ug Pb)	Water(ug Pb)	Dust(ug Pb)	Dermal(ug Pb)	Pica(ug Pb)
Postnatal	59	198	76440	1456	0	0
Toddler	1124	2659	535920	21692	0	0
Preschool	580	1821	306180	19683	0	0
Youth	1997	8073	1072680	43418	0	0
Adolescent	2002	8537	919380	21890	0	0
Young Adult	10175	31958	3218880	61312	0	0
Middle Age	10001	39871	3832080	72992	0	0
Early Retirement	4404	26366	2605680	49632	0	0
Late Retirement	1126	11092	1225980	17514	0	0
1						
Export to		-				
Enporce		<u> </u>				

Figure 30. Total exposure by age for each exposure category (*.mpa)

1e1.E	dit Model Settings					X
	Parameters	Values	Units	Descriptions	Keys	
1	Exposure Age	0.002740	years	Age at acute exposure or beginning of chronic exposure	expage	
2	Last Day	32850	days	Maximum number of days	endday	
3	Maximum Cycles	32850		Maximum number of cycles	ncycle	
4	Fetal Exposure On/Off	0		Fetal Exposure Switch (1=On 0=Off) INTERNAL	ifetal	
5	Fixed Length Detta Option	0		Use fix length detta (0=Variable >0=Fixed/TimeStep)	delt0	
6	Detta Step Lengths	Edit		Delta step lengths by age range	delta	
7	Output Step Interval	100	days	Write output to file only on these steps	iskip	
8	Acute/Chronic	2		Switch for acute or chronic (1=Acute 2=Chronic)	iacute	
9	Mode of Intake	3		Selection for mode of intake (0=Injection 1=Inhalation 2=Ingestion 3=Combination	inmode	
10	Linear/Nonlinear Model	1		(0= Linear Model 1=Nonlinear Model)	irbc	1
11	RBC Threshold Concentration	60	ug/dl	Lead concentration on RBC above which a nonlinear model is used	rboni	
12	Nonlinear Parameter 2	350		Nonlinear Parameter 2	satrat	
13	Power	1.5		Power	power	
14	Chelation On/Off	0		Chelation Switch (0=Off 1=On)	ichel	
15	Fixed Delta	1.0	days	Length for Fixed Delta	deltfix	
16	Body Size Curve (on/off)	0		Computation will use growth curve (1=on, 0=off)	bUseBodyS	i
17	Bone Computation?	Use Leggett		Bone computation to be used (i.e. Leggett or O'Flaherty)	sBoneComp	j I
18						1
19						
20						
21						
22						
2		1	1			
		Default		Save	Cancel	

Figure 31: Exposure Model options with options to change internal exposure and growth parameters.

SECTION 3: SPECIAL EXPOSURE FEATURES

There are three model features that have been incorporated into the Exposure component of the All Ages Lead Model that reflect special cases with respect to human exposure to lead. These are pica ingestion, dermal absorption, and maternal (or prenatal) exposure. The AALM permits the modeler to explore these routes of exposure with the caveat that these features are under development and have not been fully tested. Furthermore, there is little scientific information available to document the basis for parameter inputs and very little guidance can be given.

3.1 PICA EXPOSURE

This document defines pica as the habitual tendency or craving to ingest non-food items such as paint chips and soil. In the default mode, numeric values for the amount of ingested soil and paint chips are set at 0, so that no contribution from pica will be included in a study without the intentional specification of some positive value for this ingestion rate. To activate the pica mode, simply enter an estimate the amount of soil or paint chips that might be ingested during any of the life stages of the individual, and revise, if necessary, the lead concentration in that soil or dust (Figure 32). Experience will soon show that pica can, at reasonable ingestion rates, dominate the exposure of the individual.

Parameter	Postnatal	Toddler	Preschool	Youth	Adolescent	Young Adult	Middle Age	Early Retirement	Late Retirement
tart Date	08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/18/1928	08/17/1934	08/12/1955	08/05/1980	08/01/1997
uration (days)	182	1276	729	2554	2189	7664	9124	6204	2919
ead Conc. of Soil (ug/g)	200	200	200	200	200	200	200	200	200
oil Ingested (g/day)	0	0	0	0	0	0	0	0	0
ead Conc. of Paint Chip (ug/cm2)	500	500	500	500	500	500	500	500	500
aint Chip Ingested (cm2/day)	0	0	0	0	0	0	0	0	0

Figure 32: Window for the Pica component of Media Exposure.

3.2 DERMAL EXPOSURE

The frequency of lead exposure via the dermal route remains poorly defined in the literature, especially the mechanism by which it is absorbed into the body. Nevertheless, dermal absorption can be an important contributor to children playing in dusty situations, such as sand boxes. It may also occur with occupationally-exposed subjects working with fine dusts and under hot,

sweaty conditions. The default values for dermal absorption parameters are arbitrary and not based on specific scientific studies. Like pica, dermal exposure is set to zero by default values for soil and paint chip ingestion (Figure 33) and may be activated by changing these values to some positive number.

08/22/1915	02/21/1916	08/21/1919	08/20/1921	08/08/028	00474004		00 ID F H 000	
182	4.070			00/10/1320	08/17/1934	08/12/1955	08/05/1980	08/01/1997
	1276	729	2554	2189	7664	9124	6204	2919
40	40	40	40	40	40	40	40	40
0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5
) 200	200	200	200	200	200	200	200	200
	0.3 200	200 200	200 200 200	200 200 200 200 200	200 200 200 200 200 200	200 200 200 200 200 200 200 200	200 200 200 200 200 200 200 200 200 200	200 200 200 200 200 200 200 200 200 200

Figure 33: Window for the Dermal component of Media Exposure.

3.3 HISTORICAL LEAD EXPOSURES

In the case where the exposure scenario extends back in time to the early 1900s, the model default values for some parameters may not be accurate in light of high exposures, due to the use of food cans with lead soldered side seams and gasoline with tetraethyl lead. This feature, still under development, will allow the modeler to insert historical food data from 1982 to the present and actual air lead measurements from the early 1970s. This information will also provide guidance for reasonable estimates of exposure to these media for the years prior to the time accurate date were available.

3.4 MATERNAL EXPOSURE (Conceptual Design)

Maternal exposure results from the transfer of lead from the mother to the developing fetus during the nine months of fetal development. The model treats this exposure as a single source (the mother's blood) that determines the tissue concentrations of the infant at the time of birth. Thus, with this component activated, exposure begins at the formation of the embryo. Until birth, this developing fetus is modeled as a separate compartment of the mother's model. At the time of birth, a new file is created, representing the child's accumulated exposure with tissue and organ concentrations derived from maternal exposure rather than the default for the postnatal infant model.

SECTION 4: ABSORPTION MODULE

4.1 MODELING LEAD ABSORPTION

Similar to the Exposure Module, the Absorption Module has a set of parameters that represent specific body functions related to lead absorption, and these parameters can be adjusted to reflect the absorption of lead from the lungs and the digestive tract. A third absorption option would be dermal absorption when it is fully implemented. As with the Exposure module, all parameters in the Absorption module can be changed by replacing the default values with "exploratory" values that the modeler may find applicable to a specific situation. Two groups of parameters can be edited in this module, lung absorption and digestive tract absorption.

Lung absorption is derived from the Leggett model, where there are four lung compartments, each with two absorption parameters (Figure 34). The first parameter (Fraction) is the fraction of the total amount of lung lead that is deposited in the lung, and the second parameter is the amount per day that is transferred to the blood compartment.

For ingested food, the rate of movement through the gastrointestinal tract varies by more than just the age of the individual. No distinction is made in the model among the three sections of the small intestine (duodenum, jejunum, ileum) or their function. The gastro-intestinal absorption of lead in children and adults has been thoroughly described (Mushak, 1991). Issues that are of special interest for modeling purposes include the physical process of preparing the food for nutrient absorption, as well as the physiological and biochemical process of absorption. Following digestion in the stomach, absorption occurs across the epithelial lining of the intestinal mucosa, the additional folds of the submucosa, and the microscopic folds of the villi and microvilli. This gives the total absorbing surface not only a large physical area across which nutrients can pass, but a highly complex microenvironment whereby non-nutrients such as lead can pass into the blood stream with a high degree of localized variability.

The rate of movement through the gastrointestinal tract is the reciprocal of the number of days the lead would reside in the G.I. tract (Leggett, 1993). During this time, several digestive and absorption mechanisms break down the food ingested into simpler chemical and biochemical components that can be transported across the wall of the small intestine. Among these are a group of nutrient elements, including calcium, iron, and phosphorus, that are chemically similar to lead. Lead is postulated to follow these essential nutrients by a combination of three mechanisms: diffusion, pinocytosis, and facilitated transport (Leggett, 1993).

Stable isotope studies have shown, at least for adults, that fasting can greatly increase the rate of absorption (Maddeloni et al, 1998). Whereas the mechanism is not fully understood, it is highly likely that calcium plays a significant role is preventing the uptake of lead. Theoretically, if there is sufficient calcium in the digestive tract, lead will be excluded to some degree. But, at a microsite of absorption that is partially isolated by the intricate folds of the submucossa, the calcium may be absorbed during fasting periods and not replenished, leaving the lead available at higher relative concentrations. Without modeling calcium ingestion and absorption at the same time as lead ingestion and absorption, it is difficult to correctly estimate the true rate of lead absorption.



Figure 34. Schematic diagram of the Leggett Model. This model shows the routes of absorption as the Respiratory Tract (RT Tract) and Gastrointestinal Tract (GI Tract) contributing lead to the Diffusible Plasma compartment. A further discussion of this model may be found in Chapter 5.

Recent data from stable isotopic studies suggest that the GI absorption fraction in children is similar to adults (Gulson, 2005).

It is clear that much remains to be discovered with respect to the mechanisms of gastrointestinal absorption of lead. With its open architecture and flexible application of input parameters, the All Ages Lead Model is designed to facilitate continued research on such critical matters as the intake and absorption of lead, and to accommodate the results of that research through the adjustment of the key absorption parameters.

Descripti	on: Delta step lengths by	age range		Interpolation Type: Linear			Unit
	A	В	с	D	E	F	
1	Age Cut-off	1000	1900	2800	6700	21900	
2	Delta Step Lengths	1.0	1.0	1.0	1.0	1.0	
		A 110 1				0	

Figure 35: Drop down window for reviewing and revising Delta step values. These are the internal time step functions for Absorption and Biokinetics, and, by default are set at 1 day of simulation. They can be varied individually. The model also interpolates between these steps using the Linear option. The other options for Interpolation Type are Floor, Ceiling, and Index, which have not been evaluated for AALM applications, but are included for research purposes.

1c. Edit Absorption 🛛 🔀						
	Parameters	Values	Units	Descriptions	Keys	
1	Dermal Absorption Factor	0		Dermal absorption factor should be between 0 to 1 with 0 as default	kdermal	
2	Lung Compartment 1 Fraction	0.25		Fraction of Pb deposited amount assigned to lung compartment 1	r1	
3	Lung Compartment 2 Fraction	0.35		Fraction of Pb deposited amount assigned to lung compartment 2	r2	
4	Lung Compartment 3 Fraction	0.3		Fraction of Pb deposited amount assigned to lung compartment 3	r3	
5	Lung Compartment 4 Fraction	0.1		Fraction of Pb deposited amount assigned to lung compartment 4	r4	
6	Lung Compartment 1 Rate	16.6	/day	Transfer rate from lung compartment 1	br1	
7	Lung Compartment 2 Rate	5.54	/day	Transfer rate from lung compartment 2	br2	
8	Lung Compartment 3 Rate	1.66	/day	Transfer rate from lung compartment 3	br3	
9	Lung Compartment 4 Rate	0.347	/day	Transfer rate from lung compartment 4	br4	
10	GI Absorption Fraction	Edit		Gastrointestinal absorption fraction	af1	
11	Rate of Movement thru GI Tract	Edit	/day	Scales rate of movement through the gastrointestinal tract	agscal	
12	Fraction to GI	0.04		Fractional movement to the gastrointestinal tract	ciliar	
13	Stomach Transfer Rate	24	/day	Adult transfer rate from stomach	rstmc	
14	Small Intestine Transfer Rate	6	/day	Adult transfer rate from small intestine	rsic	
15	Upper Large Intestine Transfer Rate	1.85	/day	Adult transfer rate from upper large intestine	ruli	
16	Lower Large Intestine Transfer Rate	1	/day	Adult transfer rate from lower large intestine	rlli	
4		10				
					· ·	
Default Save C						

Figure 36: Window for editing Absorption Parameters.

4.2 ABSORPTION OF INHALED LEAD

This represents the amount of lead transferred from the lungs (through mucociliary lift into the stomach) and liver to the gastrointestinal tract where it is eliminated as feces. It does not include ingested lead that was not absorbed during digestion. The coefficient for the deposition fraction in feces (0.006) is smaller (slower) than urine (0.015) and is also constant for all age groups. Gastrointestinal absorption is treated as a simple percent of absorption through the small

DRAFT: DO NOT QUOTE OR CITE
intestine, with no variation given to the type of food consumed or daily rhythms. This function includes lead that has passed from the stomach to the duodenum and lead in bile excreted from the liver.

4.3 ABSORPTION OF INGESTED LEAD

10). Edit	Values for GI Absor	ption Fractio	n				
C) escriptio	on: Gastrointestinal absor	ption fraction		Interpolation	Type: Linear	▼ Unit:	
Γ		A	В	с	D	E	F	G
	1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
	2	GI Absorption Fraction	0.45	0.45	0.3	0.3	0.3	0.3
	<							>
	Import	Data Export Data	Add Colu	mn Delete Coli	umn	Plot	Save (Cancel

Figure 37: Drop down window for gastrointestinal absorption showing default values as a decimal percentage by age range (in years). The age breakpoints are not the same as the breakpoints for exposure that are assigned on the Opening Menu. The absorption values may be reset within the range of 0.0-1.0. Data are plotted on Figure 38.



Figure 38: Plot of the default values for the gastrointestinal absorption fraction. Note that the GI absorption fraction for the first three months, when the diet is largely liquid, is much higher than at age one year and higher.

Descript	ion: Scales rate of movem	ent through the g	jastrointestina	Interpolation	Type: Linear	▼ Unit	/day
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Rate of Movement thru	1.66667	1.66667 1.66667	1.66667	1.33333	1	
< Import	Data Export Data	Add Colur	mn Delete Colu	umn	Plot	Save (> Cancel

Figure 39: Default values for the gastrointestinal rate of movement function. Units are fraction of whole mass per day (e.g. from age 0-1, the child passes food at the rate of 1.6 volumes per day). Data are plotted on Figure 40.



Figure 40: Plot of the default values for the rate of gastrointestinal movement function. This function, in default form, reaches maturity at age 15 years and remains constant thereafter.

SECTION 5: BIOKINETICS MODULE

There are three groups of parameters in this module for which parameters may be edited. These are: Model Settings, Biokinetic Parameters, and Body Size. As with some absorption parameters in the previous section, there are several biokinetic parameters that may be adjusted according to age, thus giving considerable flexibility in the application of the model to a wide range of populations.

5.1 MODEL SETTINGS

This set of seventeen functions (Figure 42), shown also in Table 1) has three types of parameters: that are assigned default values, at the initialization of the model. Some of these are simply on/off switches that turn on or off an optional model feature selected by the modeler for a particular function. Other parameters for which default values have been assigned can be changed by the modeler. Two parameters, Delta Step Lengths and Bone Computation, have drop down windows that allow the modeler to select a specific value or option. The Delta step lengths are the 'inflection' points along the age range, and are not necessarily the same points as for age categories for the earlier exposure model component. This option allows for changes, for example, in cortical bone turnover to be specified for times different from the age categories related to growth and exposure.

5.1.1 PARAMETERS WITH EXPLICIT VALUES

These eight parameters have only one variable each (i.e., they are not age-dependent). This default value can be modified by the modeler. For example, the Exposure age (age at the beginning of exposure), which has a default of one day (0.00274 years) can be changed to 1, which would start the exposure at age one year.

The Last Day has a maximum value of 32850 (days), as does the Maximum Cycles. If either of these is changed, the lower value would prevail.

The Output Step Interval simplifies the output by tallying the results over each 100 day segment, and this interval can be changes as well.

The RBC (Red Blood Cell) Threshold Concentration has a default value of 60 μ g/dL, below which lead is taken up by red blood cells linearly; above that concentration, the uptake process is non-linear. One of the parameters of this nonlinearity (Parameter 2) has a value of 350, and a part of the equation is expressed to the power of 1.5. The Fixed Delta has a default value of 1.0 days, which may be changed by the modeler.

Of the seven On/Off switches, five simply turn a specific function on or off (Fetal Exposuredefault off; Fixed Length Delta option – default variable i.e. off; Acute/Chronic – default chronic exposure; Linear/Non Linear Model – default Linear; Chelation – default off; Body Size Curve – default off)⁶

	DEFAULI	
PARAMETER	VALUE	TYPE
Exposure age	0.002740	Function with default value
Last Day	32850	Function with default value
Maximum Cycles	32850	Function with default value
Fetal Exposure	0	On/Off Switch
Fixed Length Delta Option	0	On/Off Switch
Delta Step Lengths	EDIT	Drop down window with age-related variables
Output Step Interval	100	Function with default value
Acute/Chronic	2	On/Off Switch
Mode of Intake	3	On/Off Switch
Linear/Nonlinear	1	On/Off Switch
RBC Threshold Concentration	60	Function with default value
Nonlinear Parameter 2	350	Function with default value
Power	1.5	Function with default value
Chelation	0	On/Off Switch
Fixed Delta	1.0	Function with default value
Body Size Curve	0	On/Off Switch
Bone Computation	Use Leggett	Drop down window to select the Leggett or O'Flaherty method of bone lead kinetics

One switch has four possible values: Mode of Intake: 0= injection, 1=inhalation, 2=ingestion, 3=combination; the default is combination and the model has not been tested for any other mode.

Table 1. Biokinetic model settings showing the three types of editable variables (Function
with default value, On/Off Switch, or Drop Down window) common to several
model options.

The second drop down window allows the modeler to select the bone methodology of either Leggett or O'Flaherty. The Leggett method is generally considered to be anatomically based and offers the modeler several options within the complex of bone compartments. The O'Flaherty method is physiologically based and promotes the concept of lead following the physiological pathway of calcium in bone tissue (O'Flaherty et al, 1998). This feature is a relatively new addition to the model and has not been fully implemented. The schematic diagram for the O'Flaherty model appears in Figure 41.

⁶ Note that one of the earlier uses of the Leggett model involved the removal of lead from blood by chelation. This feature is rarely used, but was retained against the possibility of some future use.



Figure 41. Schematic Diagram of O'Flaherty Model. This model addresses the physiological process of lead movement in bone tissue and is currently being tested for compatibility with the All Ages Lead Model

	Parameters	Values	Units	Descriptions	Keys	1
	Exposure Age	0.002740	years	Age at acute exposure or beginning of chronic exposure	expage	ĺ
	Last Day	32850.000000	days	Maximum number of days	endday	1
	Maximum Cycles	32850.000000		Maximum number of cycles	ncycle	1
	Fetal Exposure On/Off	0		Fetal Exposure Switch (1=On 0=Off) INTERNAL	ifetal	1
	Fixed Length Delta Option	0		Use fix length detta (0=Variable >0=Fixed/TimeStep)	delt0	1
	Detta Step Lengths	Edit		Detta step lengths by age range	delta	1
	Output Step Interval	100	days	Write output to file only on these steps	iskip	1
	Acute/Chronic	2		Switch for acute or chronic (1=Acute 2=Chronic)	iacute	1
	Mode of Intake	3		Selection for mode of intake (0=Injection 1=Inhalation 2=Ingestion 3=Corr	n inmode	1
0	Linear/Nonlinear Model	1		(0= Linear Model 1=Nonlinear Model)	irbc	1
1	RBC Threshold Concentration	60	ug/dl	Lead concentration on RBC above whcih a nonlinear model is used	rboni	1
2	Nonlinear Parameter 2	350	_	Nonlinear Parameter 2	satrat	
3	Power	1.5		Power	power	
4	Chelation On/Off	0		Chelation Switch (0=Off 1=On)	ichel	
5	Fixed Delta	1.0	days	Length for Fixed Delta	deltfix	
6	Body Size Curve (on/off)	0		Computation will use growth curve (1=on, 0=off)	bUseBodySiz	έ
7	Bone Computation?	Use Leggett		Bone computation to be used (i.e. Leggett or O'Flaherty)	sBoneComput	à
8	· · · · · · · · · · · · · · · · · · ·			· · · _ · _ · _ · _ ·		1
9						1
0						1
1						1
2						1
3						1
A	1				1	1

Figure 42: Window for editing the model settings. Features other than the Delta Step Length have options than do not vary with age.

escripti	on: Delta step lengths by	age range		Interpolation 1	• L	Jnit:	
	A	В	с	D	E	F	
1	Age Cut-off	1000	1900	2800	6700	21900	
2	Delta Step Lengths	1.0	1.0	1.0	1.0	1.0	
2	Delta Step Lengths	1.0	1.0	1.0	1.0	1.0	

Figure 43: Drop down window for editing the default Delta Step lengths.

5.2 BIOKINETIC PARAMETERS

The second option in the Biokinetic Parameter Window is the option to "Edit Biokinetic Parameters," a window that lists 47 parameters, all of which have numeric values or age-related sets of numeric values. Most of these are directly or indirectly related to the arrows connecting the various compartments of the schematic model, Figure 45 (Pounds, J.G.; Leggett, R.W. (1998).

	Parameters	Values	Units	Descriptions
	Pb Decay Rate	0		Rate of decay for Pb
	Cortical Bone Turnover	Edit		Cortical bone turnover rate by age range
	Trabecular Bone Turnover	Edit		Trabecular bone turnover rate by age range
	Transfer from Cortical Surface to Blood	Edit	/day	Transfer rate from the cortical surface to blood by age range
	Transfer from Trabecular Surface to Blood	Edit	/day	Transfer rate for tabecular surface to blood by age range
	Cortical Surface to Volume Transfer	Edit		Fraction transfer from the cortical bone surface to the volume of
	Trabicular Surface to Volume Transfer			Fraction transfer from the trabecular bone surface to the volume
;	Total Transfer from Exchange Bone Volume	Edit	/day	Total transfer rate from exchange bone volume by age range
	Transfer from Exchange to Non-exchange Volume	Edit		Fraction transfer from exchange volume to non-exchange volume
0	Transfer from Liver 1	0.0693	/day	Transfer rate from liver 1
1	Transfer from Kidney 1	0.139		Transfer rate from kidney 1
2	Transfer from Bladder to Urine	Edit	/day	Transfer rate from uninary bladder to urinary tract by age range
3	Transfer from Liver 2	Edit	/day	Transfer rates from liver section 2 by age range
4	Transfer from Kidney 2	Edit	/day	Transfer rate to kidney section 2 by age range
5	Transfer from Fast Soft Tissue	2.079	/day	Transfer from fast soft tissue 1 by age range
6	Transfer from Intermediate Soft Tissue	0.00693	/day	Transfer from intermediate soft tissue 2 by age range
7	Transfer from Slow Soft Tissue	0.00038	/day	Transfer from slow soft tissue 3 by age range
8	Transfer Rates from Brain	Edit	/day	Transfer rates of lead from brain by age range
9	Deposition Fraction in Urine	0.015		Deposition fraction in urine
20	Deposition Fraction in Feces	0.006		Deposition fraction in feces
21	Deposition Fraction in Sweat	0.0035		Deposition fraction in sweat
22	Intermediate Soft Tissue to Excretion Fraction	0.4		Fraction transfer from intermediate soft tissue to excretion
23	Deposition Fraction in Bone	Edit		Deposition fraction in bone by age range
54. I	T	- 11	4	For the offered descendence of the forther descendence in the second sec

Figure 44: Window for editing biokinetic parameters. The bottom half of this window includes parameters 24 through 47 and appears in Figure 78.

5.2.1 Pb DECAY RATE.

This rate of decay for lead has a default of 0 and should be left at this value as it is not currently used. It is reserved for future use.



Figure 45. Schematic Diagram of Leggett Model. The small numbers correspond to the numbered parameters in Figure 44 and Figure 78. The central exchange compartment is *diffusible* plasma. Bone is represented as having surface (which exchanges with plasma) and volume compartments; the latter simulates slow exchange with the surface and slow return of lead to the plasma from bone resorption. Pounds, J.G.; Leggett, R.W. (1998).

5.2.2 CORTICAL BONE TURNOVER.

This is one of three parameters that regulate the transfer into and within the cortical bone tissue (Figure 44). The value for the turnover rate is age-dependent and decays to a value nearly zero by age fifteen (See Figure 46). This allows changes in the rate of cortical bone turnover to be made for the Age Ranges (Delta Steps) that are specified in Biokinetics Model Settings, line 6. Units

2. Edit V	Values for Cortical B	one Turnove	r				
Descripti	ion: Cortical bone turnove	r rate by age ran	ge	Interpolatio	on Type: Linear	▼ Ut	nit:
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Cortical Bone Turnover	0.0102	0.00822	0.00288	0.00154	0.00089	0.000512
<	_						>
Import	Data Export Data	Add Colu	ımn Delete (Column	Plot	Save	Cancel

Figure 46: Values for cortical bone turnover rate which decreases with age. Data are plotted on Figure 47.



Figure 47: Plot of default values for cortical bone turnover rate by age. This Y-axis scale appears to show values less than zero, but this is a rounding error. Above the age 25 years, this value is 0.000082

are percent fraction per day. The values decrease exponentially from age 0 to 25 years, from 0.01% at birth to 0.00008% at age 25 and older. Interpolation between age groups may be specified as linear, floor, ceiling, or index. During early development, cortical and trabecular bone turnover are modeled at nearly the same rate; by age 25, cortical turnover is about 17% of trabecular.

5.2.3. TRABECULAR BONE TURNOVER RATE.

At one year and younger, the turnover for trabecular bone is modeled at the same rate as cortical (See Figure 48 et seq.). Above one year, the trabecular bone turnover rate decreases slightly less than cortical. Changes to trabecular bone turnover rates can be made for the <u>same</u> Age Ranges or independently of the changes for cortical bone turnover rates.

7. Edit V Descriptio	Edit Values for Trabicular Surface to Volume Transfer								
	A	В	с	D	E	F	G		
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000		
2	Trabicular Surface to V	0.35	0.35	0.35	0.35	0.35	0.35		
<						·	>		
Import	Data Export Data	Add Colu	mn Delete Col	umn	Plot	Save	Cancel		

Figure 48. Age specific values for Surface to Volume in Trabecular Bone. Data are plotted on Figure 49.



Figure 49. Plot of Trabecular Surface to Volume.



Figure 50: Values for trabecular bone turnover rate by age. Data are plotted on Figure 51. Figure

51



Figure 51: Plot of default values for trabecular bone turnover rate by age.

5.2.3 TRANSFER FROM CORTICAL SURFACE TO BLOOD PLASMA.

The model treats the bone surface as the major route of transfer from cortical bone surface to both blood plasma and bone volume. In the case of adults, the transfer from the cortical bone surface is split evenly between transfer to blood plasma (Figure 52) and bone volume (Figure 56), thus the coefficient for ages over 15 is 0.5/day. For the younger age groups, this distribution differs between the two: 0.65/day to blood plasma and 0.35/day to cortical bone volume.

. Edit	Values for Transfer f	. Edit Values for Transfer from Cortical Surface to Blood											
Descript	tion: Transfer rate from the	cortical surface (to blood by a	Interpolation	Type: Linear	▼ Unit:	/day						
	A	В	с	D	E	F	G						
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000						
2	Transfer from Cortical S	0.65	0.65	0.65	0.65	0.65	0.65						
<							>						
	Import Data Export Data Add Column Delete Column Plot Save Cancel												

Figure 52: Dropdown window for editing the parameter for transfer from cortical surface to blood. Data are plotted on Figure 53.





5.2.4. TRANSFER FROM TRABECULAR SURFACE TO BLOOD PLASMA.

Trabecular surface to blood plasma transfer rates are modeled the same as cortical surface to plasma. For the younger age groups, this distribution is 0.65/day to blood plasma (Figure 54) and 0.35/day to trabecular bone volume (Figure 58). For ages over 15 is 0.5/day for each route.

. Edit \	. Edit Values for Transfer from Trabecular Surface to Blood											
Descripti	on: Transfer rate for tabe	cular surface to b	lood by age r	Interpolation	Type: Linear		/day					
	A	В	с	D	E	F	G					
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000					
2	Transfer from Trabecul	0.65	0.65	0.65	0.65	0.65	0.65					
<		1111					>					
Import	Import Data Export Data Add Column Delete Column Plot Save Cancel											

Figure 54: Drop down window for editing the parameter for the transfer rate for trabecular surface to blood. Data are plotted on Figure 55.



Figure 55: Plot of the default values for the transfer rate for trabecular surface to blood.

5.2.6 CORTICAL SURFACE TO VOLUME TRANSFER.

This is the fraction per day that is transferred from the cortical bone surface to the volume of the cortical bone and is equal to [1-(transfer to blood)]. The default values are 0.35 for the younger ages (0-15 years), then 0.5 for the older ages.

. Edit V	Edit Values for Cortical Surface to Volume Transfer											
Descripti	on: Fraction transfer from	the cortical bone	surface to th	Interpolation	Type: Linear							
	A	В	с	D	E	F	G					
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000					
2	Cortical Surface to Volu	0.35	0.35	0.35	0.35	0.35	0.35					
<					<u></u>		>					
				Connect								

Figure 56. Window for editing the value for transfer from cortical surface to volume. Data are plotted on Figure 57.



Figure 57: Plot of the default values for transfer from cortical surface to volume.

5.2.7 TRABECULAR SURFACE TO TRABECULAR VOLUME TRANSFER.

Similar to the cortical surface to volume transfer, this parameter accounts for the fraction of trabecular surface lead that is transferred trabecular volume. The age-dependent default values are the same as the cortical counterpart, 0.35 for children and 0.5 for older children and adults.

. Edit V Descriptio	Alues for Trabicula	Edit Values for Trabicular Surface to Volume Transfer Description: Fraction transfer from the trabecular bone surface to Interpolation Type: Linear Unit:										
	A B C			D	E	E F	G					
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000					
2	Trabicular Surface to V	0.35	0.35	0.35	0.35	0.35	0.35					
<		1111					>					
Import	Data Export Data	Add Colu	mn Delete Colu	umn	Plot	Save (Cancel					

Figure 58: Window for editing the values for transfer from trabecular surface to volume. Data are plotted on Figure 59.



Figure 59: Plot of the default values for transfer from trabecular surface to volume.

5.2.8 TOTAL TRANSFER FROM EXCHANGE BONE VOLUME.

Based on radium data, the removal halftime of lead from the exchangeable bone volume is 30 days, during which time 20% of the lead passes to the non-exchangeable bone tissue and 80% returns to the bone surface, and is the same for cortical and trabecular bone. Although this is an age-dependent variable, the current default value is the same for all age groups (Figure 60).

8. Edit V	Values for Total Tra ion: Total transfer rate from	nsfer from Ex	c hange Bone volume by a	Volume Interpolati	on Type: Index	▼ Ur	nit: /day
	A	В	с	D	E	F	G
1	Indexes	0.000	0.274	1.000	5.000	10.000	15.000
2	Total Transfer from Exc	0.023105	0.023105	0.023105	0.023105	0.023105	0.023105
<		1111					>
Import	Data Export Data	Add Colu	mn Delete Co	lumn	Plot	Save	Cancel

Figure 60: Window for editing values for total transfer from the exchange bone volume. The plot feature has been temporarily disabled.

5.2.9. TRANSFER FROM THE EXCHANGE TO THE NON-EXCHANGE VOLUME.

This parameter combines the cortical and trabecular bone lead transfer from the exchange volume to the non-exchange volume compartments and has the same default value (0.2) for all age ranges (Figure 61).

9. Edit	Values for Transfer	from Exchanş	ge to Non-ex	change Volu	me		
Descrip	otion: Fraction transfer from	exchange volum	e to non-exch	Interpolation	n Type: Index	👻 Uni	:
	A	В	с	D	E	F	G
1	Indexes	0.000	0.274	1.000	5.000	10.000	15.000
2	Transfer from Exchang	0.2	0.2	0.2	0.2	0.2	0.2
<							>
Impo	rt Data	Add Colu	mn Delete Co	luma	Plat	Caua	Cancel

Figure 61: Window for editing values for total transfer from the non-exchange bone volume. The plot feature has been temporarily disabled.

5.2.10. TRANSFER FROM LIVER1.

This is the fraction per day that is transferred from Liver1 to the three compartments: Gastrointestinal tract, Diffusible Plasma and Liver 2. The default value of 0.0693 is constant for all age ranges. The liver and kidney compartments are excretion pathways. The liver may contain 2-3% of the total body burden of lead, and this is treated as residing in two compartments, Liver1 and Liver2. From Liver1, lead may pass to the small intestine, back to the Diffusible Plasma, or to the Liver2 compartment (Figure 45).

5.2.11 TRANSFER FROM KIDNEY1.

Transfer from the kidney has two components: Kidney 1 from the kidney back to diffusible plasma and Kidney 2 to Kidney 1 from the bladder. Kidney 2 is discussed in Section 5.2.14. The fraction that is transferred from Kidney 1 has a transfer rate constant (0.139/day) for all ages.

5.2.12 TRANSFER FROM BLADDER TO URINE.

Lead that enters the bladder passes from the diffusible plasma to the urinary bladder and from the renal tubules. Estimates of urinary excretion removal halftime in humans *in vivo* are about 0.1 day (Leggett, 1993). The model includes two routes: Diffusible Plasma to the Urinary Path and Diffusible Plasma to Other Kidney Tissue. From the Bladder to Urine is, for all practical purposes, instantaneous.

12. Edi Descrip	t Values for Transfer	from Bladde	er to Urine	Interpolati	on Type: Linear	- L	X Jnit: /day
	A	В	c	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Transfer from Bladder t	12	12	15	11	8	5
<							>
Impor	t Data Export Data	Add Colu	mn Delete Co	lumn	Plot	Save	Cancel

Figure 62: Window for editing the values for the transfer of lead from the bladder to urine. Data are plotted on Figure 63.



Figure 63: Plot of the default values for the transfer of lead from the bladder to urine. This age-dependent variable excretes lead in urine away from the bladder. The values are units per day and are larger for younger children. They become constant at about age fifteen.

5.2.13 TRANSFER FROM LIVER2.

Of the lead leaving the Liver1 compartment, ten percent is assigned to a long term storage compartment, Liver2. The rest is split between the small intestine and diffusible plasma (45% each). This age-dependent variable becomes constant at about age 10 years. The transfer rate is calculated as 0.1×0.0693 /day, for a value of 0.00693 in the younger age groups (Figure 64).

3. Edit Descriptio	Values for Transfer	r from Liver 2 er section 2 by ag	ge range	Interpolatio	n Type: Linear	💌 Uni	t /day
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Transfer from Liver 2	0.00693	0.00693	0.00693	0.00693	0.0019	0.0019
<	-						>
Import I	Data Export Data	Add Colur	mn 🛛 Delete Colu	umn	Plot	Save	Cancel

Figure 64. Age-related coefficients for transfer of lead from Liver2. Data are plotted on Figure 65.



Figure 65. Plot of the transfer of lead from the liver 2 compartment.

5.2.14 TRANSFER FROM KIDNEY2.

This transfer for Kidney 2 represents the amount of lead that passes to the Bladder from the kidney. The numeric values are the same as for Liver 2, and they are age specific (Figure 66).

escrip	otion: Transfer rate to kidne	y section 2 by a	ge range	Interpolati	on Type: Linear	▼ Ur	nit: /day
	A	В	с	D	E	F	
	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
	Transfer from Kidney 2	0.00693	0.00693	0.00693	0.00693	0.0019	0.0019
1							3

Figure 66: Drop down window for editing the default values for the transfer of lead from the kidney. Data are plotted on Figure 67.



Figure 67: Plot of the default values for the transfer of lead from Kidney 2 (From Kidney to Bladder).

5.2.15 TRANSFER FROM FAST SOFT TISSUE.

Soft tissue consists of muscle, fat, and skin. Although this compartment represents a large fraction of the total body mass, the concentration of lead is low. This parameter is substantially larger than other values and is age-related (Figure 68). All of the soft tissue is in one compartment, but the binding capacity within that compartment is not uniform throughout. The model treats the binding capacity as three strengths, although in reality there may be a continuum from low to high. The three strengths are Tenacious Turnover, Intermediate Turnover, and Rapid Turnover. These three turnover rates, with one exception, return the lead to the Diffusible Plasma. The exception is hair, nails and skin, which are described below.

1. Edit Descriptio	Values for Depostio	n Fraction in ead in fast turnov	Fast Soft Tis ver soft tissue	s ue Interpolation	Type: Linear	✓ Unit:	
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Depostion Fraction in Fa	0.08345	0.08345	0.08345	0.08375	0.08375	0.08375
<		1111					>
			1	1			
Import I	Data Export Data	Add Colu	mn Delete Col	umn	Plot	Save	Cancel

Figure 68: Drop down window for editing the default values for the deposition fraction of lead in fast soft tissue. Data are plotted on Figure 69.



Figure 69: Plot of the default values for the deposition fraction of lead in fast soft tissue.

5.2.16 TRANSFER FROM INTERMEDIATE SOFT TISSUE.

This parameter is not age-related and is nearly three orders of magnitude lower than Fast Soft Tissue. Some of this lead is "excreted" in the form of hair, nails, and skin. Note the relative significance of the pathways for the deposition fraction of lead from the body: feces (0.006) urine (0.015, hair/skin/nails (0.4) and sweat (0.0035). The Intermediate Soft tissue excretion fraction (0.4) represents a major excretion route (Hair, Nails, and Skin), two orders of magnitude larger than sweat. This parameter is also not age-dependent.

32. Edit	Values for Depostio	n Fraction in	Intermediat	e Soft Tissue	•		X
Descripti	on: Depostion fraction of I	lead in intermedia	ate turnover s	Interpolation	n Type: Linear	▼ Un	it
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Depostion Fraction in Int	0.01	0.01	0.01	0.01	0.01	0.01
<		Ш					>
Income at	Data Europt Data	Add Calu	and Datata Cal		DI-1	Cause	Concel

Figure 70: Drop down window for editing the default values for the deposition fraction of lead in intermediate soft tissue. Data are plotted on Figure 71.



Figure 71: Plot of the default values for the deposition fraction of lead in fast soft tissue.

5.2.17 TRANSFER FROM SLOW SOFT TISSUE.

This parameter is not age-related and is among the slowest, being more than an order of magnitude slower than Intermediate Soft Tissue.

33. Edi	t Values for Depostio	n Fraction in	n Slow Soft Ti	ssue			×
Descrip	tion: Depostion fraction of I	lead in slow turno	over soft tissu	Interpolation	n Type: Linear	▼ Unit:	
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Depostion Fraction in SI	0.001	0.001	0.001	0.001	0.001	0.001
<		1111					>
Import	t Data Export Data	Add Colu	mn 🛛 Delete Co	lumn	Plot	Save	Cancel

Figure 72: Drop down window for editing the default values for the deposition fraction of lead in slow soft tissue. Data are plotted on Figure 73.



Figure 73: Plot of the default values for the deposition fraction of lead in slow soft tissue.

5.2.18 TRANSFER RATES FROM BRAIN.

The default value for this age-dependent parameter is very small (0.00095/day) and is modeled the same for all age categories. The brain is especially vulnerable to lead, possible due to the role calcium plays in nerve conduction velocity. Although it contains very little lead, brain tissue is an important component of this biokinetic model. The small amount of lead taken up by the brain is strongly retained. The assumed deposition fraction for infants (0.00045) is arbitrarily three times higher than for older children and adults (0.00015), which is based on a removal half-time of two years (Figure 74) (Leggett, 1993).

4. Edit	Values for Depostio	n Fraction in lead in the brain t	Brain by age range	Interpolation	Type: Linear	✓ Unit:	
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Depostion Fraction in Br	0.00045	0.00045	0.00045	0.00015	0.00015	0.00015
<							>
					(
Import [Data Export Data	Add Colu	mn Delete Col	umn	Plot	Save	Cancel

Figure 74: Drop down window for editing the default values for the deposition fraction of lead in the brain. Data are plotted on Figure 75.



Figure 75: Plot of the default values for the deposition fraction of lead in the brain. The scale on the X-axis have been increased by 10³ for presentation purposes to avoid the appearance of a rounding error.

5.2.19 DEPOSITION FRACTION IN URINE.

More than twice the size of the Feces compartment, this fraction is a major elimination route for lead. The coefficient is 0.015 and is constant for all age groups. This number represents the transfer rate from the diffusible plasma to the urinary bladder. The Deposition Fraction for lead in Feces is 0.006, less than half that of Urine, and represents the lead entering the digestive tract from the mucociliary lift mechanism in the Lungs and the passage of lead in bile from the liver. The Sweat deposition fraction 0.0035 represents the small amount of the lead that is eliminated by the sweat pathway, and is not age-dependent. The value is somewhat arbitrary (0.0035) in that it is assumed that 0.35% of the lead leaving the diffusible plasma does so by the sweat route.

5.2.20 DEPOSITION FRACTION IN BONE.

This age-dependent variable represents deposition to bone surfaces and varies widely in children and young adults, up to age 25, after which it reaches a low, constant value of 0.08 (unitless), about one third that of a newborn child (Figure 76). This means that, for children under 20, the history of lead exposure, as released from their bone tissue, can be as important a contributor to blood lead as their current exposure.

23. Edit	Values for Deposition	on Fraction ii	n Bone				×
Descripti	on: Deposition fraction in	bone by age rang	je	Interpolation 1	Type: Linear	▼ Unit:	
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Deposition Fraction in B	0.24	0.24	0.144	0.128	0.179	0.237
<		1111					>
Import	Data Export Data	Add Colur	mn Delete Coli	Inn	Plot	Save (Cancel

Figure 76. Deposition fraction in Bone. Data are plotted on Figure 77.



Figure 77. Plot of Deposition Fraction in Bone. Given that the deposition of lead appears to be similar to the deposition of alkaline earth elements calcium, strontium, and barium (Leggett, 1993), it follows that the spike at age 15, which probably represents a spurt in bone growth, would occur also for lead.

Parame	eters	Values	Units	Descriptions
Trabecular Bone Depositio	n Fraction	Edit		Fraction of bone deposition going to trabecular by age range
5 Deposition Fraction in Liver	r 1	0.04		Deposition fraction in liver section 1
6 Depostion Fraction from Liv	ver 1 to Liver 2	0.1		Depostion fraction of lead from liver section 1 to section 2
7 Depostion Fraction from Liv	ver to Small Intestine	0.45		Depostion fraction of lead from liver to small intestine
8 Depostion Fraction from Liv	ver to Plasma	0.45		Depostion fraction of lead from liver to plasma
9 Deposition Fraction in Kidne	ey 1	0.02		Depostion fraction of lead in kidney section 1
Depostion Fraction in Kidne	ey 2	0.0002		Depostion fraction of lead in kidney section 2
31 Depostion Fraction in Fast	Soft Tissue	Edit		Depostion fraction of lead in fast turnover soft tissue by age range
2 Depostion Fraction in Interr	nediate Soft Tissue	Edit		Depostion fraction of lead in intermediate turnover soft tissue by
33 Deposition Fraction in Slow	Soft Tissue	Edit		Depostion fraction of lead in slow turnover soft tissue by age rar
34 Deposition Fraction in Brain	1	Edit		Depostion fraction of lead in the brain by age range
35 Depostion Fraction in RBC		0.24		Depostion fraction of lead in red blood cells
36 Deposition Fraction in EVF		0.5		Depostion fraction of lead in the EVF
37 Size of EVF Relative to Pla	sma	3		Relative size of EVF to plasma
38 Transfer Rate from Plasma	1	2000	/day	Transfer rate from plasma
39 Deposition Fraction for Plas	sma Proteins	0.0004		Deposition fraction for plasma proteins
10 Rate of Loss For Plasma P	roteins	0.139	/day	Rate of loss for plasma proteins
11 Transfer Rate From RBC		Edit	/day	Plasma transfer rate from RBC by age range
12 RBC Reference Volume		22	dl	Amount of RBC in reference adult male
13 Plasma Reference Volume		30	dl	Amount of plasma in reference adult male
44 Amount of Blood		Edit	dl	Amount of blood by age range
15 Chelation Factor 1		0.4		Chelation Factor 1
6 Chelation Factor 2		1243		Chelation Factor 2
a louise reaso		4050	1	

Figure 78. Bottom of biokinetic Parameters Window. Parameter 47 is Chelation Factor 3 and has a value of 1250.

5.2.21 TRABECULAR BONE DEPOSITION FRACTION.

The fraction of deposition going to trabecular bone ranges from 0.20 to 0.56 (Figure 79), increasing with age. As the bone grows, lead is incorporated into the crystalline structure of the bone tissue and replaces calcium. Later, as new bone tissue replaces old, the lead is released from the nonexchangeable trabecular volume and passes back to the Diffusible Plasma.

4. Edit	Values for Trabecu	lar Bone Dep	osition Fracti	on			
Descripti	on: Fraction of bone depo	osition going to tra	abecular by a	Interpolation	Type: Linear	▼ Unit:	
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Trabecular Bone Deposi	0.2	0.2	0.2	0.222	0.25	0.279
<		1111					>
			la con	-1			I

Figure 79. Age-dependent variable for Trabecular Bone Deposition Fraction. Data are plotted on Figure 80.



Figure 80. Trabecular bone deposition fraction.



Figure 81: Plot of the default values for the transfer of lead from the liver.

5.2.22 DEPOSITION FRACTION IN COMPARTMENTS WITHOUT AGE VARIABILITY.

This is the fraction of lead (0.04) deposited in Liver1 from Diffusible Plasma, and does not vary with age. This Liver 1 Fraction receives 4% of the lead released by Diffusible Plasma, giving a transfer rate of 80/day and a removal halftime of ten days. Ten percent (0.1) of the Liver 1

fraction moves to the second liver component, Liver 2. Forty-five percent (0.45) of the Liver 1 fraction is deposited in the Small Intestine through the bile duct. Most of this is eliminated with feces; a small amount may be reabsorbed into the Diffusible Plasma.

All of the lead from Liver 2 is passed back to the Diffusible Plasma compartment. Most of the lead in Liver1 passes into the bile duct and then to the small intestine. The rest (0.45) passes back to the Diffusible Plasma. The fraction deposited in Kidney 1 from Diffusible Plasma (0.02) is also called the Urinary Path (UP). It has a high deposition but a shorter retention than Kidney 2. The fraction deposited in Kidney 2 from Diffusible Plasma (0.0002) is also called "Other Kidney Tissue (OKT). Although it has a lower uptake than Kidney 1, Kidney 2 has a much larger retention time. Lead passes out to the Diffusible Plasma compartment from the OKT.

5.2.23 DEPOSITION FRACTION FROM DIFFUSIBLE PLASMA TO FAST SOFT TISSUE.

This fraction deposited in Fast Soft Tissue from Diffusible Plasma varies slightly with age and has a value from 0.08345 to 0.08875 for a transfer rate of 167-177/day. (Transfer Rate = Deposition Fraction x 2000/day) (Figure 82).

31. Edit	\$1. Edit Values for Depostion Fraction in Fast Soft Tissue									
Description: Depositon fraction of lead in fast turnover soft tissue Interpolation Type						▼ Unit:				
	A	В	с	D	E	F	G			
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000			
2	Depostion Fraction in Fa	0.08345	0.08345	0.08345	0.08375	0.08375	0.08375			
<										

Figure 82. Age-dependent variable for Diffusible Plasma to Fast Soft Tissue. Data are plotted on Figure 83.



Figure 83. Plot of the default values for the Diffusible Plasma to Fast Soft Tissue parameter.

5.2.24 DEPOSITION FRACTION FROM DIFFUSIBLE PLASMA TO INTERMEDIATE SOFT TISSUE.

This is the fraction deposited in Intermediate Soft Tissue from Diffusible Plasma with a turnover rate of 25 to 300 days and includes a component of permanent elimination (Hair, Nails, and Skin). This compartment has a Deposition Fraction from 0.005, giving a transfer rate of 0.00277/day to HNS and 0.00416 to Diffusible Plasma (Figure 84).

32. Edit	2. Edit Values for Depostion Fraction in Intermediate Soft Tissue									
Descripti	on: Depostion fraction of I	lead in intermedia	ate turnover s	Interpolation	Type: Linear	▼ Unit:				
	A	В	с	D	E	F	G			
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000			
2	Depostion Fraction in Int	0.01	0.01	0.01	0.01	0.01	0.01			
<							>			
laun aut	Import Data Export Data Add Column Delete Column Plot Save Cancel									

Figure 84. Age-dependent variable for Diffusible Plasma to Intermediate Soft Tissue. Data are plotted on Figure 85.



Figure 85. Plot of the default values for the diffusible Plasma to Intermediate Soft Tissue parameter.

5.2.25 DEPOSITION FRACTION FROM DIFFUSABLE PLASMA TO SLOW SOFT TISSUE.

A small amount of lead is transferred from Diffusible Plasma to Slow Soft Tissue with a turnover rate from 1500 to 10000 days in Slow Soft Tissue. The removal half time is 5 years and transfer rate 2/day in, and 0.00038/day out of the compartment.

3. Edit Values for Depostion Fraction in Slow Soft Tissue									
Descripti	on: Depostion fraction of	lead in slow turno	over soft tissu	Interpolation	Type: Linear	▼ Unit			
	A	В	с	D	E	F	G		
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000		
2	Depostion Fraction in SI	0.001	0.001	0.001	0.001	0.001	0.001		
<							>		
Import	Import Data Export Data Add Column Delete Column Plot Save Cancel								

Figure 86. Age-dependent variable for Diffusible Plasma to Slow Soft Tissue. Data are plotted on Figure 87.



Figure 87. Plot of the default values for the Slow Soft Tissue Parameter.

5.2.26 DEPOSITION FRACTION FROM DIFFUSIBLE PLASMA TO BRAIN.

This is the fraction deposited in Brain from Diffusible Plasma, and is age-dependent.

escripti	on: Depostion fraction of I	lead in the brain l	by age range	Interpolation	Type: Linear	▼ Un	it
	A	В	с	D	E	F	
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Depostion Fraction in Br	0.00045	0.00045	0.00045	0.00015	0.00015	0.0001:
<1							

Figure 88. Age-dependent variable for the Deposition Fraction from Diffusible Plasma to Brain. Data are plotted on Figure 89.



Figure 89. Plot of the default values for the Deposition Fraction from Diffusible Plasma to Brain parameter.

5.2.27 COMPONENTS OF THE CIRCULATORY SYSTEM.

The fraction of lead in Bound Plasma that is transferred to Red Blood Cells is the Deposition Fraction. Red blood cells have a high affinity for lead, and therefore a high percentage of the lead in blood is bound to the red blood cells. These cells, which have a half life of about 120 days may also serve as a mechanism for eliminating lead from the system. The fraction of lead that is deposited in Extra Vascular Fluids from Diffusible Plasma and Red Blood Cells is 0.5, or 50%. The Ratio of Total Extravascular Fluid mass to mass of blood Plasma is 3. The transfer of lead from diffusible plasma is on the order of 1-2/min, or 2000/day. The deposition fraction for Plasma Proteins is extremely small, 0.0004. The rate of loss for Plasma Proteins is 0.139/day; approximately 14 % of the lead in plasma proteins is transferred to other compartments per day. The amount of Red Blood Cells in the Reference Adult Male is 22 deciliters. The amount of Blood Plasma in the Reference Adult Male is 30 deciliters.

5.2.28 TRANSFER RATE FROM RBC TO PLASMA.

Age-dependent variable that represents the daily transfer of lead from Red Blood Cells to Plasma are shown in Figure 90.

41. Edit Descripti	Values for Transfer	Rate From R	BC	Interpolation	Type: Linear	Unit:	/day
	A	В	с	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Transfer Rate From RB	0.462	0.462	0.462	0.277	0.139	0.139
<						>	
Import	Data Export Data	Add Colu	umn	Plot	Save	Cancel	

Figure 90. Age-dependent variable for the Transfer Rate from RBC. Data are plotted on Figure 91.



Figure 91. Plot of the default values for the Transfer Rate for RBC parameter.

5.4.29 AMOUNT OF BLOOD.

Total volume, in deciliters, of blood, which is age-dependent, is shown in Figure 92.

44. Edi	t Values for Amount	t of Blood					
Descript	tion: Amount of blood by	age range		Interpolatio	on Type: Linear	• Ut	nit: dl
	A	В	c	D	E	F	G
1	Age Range	0.000	0.274	1.000	5.000	10.000	15.000
2	Amount of Blood	3	5.5	8	17	30	50
<		1111					>
	- 1 1		1				(
Import	t Data Export Data	Add Col	umn Delete Co	olumn	Plot	Save	Cancel

Figure 92. Age-dependent variable for total volume of blood (deciliters) Data are plotted on Figure 93.



Figure 93. Plot of the default values for the Amount of Blood parameter.

5.2.30 CHELATION FACTORS 1, 2, and 3.

These chelation factors are not currently used. They were a part of the original Leggett model because chelation data were used at that time. The options are included against the possibility they might have some future use.

SECTION 6. CONTROLLING THE MODEL OUTPUT

6.1. A SIMPLE CONFIDENCE TEST

At this point, it is helpful to begin developing your familiarity with the model. First, observe the model output in the case where no changes have been made to the model defaults. Restart the model and select the option to "Create and new study." Click on "Youth" to populate the age ranges through 13 years of age. Click Run, then OK to save the file; specify a new file named First Flight, and save this file; click "Exposure" and "Absorption," then "Continue" "Start" and "Continue" to run the model with all default settings. Select "Plot" to produce the graph shown in Figure 94. This is a plot of the predicted blood lead concentration trough 4700 days, or thirteen years. Note that when only one parameter (Blood Pb) is plotted, the X and Y values are given as well.



Figure 94. Output with all default settings

Now return to the opening window by closing out two windows in succession, and change the exposure by selecting Media Exposure and clicking Edit. Click edit again to get to the parameter array for Air. Click the Dust tab and change the lead concentration of dust ingested by a toddler, preschooler and youth from the default 200 μ g Pb/g to 2000 μ g Pb/g. Click OK and run the model as before. You should get:



Figure 95. Total Blood lead with elevated dust lead exposure.

Note that the Y-axis scale has changed to accommodate the range of values. While it appears that the effect of the increased exposure lasts only through 1900 days, the upper bound for the youth age range, a closer look will show that the blood lead concentrations will remain elevated throughout the next age range, adolescent, as the lead is recycled through other body tissues, including bone tissue, returning sometime later to the blood. Indeed, even at age nineteen, six years after the secession of the elevated lead exposure, the blood lead for the elevated exposure child is still about 15% higher than the child with the lower exposure.

Take a moment to look at other body tissues to become familiar with the way the model distributes the lead over the various compartments throughout the lifetime of the individual. Remember that, as an experienced modeler, you will become adept at tweaking the biokinetic parameters according to information that you may have from your own studies, or to test the reality of data published by others.

After editing the biokinetic parameters, close this option [X] and click RUN from the Primary window. The Model Run window appears next and offers you a choice of Exposure, Absorption, and Biokinetics. Click on any of these for which you want to see the output. Clicking on Absorption automatically triggers Biokinetics. Click CONTINUE and the Exposure Model Processing window appears. Click START. The processing time depends on the complexity of your inputs. The three options for historical diet, air and dust lead concentrations are under development and are not functional at this time.

There is an option on the Exposure Model Processing window to select the exposure model output file to view or bypass this step. This is fully optional and is if value if you are troubleshooting the model. The *.mod file gives you a view of the stepwise values you have selected. This file can contain over 32,000 lines and may take a few moments to generate. The *.mpa file gives you the total exposure μ Pb for each age range. Both of these files may be exported to a spreadsheet.

The next window allows you to view the Exposure model output files using different extensions for the same file name that you have designated.

Filename.mpa (Exposure Model Output by Age) gives the total amount of lead to which the modeled person was exposed by the six different routes. All of these should be listed by the age category.

6.2. BIOKINETIC OUTPUTS

6.2.1 TABLES

CLOSE on the Select window runs the model. The outputs are saved to the same location as that for the input files. In the example below, the output is for 100 day intervals as stipulated in the Primary window, Simulation time step 100 day(s).

An explanation of the various outputs and units are given in Table 2.

6.1.2 PLOTTING

Click PLOT on the window generates a default color plot of several parameters versus days.

PRINT provides a copy of the chart.



Figure 96. Using the Built in Graphic Feature.

6.1.3 EDITING THE CHART

Click EDIT allows for several options.

In the Editing Window, SERIES, ticking (or unticking) a box will automatically adjust the chart.

a. Axis: (i) Unticking Automatic allows changes to the scales. For example, the Maximum and Minimum ranges can be changed for the Y-(Left) axis and the number of decimal places.

(ii) <u>NOTE:</u> The number of decimal places **cannot** be changed on the X- (Bottom) axis.

(iii) Titles can be added to the Y-axis and also the X-axis. However, the default title for X-axis ('Days' in blue color) remains. This can be removed by going into Editing Titles, scrolling down from TITLES to FOOT and deleting 'Days' (Fig. 39).

(iv) The frequency of Ticks can be changed.

(v) The main title can be edited in TITLES. For example, change Font size and Back Color. However, even though the title states it is black, the Black color needs to be clicked to achieve this (otherwise it stays yellow).

(v) The background blue can be changed to white in PANEL by making START and END colors white.
(vii) To change output style for each parameter, for example, YCORT, go to EDITING, CHART, SERIES, click on the color square, in Color Window change to desired color, click OK, click on BORDER, in the window Border Color Editor, change Style, Click OK.

Change other parameters in the same way. PRINT chart.

7.0 REFERENCES

- Leggett, R.W. (1992b) A generic age-specific biokinetic model for calcium-like elements Radiat. Prot. Dosim. 41: 183-198.
- Leggett, R.W. (1993a) An age-specific kinetic model of lead metabolism in humans. Environ. Health Perspect. 101: 598-615.
- Mickle, M.H. (1998) Structure, use, and validation of the IEUBK model. Environ. Health Perspect. Suppl. 106: 1531-1534.
- Mushak, P. (1991) Gastro-intestinal absorption of lead in children and adults: overview of biological and biophysico-chemical aspects. Chem. Speciat. Bioavail. 3: 87-104.
- O'Flaherty, E.J. (1998) Physiologically based kinetic model for lead in children and adults. Environ. Health Perspect. Suppl. 106: 1495-1503.
- O'Flaherty, E.J.; Inskip, M.J.; Franklin, C.A.; Durbin, P.W.; Manton, W.I.; Baccanale, C.L. (1998) Evaluation and modification of a physiologically based model of lead kinetics using data from a sequential isotope study in Cynomolgus monkeys. Toxicol. Appl. Pharmacol. 149: 1-16.
- Oreskes, N. (1998) Evaluation (Not Validation) of Quantitative Models. Environ. Health Perspect. 106(Suppl 6): 1453-1460.
- Pounds, J.G.; Leggett, R.W. (1998) The ICRP age-specific biokinetic model for lead: validations, empirical comparisons, and explorations. Environ. Health Perspect. Suppl. 106: 1505-1511.
- U.S. EPA (1994a) Guidance Manual for the Integrated Exposure Uptake Biokinetic Model for Lead in Children. U.S. Environmental Protection Agency. EPA/540/R-93/081, PB93-963510.
- U.S. EPA (1994b) Technical Support Document: Parameters and Equations Used in Integrated Exposure Uptake Biokinetic Model for Lead in Children (v0.99d). U.S. Environmental Protection Agency EPA/540/R-94/040, PB94-963505.
- White, P.D.; Van Leeuwen, P.; Davis, B.D.; Maddaloni, M.; Hogan, K.A.; Marcus, A.H.; Elias, R.W. (1998) The conceptual structure of the integrated exposure uptake biokinetic model for lead in children. Environ. Health Perspect. 106(Suppl 6): 1513-1530.